BAIRD Biopsy

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ENGLISH

INSTRUCTIONS FOR USE

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Device Description

The ENCOR ENSPIRE® Breast Biopsy System provides control operations for specialized biopsy instruments intended to acquire tissue samples of suspected breast abnormalities. The ENCOR ENSPIRE® Breast Biopsy System is designed to recognize individual Biopsy Drivers, provide a user interface, and to accept inputs when used for specific procedures. Each type of Biopsy Driver performs unique functions.

The ENCOR ENSPIRE® Breast Biopsy System may be used with ENCOR® and ENCOR® MRI Drivers, ENCOR® Foot Pedals, and ENCOR® Probes. Reference **Figures 1, 2 and 3.** ENCOR® Biopsy Probes, Tubing Cassettes, Canisters, Adapters, and Needle Guides are available separately.


Figure 1 – ENCOR® Biopsy Probes

Figure 2 – ENCOR® Driver and Probe

Figure 3 – ENCOR® Foot Pedal

Indications for Use

The ENCOR ENSPIRE® Breast Biopsy System is indicated to provide breast tissue samples for diagnostic sampling of breast abnormalities.

- It is intended to provide breast tissue for histologic examination with partial or complete removal of the imaged abnormality.
- It is intended to provide breast tissue for histologic examination with partial removal of a palpable abnormality.

The extent of a histologic abnormality cannot always be readily determined from palpation or imaged appearance. Therefore, the extent of removal of the palpated or imaged evidence of an abnormality does not predict the extent of removal of a histologic abnormality, e.g., malignancy. When the sampled abnormality is not histologically benign, it is essential that the tissue margins be examined for completeness of removal using standard surgical procedures.

In instances when a patient presents with a palpable abnormality that has been classified as benign through clinical and/or radiological criteria (e.g. fibroadenoma, fibrocystic lesion), the ENCOR ENSPIRE® Breast Biopsy System may also be used to partially remove such palpable lesions. Whenever breast tissue is removed, histological evaluation of the tissue is the standard of care. When the sampled abnormality is not histologically benign, it is essential that the tissue margins be examined for completeness of removal using standard surgical procedures.

Contraindications

1. This device is not intended for use except as indicated.
2. The ENCOR ENSPIRE® Breast Biopsy System is contraindicated for those patients where, in the physician's judgment, there is an increased risk of complications associated with percutaneous removal of tissue samples.

Warnings

1. The ENCOR ENSPIRE® Breast Biopsy System must be properly grounded to ensure patient safety. The system is supplied with a medical grade Power Cord with AC plug. Do not connect the included Power Cord to extension cords or three-prong to two-prong Adapters. To avoid the risk of electric shock, this equipment must only be connected to supply mains with protective earth.
2. To minimize interference with other equipment, cables should be positioned in such a manner to prevent contact with other cables.
3. Use of accessories not compatible with the ENCOR ENSPIRE® Breast Biopsy System may create potentially hazardous conditions.
4. Only use ENCOR® and ENCOR® MRI Drivers with script version 1.19 or greater with the ENCOR ENSPIRE® Breast Biopsy System. The system is not compatible with earlier driver scripts. The script version is identified on the touch screen display during system initialization.
5. The ENCOR ENSPIRE® Breast Biopsy System console may not be placed in an MRI suite. Place the console outside of the MRI suite and use the appropriate ENCOR® MRI accessories when performing a biopsy under MRI guidance.
6. No modification of this equipment is allowed. Do not remove the ENCOR ENSPIRE® Breast Biopsy System housing. Removal of the housing may cause electrical shock.
7. The ENCOR ENSPIRE® Breast Biopsy System is not classified as an AP or APG classified device. The system is not suitable for use in the presence of flammable anesthetic.
8. Do not use in the presence of infection.
9. After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practice and applicable local, state, and federal laws and regulations.

Precautions

1. This device should only be used by physicians trained in percutaneous biopsy procedures.
2. Carefully inspect the device prior to use to verify that the device has not been damaged. Do not use if product damage is evident and/or needle is bent.

3. Locate the ENCOR ENSPIRE® Breast Biopsy System as far as possible from other electronic equipment to minimize interference or degradation of performance of the ENCOR ENSPIRE Breast Biopsy System.
4. Inspect accessories and cords for breaks, cracks, nicks, or other damage before every use. If damaged, do not use. Failure to observe this precaution may result in injury or electrical shock to the patient or to the operator.
5. Inspect tubing connections to the Vacuum Canister and the Vacuum Tubing Cassette to ensure proper vacuum levels are achieved and maintained during use.
6. Inspect the Vacuum Canister to ensure the lid is secure and that no damage has occurred during shipping or installation. A heavily scratched canister can break during use.
7. Do not leave the ENCOR ENSPIRE® Breast Biopsy System powered on overnight. Damage may occur to the Vacuum or the Vacuum and Rinse Tubing Cassette.
8. Connect the Power Cord to a hospital grade wall outlet having the correct voltage or product damage may result.
9. Patients who may have a bleeding disorder or who are receiving anticoagulant therapy may be at increased risk of complications.
10. As with any biopsy instrument, there is a potential for infection.
11. All breast biopsies should be performed under imaging guidance to confirm the Probe position relative to the target region to be sampled and to help mitigate the occurrence of a false negative biopsy.
12. When performing a biopsy with ENCOR® and ENCOR® MRI Probes, the orientation of the sample notch is dictated by the image guidance selected. Prior to initiating the procedure, confirm that the sample notch orientation is correct for the image guidance being used.
13. Ensure that the ENCOR ENSPIRE® Breast Biopsy System is positioned in such a way that the Power Cord and retainer are accessible. In the event that the system Power Switch is inoperable, release retainer and remove cord to shut off system power.

Potential Complications

Potential complications may include, but are not limited to hematoma, hemorrhage, infection, adjacent tissue injury, pain, allergic reaction, and tissue adherence to the Biopsy Probe during removal from the breast (as per routine biopsy procedures, it may be necessary to cut tissue adhering to the stylet or coaxial while removing it from the breast).

Equipment Required

The following equipment is required for a biopsy procedure:

- Appropriate imaging modality and accessories
- ENCOR® or ENCOR® MRI Driver
- ENCOR® or ENCOR® MRI Probe
- ENCOR ENSPIRE® Vacuum Tubing Cassette or ENCOR ENSPIRE® Vacuum and Rinse Tubing Cassette
- Vacuum Canister
- Saline (optional)
- Tissue Marker (optional)
- Surgical gloves and drapes
- Local anesthetic
- Scalpel
- Other equipment as necessary

Directions for Use

Unpacking Instructions

Thoroughly inspect the packaging materials and the ENCOR ENSPIRE® Breast Biopsy System for any signs of transit damage. Do not use the unit if it appears to have any damage. Contact Bard Customer Service for instructions on how to report shipping damage.

The ENCOR ENSPIRE® Breast Biopsy System is delivered assembled with the exception of the touch screen which is detached and packaged separately. Carefully remove the main body and touch screen from the packaging. Unlock the release knob on the touch screen as illustrated in **Figure 4**. Align the connecting track on the back of the display screen and slide into place on the arm of the main body as illustrated in **Figure 5**. Attach the power, ground and communication cables from the main body to the touch screen. Attach the ground cable to the ground stud and hand tighten the supplied nut. Use the supplied flat head screwdriver to tighten the screws on the communication cable. Reference **Figure 6** for an illustration of proper communication and power cable installation locations.


Figure 4 – Release knob on the touch screen

Figure 5 – Attaching the touch screen to the main body


Figure 6 – Connecting the touch screen cables
Note: Actual connector style may differ from illustration

Procedure Setup

Prior to transporting the system, unlock the system casters and lower the tray and monitor to their lowest positions. Align the tray and monitor to face the front of the system. The monitor and tray should be positioned as illustrated on the cover of this Instructions for Use. When transporting the system, push or pull using the handles located on the front and back of the tray. Do not transport the system by pushing or pulling sideways. To transport the system outside of normal hospital conditions (e.g. outdoors or to another facility), use Bard supplied packaging.

1. Position the ENCOR ENSPIRE® Breast Biopsy System as required for safe use and plug the Power Cord (located on the rear of the console) into a hospital grade wall outlet.
 - a. The system includes casters with locking brakes. After the system is positioned, lock the casters as illustrated in **Figure 7**.

- b. The touch screen height may be adjusted by releasing the locking lever on the main post and pressing the release button on the top of main body pole as illustrated in **Figure 8**. The touch screen tilt may be adjusted by loosening the lock behind the touch screen and moving the screen to the desired position as illustrated in **Figure 9**.
- c. The adjustable tray can accommodate two Drivers as well as Driver cables and Probe tubing. The position of the tray may be adjusted by releasing the locks on the tray and rotating, lifting, or lowering the tray to the desired height as illustrated in **Figure 10**.

Figure 7 – Locking the system casters


Figure 8 – Adjusting the touch screen height


Figure 9 – Adjusting the touch screen tilt


Figure 10 – Adjusting the tray position

2. Install the Vacuum Canister and Tubing Cassette per the directions provided in the ENCOR ENSPIRE® Vacuum and Rinse Tubing Cassette Instructions for Use.
3. Turn on the main power via the Power Switch on the rear of the console as illustrated in **Figure 11**.

Figure 11 – Main Power Switch

4. Turn on the standby power via the switch on the front of the system as illustrated in **Figure 12**.

Figure 12 – Standby Power Switch

Figure 13 – Driver Connectors

5. Connect the appropriate Driver to the rectangular connector on the rear of the console illustrated in **Figure 13**. The active Driver is connected by a green LED light above the connector. Reference the applicable Driver Instructions for Use for further information on Driver setup. If desired, connect the appropriate Foot Pedal to the connector on the rear of the console as illustrated in **Figure 14**.

Figure 14 – Foot pedal connector

6. Follow the touch screen prompts. Install and calibrate the Biopsy Probe. Reference the applicable Probe Instructions for Use for further information on Probe setup.

Operating the System via the Touch Screen

7. Complete the biopsy procedure per the 'Directions for Use' section in the Driver Instructions for Use.
8. The ENCOR ENSPIRE® Breast Biopsy System status is displayed on the touch screen as shown in **Figure 15**. The sample acquisition options may be selected via this screen.
 - a. To change the sampling pattern, place a finger on desired clock position within the circle on the touch screen illustrated in **Figure 15**.
 - b. To change the sampling from full to half, touch the "Half" (sample) button on the touch screen as illustrated in **Figure 15**. Touch the "Full" (sample) button to return to full sampling.
 - c. To change the Sample Mode from normal tissue to dense tissue, touch the "Dense" (tissue) button on the touch screen as illustrated in **Figure 15**. Touch the "Normal" (tissue) button to return to normal tissue sampling.

Half (Sample button)

Full (Sample button)

Dense (Tissue button)

Normal (Tissue button)

9. During the procedure, the ENCOR ENSPIRE® Breast Biopsy System has an optional Rinse Mode when using the Vacuum and Rinse Tubing Cassette and saline. This mode allows the operator to apply a saline rinse to samples collected in the sample trap. Initiate the rinse by pressing the "Rinse" button on the touch screen.

Rinse (button)

10. If needed, the ENCOR ENSPIRE® Breast Biopsy System allows the operator to administer additional anesthesia through the ENCOR® Biopsy Probe. Enter Anesthetic Mode by pressing the "Anesthetic" tab on the touch screen as illustrated in **Figures 15 and 16**. Select the desired anesthetic option to sweep the sample notch through the desired pattern while administering anesthesia. Leave Anesthetic Mode by pressing the "Sample" button on either the handswitch or the Foot Pedal.


Figure 15 – The touch screen – Sample Tab

Figure 16 – The touch screen – Anesthetic Tab

11. Following the procedure, the ENCOR ENSPIRE® Breast Biopsy System allows the operator to deploy a marker through the Probe. Enter Marker Mode by pressing the "Marker" tab on the touch screen as illustrated in Figure 17. Deploy the marker per the marker Instructions for Use. Rotate the sample notch 180° by pressing the "180°" button on the touch screen as illustrated in Figure 17. Exit Marker Mode by pressing the "Sample" button on either the handswitch or Foot Pedal.

180° button


Figure 17 – The touch screen – Marker Tab

After Use

12. Remove disposables and discard per your facility's procedures. Clean and disinfect the system hardware per the 'Cleaning and Maintenance' section.
Note: The system may be left on if additional procedures are to be performed; it is not recommended to leave the main power on overnight. It is not recommended to unplug the Driver after each use. The Driver and Foot Pedal may be left connected to the ENCOR ENSPIRE® Breast Biopsy System if all procedures are to be performed with the same Driver.
13. To set the system to standby, press and hold the switch on the front of the console (illustrated in Figure 12) until the system shuts down. To turn off the main power, press the main Power Switch located on the rear of the console (illustrated in Figure 11).

Cleaning and Maintenance

To clean the ENCOR ENSPIRE® Breast Biopsy System, turn off the power and disconnect the Power Cord from the AC outlet. Leave all other cables connected. Clean all exposed surfaces of the ENCOR ENSPIRE® Breast Biopsy System, cables, and Foot Pedal with a soft, lint-free cloth dampened with Dispatch or a similar cleaning solution or disinfectant. Follow the procedures approved by your facility or use a validated infection control procedure. Always follow the instructions provided by the manufacturer of the cleaner or disinfectant used.

Do not use abrasive cleaners or spray any fluids directly onto any part of the ENCOR ENSPIRE® Breast Biopsy System or cables. Spraying the system may cause it to malfunction and will void the warranty.

Do not submerge any component of the ENCOR ENSPIRE® Breast Biopsy System in fluids. Submerging in fluids may cause the system to malfunction and will void the warranty. Sterilization and exposure to fluids may damage the electrical components of the device. If the system is inappropriately cleaned, it may cause the system to malfunction and will void the warranty. Do not autoclave any components of the ENCOR ENSPIRE® Breast Biopsy System. Do not heat over 54°C (129°F).

Disposal of Equipment

After following the cleaning recommendations above, there are no biohazard risks involved with the disposal of the ENCOR ENSPIRE® Breast Biopsy System or any of the reusable accessories at the end of the ENCOR ENSPIRE® Breast Biopsy System's useable life per Waste of Electrical and Electronic Equipment (WEEE) directive [Directive 2002/96/EC]. Dispose of any single use device per the Instructions for Use for the disposable device.

Do not dispose of the ENCOR ENSPIRE® Breast Biopsy System by placing into trash receptacles.

Routine Maintenance

No calibrations or adjustments are necessary for the ENCOR ENSPIRE® Breast Biopsy System outside of normal use. In addition, the following routine maintenance activities may be completed:

- Electrical safety testing should be performed at intervals no greater than 12 months, using a standard medical safety analyzer.
- Regularly inspect all cables and connectors for possible wear or damage.

Odor Control Filter Replacement

Turn the Power Switch to off. Remove the old Odor Control Filter by rotating it counter-clockwise and discard. Install the new Odor Control Filter (F3000) by rotating clockwise until snug. The location of the Odor Control Filter is identified in Figure 18.


Figure 18 – Odor Control Filter

Fuse Replacement

The Fuse Holder is accessible after the Power Cord is removed. Unplug the ENCOR ENSPIRE® Breast Biopsy System from the wall outlet, and remove the Power Cord from the power inlet connector on the rear panel of the console. Using a small slotted screwdriver, press the release tab located near the center of the Fuse Holder. Pull the Fuse Holder assembly outward to remove it. Install the Fuse Holder. Remove and replace any blown fuses in the Fuse Holder assembly and snap it back into place. The Fuse Holder orientation determines the voltage. Rotate the fuse 180 degrees for 230V orientation. Replace the Fuse Holder assembly and reconnect the Power Cord to the power inlet connector.

Recommended Fuse for Replacement:

Series: 326P (3AB)
Rated: 250V, 10 Amp, Slo-Blo
Certification: UL 248-14
(For example; Littelfuse P/N: 0326010.HXP)

Troubleshooting

The ENCOR ENSPIRE® Breast Biopsy System will display any alert conditions that occur. If the system can identify the cause of the alert, the information is displayed on the touch screen.

Depending on the cause of the alert, different actions may be required to clear the alert. Certain alert conditions, when corrected, allow the system to recover automatically and return to a ready state. Other alerts require pressing a button on the touch screen to clear. If a Driver is connected and is the cause of the alert, follow the Driver Instructions for Use to clear the alert.

When no other actions are indicated, clear the alert condition by turning off the main power using the Power Switch on the rear of the system (illustrated in Figure 11). Wait 20 seconds, turn the main power back on, and restart the ENCOR ENSPIRE® Breast Biopsy System with the standby button on the front of the console (illustrated in Figure 12).

If any alert condition persists after correcting known problems and clearing the alert condition from the screen, contact Bard for service.

Table 1: Error Status and Troubleshooting		
Error Screen Text	Error Number(s)	Troubleshooting
System Error	001 – 002, 005 – 071, 073, 103 – 118, 120 – 121, 123 – 124, 126 – 128, 136, 139	A system error occurred. Contact Bard for service.
System Error	003 – 004, 072, 074 – 102, 119, 122, 125, 129, 131 – 135, 137 – 138, 140 – 150, 153 – 163	A system error occurred. Restart the system. If the problem persists, contact Bard for service.
Failure Reading Handpiece Memory	130	An error with the Driver/handpiece occurred. Restart the system. If the problem persists, contact Bard for service.
System Error	151, 152	A button on the Foot Pedal is stuck. Check the Foot Pedal and restart the system. If problem persists, contact Bard for service.
Control Module API Initialization Error	1001, 2001	An error occurred during initialization. Restart the system. If the problem persists, contact Bard for service.
Failed Vacuum – Initialization Error	1002, 2002	An error occurred during initialization. Restart the system. If the problem persists, contact Bard for service.
Invalid Gearbox I.D./Driver Type	1003	An error with the Driver/handpiece occurred. Restart the system. If the problem persists, contact Bard for service.
Probe Not Removed	1004, 2003	The Probe is not removed. Remove the Probe to continue.
Driver Button or Contact Stuck during Initialization	1005, 2004	A button on the Driver is stuck. Check the Driver and restart the system. If problem persists, contact Bard for service.
Driver Error – Initialization Failed	1006 – 1007, 2005 – 2006	An error with the Driver occurred. Disconnect the Driver to continue. Connect Driver to retry. If the problem persists, contact Bard for service.
Driver Error – Control Module Software Failure	1008	The system software is out of date. Contact Bard for service to update the system software.
Driver Not Connected	1009, 2007	Driver is not connected. Connect Driver to continue.
System Power Down Detected	1010, 2008	The system is powering down. Please wait.
Calibration Failed	1011, 2009 – 2010	Calibration failed. Re-calibrate Probe or remove Probe to continue.
Probe Error - Wrong Probe inserted	1012	MRI Probe must be used with MRI Driver. Non-MRI Probe must be used with Non-MRI Driver. Remove Probe to continue.
Restoring Vacuum	1013, 2012	Restoring Vacuum. Please wait.
Vacuum Failure Detected	1014, 2013	Vacuum failure detected. Check secure tubing connections. Check for pinched tubing. Ensure Tubing Cartridge is loaded properly. Press "Reset" to retry.
Probe Error - Hold Current Out of Range	1015	A Probe error occurred. Remove Probe to continue.
Probe Removed	1016	Probe Removed. Please wait.
Driver Error – Unrecoverable Error during Stroke	1017	Driver Error Occurred. Remove Probe or Disconnect Driver to continue.
Indexing Error	1018	Indexing Error occurred. Press "Reset" to continue.
Driver Error – Tool State Integrity Failure	1019	An error with the Driver occurred. Disconnect the Driver to continue. Connect Driver to retry. If the problem persists, contact Bard for service.
Vacuum Failure Detected	1020	Vacuum failure detected. Restart System. If problem persists, connect a different Driver.
Error Opening Cutter	1021	An error with the Probe occurred. Remove the Probe or disconnect the Driver to continue.
Cutter Not Centered	1022	An error with the Probe occurred. Remove the Probe or disconnect the Driver to continue.
Probe Unlatched	2011	The Probe is unlatched. Insert the Probe and press "Reset".
Vacuum Failure Detected	2014, 2017	A Vacuum failure occurred. Restart the system. If problem persists, contact Bard for service.
Driver Failure– Sample/Open Failed	2015	An error with the Driver occurred. Remove the Probe to continue.
Driver Error – Tool State Integrity Failure	2016	An error with the Driver occurred. Disconnect the Driver to continue.
UI Application Fault Detected	3000, 3001, 3004	An error occurred. Restart the system. If problem persists, contact Bard for service.
System Communication Fault Detected	3002	A communication error occurred. Restart the system. If problem persists, contact Bard for service.
Incompatible Driver Script Version	3003	The Driver script must be updated. Follow the onscreen instructions to initiate the update process, or contact Bard for service.
System detected UI Communications Fault	3005	A communication error occurred. Press Reset to continue. If problem persists, contact Bard for service.
Driver Script Update Recommended	3006	System detected a Driver script update is available. User may update via on-screen prompts or skip.

Specifications

Table 2: Specifications	
Environmental Conditions	Operate the device in ambient conditions (between 15°C to 25°C, 30% to 75% relative humidity, and 700 hPa to 1060 hPa atmospheric pressure). Transport and store in a cool dry place between -20°C to +60°C, 10% to 90% relative humidity, and 500 hPa to 1060 hPa atmospheric pressure.
Ingress of Water	Not protected against the ingress of water. IPX0.
Flammability	Equipment not suitable for use in the presence of flammable anesthetic.
Dimensions	16 in wide x 26 in deep x 46 in (retracted) to 57 in (extended) high (41 cm x 66 cm x 117 cm retracted, 145 cm extended)
Weight (Approximate)	105 lbs (48 kgs)
System Power Ratings	Maximum: 100-120 VAC / 127 VAC, 600VA / 750VA, 50/60 Hz, 3-wire, grounded with removable Power Cord 220-240 VAC, 600VA, 50/60 Hz, 3-wire, grounded with removable Power Cord
Electrical Conformity	This medical equipment has passed all required testing for electric shock, fire, and mechanical hazards in accordance with UL60601-1, IEC/EN 60601-1, CAN/CSA C22.2 No 601-1.
Classification	Class I, Type BF equipment

Table 3: Nominal Ratings - Standby	
Nominal Ratings	Standby
Voltage/Frequency	115 VAC ~ 50/60 Hz 230 VAC ~ 50/60 Hz
Current	4
VA (Calculated)	600
Watts	450
Power Factor	0.75

Table 4: Nominal Ratings - Operating	
Nominal Ratings	Operating
Voltage/Frequency	115 VAC ~ 50/60 Hz 230 VAC ~ 50/60 Hz
Current	4
VA (Calculated)	600
Watts	450
Power Factor	0.75

Electromagnetic Emissions and Immunity

Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided below. Portable and Mobile RF Communications Equipment can affect Medical Electrical Equipment.

Guidance and Manufacturer's Declaration – Emissions

All ME Equipment and ME Systems

Table 5: Guidance and Manufacturer's Declaration - Emissions		
The ENCOR ENSPIRE® Breast Biopsy System is intended for use in the electromagnetic environment specified below. The customer or user of the ENCOR ENSPIRE® Breast Biopsy System should ensure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF Emissions CISPR 11	Complies Group 1 Class A	The ENCOR ENSPIRE® Breast Biopsy System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Harmonics IEC 61000-3-2	Class A - Complies	The ENCOR ENSPIRE® Breast Biopsy System is suitable for use in all establishments, other than domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Flicker IEC 61000-3-3	Complies	

Guidance and Manufacturer's Declaration – Immunity

All ME Equipment and ME Systems

Table 6: Guidance and Manufacturer's Declaration – Immunity			
The ENCOR ENSPIRE® Breast Biopsy System is intended for use in the electromagnetic environment specified below. The customer or user of the ENCOR ENSPIRE® Breast Biopsy System should ensure that it is used in such an environment.			
Immunity Test	EN/IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
ESD EN/IEC 61000-4-2	± 2 kV, ± 4 kV, ± 6 kV, ± 8 kV Contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV Air	Complies	Floors should be wood, concrete or ceramic tile. If floors are synthetic, the relative humidity should be at least 30%.
EFT EN/IEC 61000-4-4	± 2 kV, 100 kHz PRF Mains ± 1 kV, 100 kHz PRF I/O's	Complies	Mains power quality should be that of a typical commercial or hospital environment.
Surge EN/IEC 61000-4-5	± 0.5 kV, ± 1 kV, ± 2 kV	Complies	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips/Dropout EN/IEC 61000-4-11	100% Dip for 0.5 Cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°	Complies	Mains power quality should be that of a typical commercial or hospital environment. If the user of the ENCOR ENSPIRE® Breast Biopsy System requires continued operation during power mains interruptions, it is recommended that the ENCOR ENSPIRE® Breast Biopsy System be powered from an uninterruptible power supply or battery.
	100% Dip for 1 Cycle	Complies	
	30% Dip for 25/30 Cycles	Complies	
	100% Dip for 250/300 Cycles	Complies	
Power Frequency 50/60 Hz Magnetic Field EN/IEC 61000-4-8	30 A/m	Complies	Power frequency magnetic fields should be that of a typical commercial or hospital environment.

Guidance and Manufacturer's Declaration – Immunity

All ME Equipment and ME Systems that are NOT Life-supporting

Table 7: Guidance and Manufacturer's Declaration – Immunity			
The ENCOR ENSPIRE® Breast Biopsy System is intended for use in the electromagnetic environment specified below. The customer or user of the ENCOR ENSPIRE® Breast Biopsy System should ensure that it is used in such an environment.			
Immunity Test	EN/IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Conducted RF EN/IEC 61000-4-6	3 Vrms 150 kHz – 80 MHz	Complies	Portable and mobile communications equipment should be separated from the ENCOR ENSPIRE® Breast Biopsy System by no less than the distances calculated/ listed below:
			$D = (3.5/V1)(\text{Sqrt } P)$
			$D = (3.5/E1)(\text{Sqrt } P)$ 80 MHz to 800 MHz
Radiated RF EN/IEC 61000-4-3	3 V/m 80 MHz – 2.7 GHz 80% AM at 1 kHz	Complies	$D = (7/E1)(\text{Sqrt } P)$ 800 MHz to 2.7 GHz
			where P is the max power in watts and D is the recommended separation distance in meters.
			Field strengths from fixed transmitters, as determined by an electromagnetic site survey, should be less than the compliance levels (V1 and E1). Interference may occur in the vicinity of equipment containing a transmitter.

Recommended Separation Distances between Portable and Mobile RF Communications Equipment and the ENCOR ENSPIRE® Breast Biopsy System Equipment and Systems that are NOT Life-supporting

Table 8: Recommended Separations Distances for the ENCOR ENSPIRE® Breast Biopsy System			
The ENCOR ENSPIRE® Breast Biopsy System is intended for use in the electromagnetic environment in which radiated disturbances are controlled. The customer or user of the ENCOR ENSPIRE® Breast Biopsy System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF Communications Equipment and the ENCOR ENSPIRE® Breast Biopsy System as recommended below, according to the maximum output power of the communications equipment.			
Max Output Power (Watts)	Separation (m) 150 kHz to 80 MHz $D = (3.5/V1)(\text{Sqrt } P)$	Separation (m) 80 MHz to 800 MHz $D = (3.5/E1)(\text{Sqrt } P)$	Separation (m) 800 MHz to 2.7 GHz $D = (7/E1)(\text{Sqrt } P)$
0.01	.1166	.1166	.2333
0.1	.3689	.3689	.7378
1	1.1666	1.1666	2.3333
10	3.6893	3.6893	7.3786
100	11.6666	11.6666	23.3333

Table 9: Cable Length	
Cable	Length (Minimum)
AC Power Cord	12 ft (US) / 2.5 m (International)

Equipotential Post

The purpose of the Industry Standard Equipotential Post on the rear of the ENCOR ENSPIRE® Breast Biopsy System is to provide a connection point to reduce the possibility of voltage potentials between touchable conductive parts of all the different Medical Equipment Systems within the Patient and Operator area.

To connect the ENCOR ENSPIRE® Breast Biopsy System to an Equipotential Network, snap the Equipotential Network Socket over the Equipotential Post on the rear of the ENCOR ENSPIRE® Breast Biopsy System base unit. To remove The ENCOR ENSPIRE® Breast Biopsy System from the Equipotential Network, pull the equipotential Network Socket off the Equipotential Post.

How Supplied/Stored

The ENCOR ENSPIRE® Breast Biopsy System is supplied non-sterile and should be cleaned prior to first use per the 'Cleaning and Maintenance' section.

- ENCOR® and ENCOR® MRI Drivers are sold separately. The Drivers are supplied non-sterile. Reference the Instructions for Use included with each Driver for Cleaning and Maintenance instructions.
- The Vacuum Canister, Vacuum Tubing Cassette, Vacuum and Rinse Tubing Cassette, and ENCOR® and ENCOR® MRI Probes are sold separately. The Vacuum Canister, Probes, and Tubing Cassettes are supplied sterile for single patient use only.

Warranty

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

Extended warranty service agreements are available. Consult a Bard representative for details regarding terms and conditions.

An issue or revision date and a revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular, Inc. to see if additional product information is available.

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Copyright © 2019 C.R. Bard, Inc. All rights reserved.
Copyright © 2019 C.R. Bard, Inc. Tous droits réservés.
Copyright © 2019 C.R. Bard, Inc. Alle Rechte vorbehalten.
Copyright © 2019 C.R. Bard, Inc. Tutti i diritti riservati.
Copyright © 2019 C.R. Bard, Inc. Todos los derechos reservados.
Copyright © 2019 C.R. Bard, Inc. Alle rechten voorbehouden.