TERUMOBCT

T-RAC III

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02/2020

MỤC LỤC

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
Windows Security

iexplore

The server 192.168.1.253 is asking for your user name and password. The server reports that it is from .

Warning: Your user name and password will be sent using basic authentication on a connection that isn't secure.

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The screenshot shows the B&B Electronics web interface for an Ethernet Router. The main menu includes Status, Configuration, Certificates, Network, and Maintenance. The Configuration menu is selected. On the left, a sidebar lists options like Express Setup, WLAN Settings, WLAN Security Settings, Network Settings, Serial Port Settings, Serial Port 2 Settings, Connection Settings, Ethernet Settings, Event Settings, Advanced Settings, Upload Configuration File, List Configuration Files, and Delete Configuration File. The Active Configuration and Clear Configuration buttons are also visible. The right panel displays the 'WLAN Parameters' configuration page. The 'Current Values' section includes fields for Radio Startup Mode (On), WLAN Connection Type (Infrastructure), SSID (linksys), WLAN TX Power (15 dBm), Maximum Wireless Data Rate (Auto), Specify Specific Data Rates (0x00000000), Use Fixed Data Rate (Disabled), WLAN Region (United States (US)), Preferred Radio Band (Auto), WLAN Maximum Retries (13), and WLAN RTS Threshold (0). At the bottom of this panel are 'Commit', 'Cancel', and 'Defaults' buttons, with 'Commit' being highlighted by a yellow oval. The status bar at the bottom of the browser window shows 'http://192.168.1.253/ Airborne0D6963 - Airborne ...'.

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The screenshot shows the B&B Electronics web interface after committing configuration changes. The main menu and sidebar are identical to the previous screenshot. The right panel now displays a large bold message: 'Configuration changes committed successfully'. Below this message are two buttons: 'Reload' and 'Restart', with 'Reload' being highlighted by a yellow oval. The status bar at the bottom of the browser window shows 'http://192.168.1.253/ Airborne0D6963 - Airborne ...'.

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



Low Frequency

OR

***Medium Frequency
(Optimum Settings)**

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
TERUMOBCT

T-RAC II

TERUMO RECORDING & AUTOMATIC BLOOD COLLECTOR

INSTRUCTIONS FOR USE

02/2020

CONTENTS

1	Introduction	4
2	Symbols and definitions	4
2.1	Symbols	4
2.2	Definitions.....	5
3	Important safety information.....	5
4	Specifications and accessories	7
4.1	Specifications.....	7
4.2	Accessories	8
5	Operating principles and functions.....	9
5.1	Battery	10
5.1.1	Emergency battery (only on T-RAC II Jr)	11
5.1.2	High-energy batteries.....	11
5.2	Clamp	13
5.3	On/off button	16
5.4	Tray.....	16
5.5	Touch screen.....	16
5.6	Barcode reader.....	16
5.7	Remote control.....	17
5.7.1	Multifunctional Pole	18
5.7.2	Desktop Remote Control	19
5.8	Tube Holder.....	19
5.9	TSM Battery Unit Holder	19
5.10	TOMEs connection	20
5.10.1	TOMEs Mobile Collection Tool (MCT)	22
5.11	USB flash drive	22
5.12	Transport Case	23
5.12.1	Installation.....	24
5.12.2	Battery charging	25
5.13	T-RAC II Quad charger	25
6	Installation and operation instructions	26
6.1	Installation – first time use.....	26
6.2	Operating instructions	27
6.2.1	Startup.....	27
6.2.2	Idle state	28
6.2.3	Donation cycle.....	29
6.2.4	Events	34
7	Settings.....	41
7.1	Device settings	45
7.2	Programs	46
7.2.1	Program settings	48
7.2.2	Program flow.....	55
7.2.3	Barcode rules	56
7.3	Communication settings	57
8	Alarms and error messages	58
8.1	Donation alarms	58
8.2	System alarms	60
8.3	Event warnings	63
9	Troubleshooting	65
10	Cleaning and maintenance	65
10.1	General cleaning	65
10.2	Preventive maintenance.....	66
10.2.1	Check of the weighing function	66
10.2.2	Check of the alarm function	66
10.2.3	General maintenance	66
11	High-energy battery shipping information	67
12	Disposal (end of life).....	68
13	T-RAC II Wifi setup	69
14	Symbol LS-2208 barcode reader set-up	74

1 Introduction

Intended use

T-RAC II is a blood collection weighing and mixing device with advanced safety and usability features. T-RAC II controls, by continuously measuring the weight, the collection of blood during donation. T-RAC II functions with different types and brands of blood bags. It mixes the collected blood to prevent clotting. The device is intended to be operated by nurses experienced in blood collection, however only under supervision of a doctor present during donation. The device can be transported for usage at different places.



T-RAC II shall be used in a quiet, clean and well-lit room, which is free from heavily vibrating devices, strongly electromagnetic radiating devices and noisy devices.

T-RAC II Ethernet and Wifi can be connected to the following external devices: a central computer system, a barcode reader, a remote control and a USB flash drive. This enables verification of barcodes, remote programming of T-RAC II parameters and collection of donation data.

Throughout this manual, information that only applies to T-RAC II Ethernet and Wifi will be printed in *italics*. This information doesn't apply to T-RAC II Jr.

2 Symbols and definitions

2.1 Symbols

REF	Catalogue number		This way up
SN	Serial number		Fragile, handle with care
LOT	Batch code		Keep dry
	Class II equipment		Temperature limit
	Humidity limitation.		<i>Non-ionizing radiation</i>
	Consult instructions for use		On/off button (stand-by)
	Manufacturer		<i>Ethernet connector (at the back)</i>
	Direct current		<i>USB connector for the barcode reader (at the back)</i>
EC REP	Authorized representative in the European Community		<i>USB connector for the USB flash drive</i>
	<i>Ethernet port for wifi</i>		<i>Connector for the Desktop Remote Control (at the back)</i>

This product complies with the requirements of the Medical Device Directive 93/42/EEC.

Separate collection for electrical and electronic equipment and for batteries (European symbol)

2.2 Definitions

WARNING: Informs the operator of a hazard or an unsafe practice that could result in personal injury, affect the operator's health or contaminate the environment.

NOTE: Highlights essential information.

TOMEs (Terumo Operational Medical Equipment software): dedicated system for data-exchange with several Terumo transfusion devices. More information can be found in chapter 5.10.

BBIS (Blood Bank Information System): software used in the blood bank for administration of donors (and donations).

TSM (T-SEAL MOBILE): Portable Battery sealer which is a fully automatic system for sealing PVC (polyvinyl chloride) or EVA (ethyl vinyl acetate) tubes, especially for tubes in blood pack systems.

3 Important safety information

If the equipment is not used in accordance with the instruction as described in this manual, the manufacturer will not guarantee the safety, reliability nor the performance of the device and shall not be liable.

These Instructions For Use must be followed whether for use or maintenance of the device. All users must be trained according to these instructions. The presence of a doctor approved by the transfusion center is obligatory during blood collection.

T-RAC II should be used in an indoor place which is quiet, clean and well-lit. The place should be free of heavily vibrating devices, strongly electromagnetic radiating devices and noisy devices. T-RAC II is not to be used in oxygen-rich environments (e.g. operating room).

Portable and mobile RF (radio frequency)communications equipment can affect medical electrical equipment.

When transporting or storing T-RAC II or its accessories, place them in the dedicated Transport Case or original shipping carton. Make sure that no liquid can fall on the case or shipping carton and that the storage conditions are met (chapter 4.1).

The essential performance of T-RAC II is: weighing, agitating and displaying information on the screen.

T-RAC II Wifi includes an RF transmitter which works on frequency bands [2400MHz - 2483,5MHz], [5150MHz - 5350MHz] and [5725MHz - 5825MHz]. It is compatible with IEEE 802.11 standard 'a', 'b', 'g' and 'n' and emits less power than 100mW.

WARNINGS

- Do not touch the connectors of the battery *or Multifunctional Pole* with bare hands in order to prevent damage to the device because of electrostatic discharge.
- Only the power supply which is delivered with T-RAC II should be used. The power supply is a part of this equipment.
- No modification of this equipment is allowed.
- T-RAC II is a measuring device. A precision strain gauge sensor measures the volume donated. The integral device should be handled with the greatest possible care. Whilst T-RAC II is capable of weighing objects, its use in commercial transactions is forbidden.
- The volume stated on the primary blood bag must not exceed 650ml.
- Do not use T-RAC II near devices which emit strong electromagnetical radiation. This may cause electromagnetical interference. Portable and mobile RF communications equipment can affect medical electrical equipment.
- When connecting the T-RAC II to an ME system, ensure safe operation by checking the specifications with the ME system manufacturer. The responsible organization is responsible for verifying that the system conforms to IEC 60601-1 and IEC 60601-1-2.
- The T-RAC II should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, observe the T-RAC II closely to verify normal operation in the configuration in which it will be used.
- According to IEC 60601-1-2 you should only use those accessories, transducers or cables specified in the manual, because the use of non-specified items may result in increased emissions or decreased immunity of the T-RAC II.
- When transporting the T-RAC II you should use the original shipping carton or the Transport Case especially made for this purpose.
- *If T-RAC II devices in the same area are running different programs, their alarm systems can react differently if the user has configured so. This can lead to a potential hazard.*
- *Make sure that the High-Energy battery is fully charged before usage (charging time of 3 hours).*

T-RAC II will lose donation data (including information from the alarms) of the ongoing donation if both the battery power and the mains power fall out during the donation process.

4 Specifications and accessories

4.1 Specifications

Product code	T-RAC II Junior: REF code: ME-TRAC2JR T-RAC II Ethernet: REF code: ME-TRAC2 T-RAC II WiFi: REF code: ME-TRAC2W	Order code: 2METRAC2JR Order code: 2METRAC2 Order code: 2METRAC2W
Size (W x H x D)	192x170x493mm	
Weight	T-RAC II Junior: 3.7kg T-RAC II Ethernet: 4.0kg T-RAC II WiFi: 4.1kg	
Voltage input	100 – 240 V AC	
Voltage output of power supply	12V DC	
Voltage input frequency of power supply	47–63 Hz	
Voltage input of T-RAC II	12V DC	
Power	57VA (100V inlet power) – 74.4VA (240V inlet power)	
Emergency battery	Type: NiMH Nominal voltage: 7.2V Nominal power: 10Wh	
High-energy battery	Type: Li-ion Nominal voltage: 7.5V Nominal power: 51Wh	
High-energy battery XLR	Type: Li-ion Nominal voltage: 7.3V Nominal power: 49.6Wh	
Range of balance	0 – 1000g (1g step)	
Weighing accuracy of the balance	± 3g for 0 – 500g; ±1% for 501 – 1000g	
Display range of flow rate	0 – 999 ml/min (1ml step)	
Display range of volume delivered	0 – 650 ml (1ml step)	
Display range of donation time	00:00 – 30:00 (mm:ss)	
Operating conditions	Ambient temperature: from 10°C to 40°C (50 up to 104°F) Relative humidity: 30% - 85% (non-condensing)	
Maximum altitude	2000m	
Storage conditions	Ambient temperature: from -15°C to 40°C (5 up to 104°F) Relative humidity: 10% - 90% (non-condensing)	
Pollution degree	2	
Classification (IEC 60601-1)	Class II and internally powered equipment, continuous operation, IPX0	
Noise level (continuous)	Agitating : 47dB Non-agitating : No noise Alarms : 55dB – 70dB	
Communication ports (only on T-RAC II Ethernet and WiFi)	1 x USB at the front for the USB flash drive 1x USB at the back for the barcode reader 1 x connector for Desktop Remote Control 2 x connectors for Multifunctional Pole 1 x RJ-45 for Ethernet connection 1 x RJ-45 for wifi connection (only on T-RAC II WiFi)	
Compliance to directives	Medical Device Directive 93/42/EEC	
EMC Classification	Group 1, Class B	
Manufacturer	Terumo BCT Inc., 10811 W. Collins Ave., Lakewood, CO 80215, USA	

4.2 Accessories

List of accessories		
T-RAC II Junior	<i>T-RAC II Ethernet and Wifi</i>	
D	D	AC power cord and power supply; Instructions for use (English); Instructions for use other languages: CD-ROM; Warranty Card
D*	X	Emergency battery, built-in in T-RAC II Junior (REF code: ME+05TRC225 – Order code: 8ME05TRC225)
O*	D	High-energy battery (REF code: ME+05TRC226 – Order code: 8ME05TRC226)
O*	D	High-energy battery XLR (REF code: ME+05TRC320 – Order code: 8ME05TRC320)
X	D	<i>USB Flash drive</i> (REF code: ME+05TRC227 – Order code: 8ME05TRC227)
X	D	<i>Barcode reader + cable</i> (REF code: ME+05TRC228 – Order code : 8ME05TRC228)
O	O	Transport Case (REF code: ME-TRAC203 – Order code : 2METRAC203)
X	O	<i>Multifunctional Pole including long antenna 25 cm with barcode reader holder and TSM Hand Unit holder</i> (REF code: ME-TRAC204 – Order code : 2METRAC204)
X	O	<i>Short antenna 12.5cm for Multifunctional Pole</i> (REF code: ME+05TRC217 – Order code : 8ME05TRC217)
X	O	<i>Tube Holder (for Multifunctional Pole: 2 parts / packaging)</i> (REF code: ME-TRAC205 – Order code: 2METRAC205)
X	O	<i>Desktop Remote Control</i> (REF code: ME-TRAC202 – Order code : 2METRAC202)
X	O	<i>TSM Battery Unit Holder (for attaching TSM Battery Unit to T-RAC II)</i> (REF code: ME-TRAC206 – Order code : 2METRAC206)
X	O	<i>TOMEs</i> (REF code: ME-TOMES – Order code : 2METOMES (only until TOMEs version 6.2 and for China) REF/order code: 8METOMESHW (Hardware only, no TOMEs preinstalled))
X	O	<i>TOMEs software license for T-RAC II</i> (REF/order code: 2METOMESTRA01 TOMEs software license for T-RAC for 1 (one) connection REF/order code: 2METOMESTRA05 TOMEs software license for T-RAC for 5 (five) connections REF/order code: 2METOMESTRA10 TOMEs software license for T-RAC for 10 (ten) connections)
X	O	<i>TOMEs Installer</i> (REF code: ME-TOMESINST – Order code : 2METOMESINST)
X	O	<i>T-RAC II Quad Charger</i> (REF/order code: 2METRAC207)

(D = Included by **default**, O = optional accessory, X = not available for this device)

[WARNING] Only the accessories listed above may be used with T-RAC II.


[NOTE] Both 2METOMES / 8METOMESHW and 2METOMESTRA01 / 2METOMESTRA05 / 2METOMESTRA10 are needed in order to transfer data between T-RAC II and TOMEs.

[NOTE] * Replacing the emergency battery on T-RAC II Junior needs to be done by a service engineer.

5 Operating principles and functions

T-RAC II (Terumo Recording and Automatic Collector) is designed to collect blood in different types and brands of blood bag systems. During donation, it shakes the blood bags in order to mix the blood with anti-coagulant. At the same time, the weight of the collected blood is measured. At the end of the donation, the clamp closes to stop the blood flow. T-RAC II shows detailed information about the donation on the screen and activates alarms in case of problems.

T-RAC II features a data-acquisition system. A barcode reader allows the operator to collect more information about the donation. This data can be transferred to a central computer system called TOMEs over a network or by using the USB flash drive. TOMEs can also be connected to a blood bank management system and makes it possible to set device or donation parameters on T-RAC II devices.

Desktop Remote Control (chapter 5.7.2)

Transport Case (chapter 5.12)

T-RAC II Quad Charger (chapter 5.13)

5.1 Battery

Three types of battery exist for T-RAC II:

- Emergency battery, as part of T-RAC II Junior.
- High-energy battery, delivered with T-RAC II Ethernet and Wifi. This battery can also be ordered separately for T-RAC II Junior. However this version will become obsolete in the future and will be replaced by its XLR version.
- High-energy battery XLR, delivered with T-RAC II Ethernet and Wifi. This battery can also be ordered separately for T-RAC II Junior.
- The batteries are charged automatically when T-RAC II is connected to the mains, regardless whether T-RAC II is switched on or off. Make sure that the High-Energy batteries are fully charged before usage (charging time of 3 hours).

WARNINGS

- Never open the battery or unscrew its screws. The battery must be protected from extreme temperatures, extreme humidity or fire. Replace the high-energy battery in case it's broken or contact Terumo BCT for replacing the emergency battery (T-RAC II Jr.).
- Batteries need special recycling. Do not throw them away with domestic waste or send a broken battery back to Terumo BCT (chapter 12).

NOTES

- There is also a small non-rechargeable lithium battery inside the device which allows to keep time and date. Normally, this battery shouldn't be replaced during the life time of T-RAC II. If time and date aren't kept correctly by the device, please contact your Terumo BCT dealer in order to replace this battery.

5.1.1 Emergency battery (only on T-RAC II Jr)

If the mains power is lost, T-RAC II Junior will detect this. The system will automatically switch to battery power. In case of an ongoing donation, the battery will allow the user to safely finish the donation. This battery isn't meant to replace the power supply, it must only be used in emergencies. No donations can be started when T-RAC II Junior is not connected to the mains.

The battery's state can be monitored on the touch screen:

- Emergency battery status OK. A donation can be ended in a safe way.

- Emergency battery empty or broken.

Battery technology: NiMH (Nickel-metal hydride)

Terumo BCT advises to have the battery replaced each 4 years or if the battery icon is red when working on mains power.


[WARNING] The emergency battery is part of T-RAC II Junior. Do not try to remove or replace the emergency battery. Contact Terumo BCT for this kind of servicing.

5.1.2 High-energy batteries


The battery provides T-RAC II with a fully independent power source, especially useful for mobile donations. When running on the battery, T-RAC II will work the same as if it was connected to the mains.

The battery is inserted at the front of T-RAC II. Switch off the device before (re)placing the battery, unless it is connected to the mains. If the device is connected to the mains, it does not need to be switched off for (re)placing the batteries.

To place the battery:


1. Push the battery forward into the battery compartment.


2. The battery will 'click' into place.

To remove the battery:

1. Press down the battery handle.

2. Pull-out the battery pack.

When the mains are disconnected, the battery can also be charged while T-RAC II is in its Transport Case. More information can be found in chapter 5.12.

The battery's state can be monitored on the touch screen:

- Battery level 0% - 9%, battery broken or not present. When running only on batteries, it will be impossible to start a donation.
- Battery level 10% - 24%
- Battery level 25% - 39%
- Battery level 40% – 54%
- Battery level 55% – 69%
- Battery level 70% – 84%
- Battery charging

Because of the technology used in this high-energy battery, no specific battery service is required. The battery can be charged and discharged regardless the battery level.

Battery technology: Li-ion (Lithium-ion)

Recharging typically takes 3 hours.

With a new battery, following amounts of donations can be performed (tested for donations taking 8 minutes of full-speed agitation, without events). For the High-energy battery:

- *±100 donations on T-RAC II Junior.*
- *±74 donations on T-RAC II Ethernet.*
- *±54 donations on T-RAC II Wifi when the setting 'Enable LAN' is set to 'Yes'.*

For the High-energy battery XLR:

- *±99 donations on T-RAC II Junior.*
- *±56 donations on T-RAC II Ethernet.*
- *±39 donations on T-RAC II Wifi when the setting 'Enable LAN' is set to 'Yes'.*

It depends on the settings of the device (e.g. Bcklight lev. and Mix freq.) and on the way the device is used (e.g. switching off the device after donation or not) how much donations can be done.

Because batteries age, this value can be lower after some time. If the battery power isn't sufficient anymore, replace it with a new one.

[WARNINGS]

- Please refer to chapter 11 for the shipping information of the high-energy batteries!
- If no battery is present, an interruption of the mains current will cause T-RAC II to shut down immediately. If that happens during donation, the clamp will not close. *T-RAC II will lose donation data (including information from the alarms) of the ongoing donation.*
- To prevent short-circuits, never connect the battery's connector to anything else than the battery compartment in T-RAC II. Never connect the battery directly to the Multifunctional Pole. This will cause permanent damage to both.
- Only insert T-RAC II batteries in the battery compartment.
- When the high-energy battery is not likely to be used for some time, take it out to prevent leakage.
- When the battery is discharged, it needs to be recharged immediately. Never store batteries which are completely discharged.

[NOTE]

- When running only on the battery and the battery is nearly empty, T-RAC II can refuse to start a donation. An error message will be given.

5.2 Clamp

The clamp controls the blood flow. The donor tube between the venipuncture needle and the primary blood bag should be inserted in the clamp. The clamp has three positions which can be monitored on the touch screen:


- Open. In this state, the tube can be inserted in and removed out of the clamp. This state is used when the device is idle.

- Locked. The tube cannot be removed from the clamp, but the tube can still move axial. Blood can flow through the tube. This state is used while donating in order to prevent the tube from leaving the clamp.


- Closed. The tube is squeezed by the clamp. The tube cannot move and blood cannot flow. This state is used before and after donation, and during some alarms.


After pushing the clamp head backwards after tube insertion (see step 2 in chapter 5.2) and then releasing it, the clamp will change from the open to the locked state. This facilitates handling because it allows the user to feel immediately whether the tube is correctly locked in place.


To remove the clamp in case of loss of power:

Make sure the device is switched-off.
Rotate the clamp head 20° counter clockwise.


Pull-up the clamp head.
The tube can now be removed out of the clamp head.


Make sure that the inner black piece (in the center of this top view) does not rotate.
Do not switch-on T-RAC II while the clamp head is removed.

To put the clamp back in position:


Make sure that the high part of the inner black piece is pointing towards the front, on a 20° angle, as shown in this top view.

Under an angle of 20°, put the clamp head back in position.

Without forcing, the clamp should slide completely down, as shown in the picture.

If this doesn't work the first time, make sure that both the inner black piece and the clamp head are under an angle of 20°. Then, try inserting again.

Without forcing, rotate the clamp head 20° clockwise.

If the clamp head points exactly towards the front of the device, the position is correct. The device can be switched-on again.

WARNINGS

- Always switch off the device before performing the clamp emergency exit.
- The clamp must be put back in place before supplying power to the device.
- At the moment of locking the clamp:
 - Make sure that the tube is inserted completely to the back of the clamp.
 - Make sure that the tube is not bent. It needs to be straight.
- After locking the clamp, always make sure that the tube is correctly locked inside the clamp. It must be possible to move the tube from left to right and it must be impossible to take the tube out of the clamp. If the tube is stuck in any way:
 1. Open the clamp by pressing 'Open clamp' and 'Yes'.
 2. Insert the tube correctly.
 3. Lock the clamp again.
 4. Recheck that the tube is correctly locked.
- Do not insert sharp objects in the clamp.

5.3 On/off button

Push the button once in order to switch on the device. Push the button and hold for 3 seconds in order to switch it off. The button features a bi-color LED:

- Off: T-RAC II is switched off and is not connected to the mains.
- Orange: T-RAC II is switched off, but is connected to the mains. The battery, if present, is charging or is fully charged.
- Green: T-RAC II is switched on.

5.4 Tray

The tray has been tested to fit common blood bag systems. The primary bag of the system (in which the blood is collected) should lay on top of the other bags and the bags should be placed in the middle of the tray. The tray features an elastic band as filter holder. Place the filter behind this elastic band. The tray is connected by a magnet to the main device. To remove it, push one of its sides down so the tray will come loose.

WARNINGS

- Make sure the tray is in place when the device is agitating. Never touch the area of the moving magnet. This may cause injury.
- Do not touch the tray at startup, at the screen 'initializing donation' or during donation. Otherwise, wrong donation data can be displayed and recorded.
- The volume stated on the primary blood bag must not exceed 650ml.
- Do not put heavy objects (more than 2kg) in the tray as this can damage the balance.
- The blood bag system should be placed in the center of the tray.

5.5 Touch screen

T-RAC II is equipped with a touch screen. The touch screen visualizes all information about the donation process and guides the user through the complete procedure with pictograms and illustrated instructions. It allows the user to set device and donation parameters and to give other inputs about a donation.

In case the touch screen doesn't respond correctly anymore, it needs to be calibrated. See chapter 9 for more information.

5.6 Barcode reader

The barcode reader can be connected to the USB slot at the back of T-RAC II, which is marked as follows.


In order to set-up the barcode reader to work with T-RAC II, scan the codes mentioned in chapter 14

5.7 Remote control

The remote control is designed to decrease the number of times the user needs to bend down. On top of the remote control, an antenna gives visual information about the state of the device.

The remote control exists in two versions which can be ordered separately: as part of the Multifunctional Pole or as a Desktop Remote Control.

- Four buttons allow the user to operate the touch screen. The pictograms on the buttons can also be found on the touch screen. Their function is the same. Buttons whose pictograms aren't visible on the touch screen, have no function.

In this case, the right button has no function since its pictogram is not shown on the touch screen.


- *The antenna is colored as follows:*

Red flashing slow (0.7Hz)	Alarm indicating the end of donation (medium priority alarm).
Red flashing fast (1.5Hz)	Alarm/error during a donation (medium priority alarm).
White	Donation process started (including processing events).
White flashing	Alarm / error not related to a donation (low priority alarm).
Yellow	Donation in progress, collected volume 0-25%.
Green	Donation in progress, collected volume 25-50%.
Blue	Donation in progress, collected volume 50-75%.
Purple	Donation in progress, collected volume 75-100%.


The color during donation changes gradually as follows. If the antenna flashes during the donation stage in one of these colors, a low flow rate pre-alarm is triggered.


5.7.1 Multifunctional Pole

The Multifunctional Pole can be put on two locations in the device. It can be removed or put in place even when the device is active.

On the pole, the antenna needs to be inserted manually in the dedicated hole. Push the antenna down until it blocks.

The pole features a barcode reader holder which can be rotated according to the user's preference. The pole also features a holder for the Hand Unit of T-SEAL Mobile. Both brackets of the holder can be rotated separately and according to the user's preference (2METRAC204).

Length of the pole: 420mm (excl. antenna); 630mm (incl. antenna).

[WARNING] Never connect the Multifunctional Pole to T-RAC II at the same time that a Desktop Remote Control is connected. Never connect more than one Multifunctional Pole to one T-RAC II device. This can permanently damage the device.

[WARNING] Never connect the Multifunctional Pole directly to the battery. This will permanently damage both.

5.7.2 Desktop Remote Control


The Desktop Remote Control can be used on a table near the device. Connect it to the dedicated slot at the back of the device, which is marked as follows.

The antenna on this remote control cannot be removed.



5.8 Tube Holder

A Tube Holder can be connected to the back bracket of the barcode reader holder on the Multifunctional Pole. This holder can be used for inserting the sample blood tubes before, during and after donation.

5.9 TSM Battery Unit Holder

A TSM Battery Unit Holder can be connected to the back end of the T-RAC II for inserting the Battery Unit of the TSM.


5.10 TOMEs connection

A TOMEs network connects one or more Terumo BCT transfusion devices with TOMEs, the dedicated central computing system. TOMEs can also be connected to a customer's network. TOMEs can be accessed via a web browser on the computers which are connected to the customer's network. Via the customer's network, TOMEs can also be connected to a Blood Bank Information System (BBIS) for bi-directional communication.

TOMEs adds following functionality to T-RAC II:

- Collecting and viewing donation data (this is done automatically when a network connection is present).
- Configuring multiple T-RAC II devices at once.
- Monitoring of the T-RAC II devices.
- Sending donation-specific settings to T-RAC II, based on information from a BBIS (bi-directional communication).

To connect T-RAC II to the network:

1. Make sure that the correct communication settings are filled-in (see chapter 7.3), according to your network.
 2. Make sure that the Ethernet functionality of T-RAC II is enabled.
 3. For wired Ethernet connection: connect the Ethernet cable (CAT5 UTP or better) to the port indicated as follows. **Wired Ethernet connection is possible both on T-RAC II Ethernet and T-RAC II Wifi.**
-
4. For wireless Ethernet connection (only on T-RAC II Wifi), follow the steps as described in chapter 13 .

If the network connection is enabled in the T-RAC II settings (chapter 7.2.1), the status of the network connection can be monitored on the touch screen:

-  Ethernet/Wifi disabled by the program setting 'Enable LAN' = No.
-  T-RAC II Ethernet/Wifi connected to TOMEs
-  T-RAC II Ethernet/Wifi disconnected

Pressing the pictogram will show a screen in which the Ethernet functionality can be quickly enabled / disabled.

[NOTE]

- When the Ethernet pictogram is green and blinking, data is being sent.

When the pictogram is red and blinking, the internal memory of T-RAC II is full. T-RAC II won't be able to save any more information from donations until the network connection is re-established or until the USB flash drive is inserted.

T-RAC II – TOMEs Compatibility Table:

Compatible?	T-RAC II firmware version					
	1.55	3.25	5.35	6.25	6.55	7.05
TOMEs version 3.8	Yes	Yes	No	No	No	No
TOMEs version 4.4	Yes	Yes	Yes	No	No	No
TOMEs version 5.0	Yes	Yes	Yes	No	No	No
TOMEs version 5.1	Yes	Yes	Yes	No	No	No
TOMEs version 5.2	Yes	Yes	Yes	Yes	Yes	Yes
TOMEs version 6.0	No	No	Yes	Yes	Yes	Yes
TOMEs version 6.1	No	No	Yes	Yes	Yes	Yes
TOMEs version 6.2	No	No	Yes	Yes	Yes	Yes
TOMEs version 6.3	No	No	Yes	Yes	Yes	Yes
TOMEs version 7.0	No	No	Yes	Yes	Yes	Yes
TOMEs version 7.1	No	No	Yes	Yes	Yes	Yes
TOMEs version 7.2	No	No	Yes	Yes	Yes	Yes

[NOTE] The network functions of T-RAC II can only be used with TOMEs or TOMEs Mobile Collection Tool. Following optional accessories are required:

- TOMEs hardware (2METOMES/8METOMESHW) or TOMEs Installer (2METOMEsINST / ME-TOMESINST)
- TOMEs software license for T-RAC II (2METOMESTRA01 / 2METOMESTRA05 / 2METOMESTRA10)

[NOTE] Please refer to the help-file of the T-RAC II TOMEs software module for more information.

[NOTE] Check that the network connection is stable in order to make sure that there is no delay in sending donation data to TOMEs or TOMEs Mobile Collection Tool.

5.10.1 TOMEs Mobile Collection Tool (MCT)

TOMEs Mobile Collection Tool (MCT) is a program which can be downloaded from TOMEs (starting from TOMEs version 6.1). It can be run from a computer without installation. The purpose of TOMEs Mobile Collection Tool is to use it on mobile donation sites. A small mobile network is set up with a pc running the TOMEs Mobile Collection Tool. T-RAC II devices can connect to the PC just as they can connect to TOMEs. The tool then saves all data from all connected T-RAC II devices on the computer running the tool. Alternatively, data from T-RAC II USB flash drives can be imported. When the computer reconnects with the TOMEs network, the tool automatically synchronizes all saved data to TOMEs.

[NOTE] Please refer to the TOMEs software module IFU for more information.

5.11 USB flash drive

The USB flash drive offers a means of data transfer to and from TOMEs or to TOMEs Mobile Collection Tool without the need for a TOMEs/MCT network. The drive can be inserted in the USB slot at the front of T-RAC II when T-RAC II is switched off.

Two types of data can be transferred with the drive:

- Donation data from T-RAC II to TOMEs or TOMEs Mobile Collection Tool. This is explained below.
- Programs from T-RAC II to other T-RAC II devices, to TOMEs and vice versa. See chapter 7.2.

Copying donation data from T-RAC II

T-RAC II features an internal memory to temporarily save donation data (information from the blood collection). When a USB drive is connected, all donation data present in the T-RAC II device are automatically copied to the USB drive if the device is not connected to the network. Also, after a new donation is performed on T-RAC II, its donation data will be saved to the USB flash drive. The USB flash drive may only be inserted or removed from T-RAC II when T-RAC II is switched off. Data which is copied to the drive, is automatically erased from the internal memory of T-RAC II. Terumo BCT suggests two ways of using the USB drive:

- **One USB flash drive for all T-RAC II devices**

When donations on all devices are finished, the USB drive is inserted one by one on each T-RAC II. All donation data will be on one stick, making it easier to copy the data to TOMEs or TOMEs Mobile Collection Tool.

- **One USB flash drive per T-RAC II**

In each T-RAC II, a USB drive is inserted. In this way, all drives can be collected at the end of all donations. This makes it easier to copy data from the T-RAC II devices to the USB drive.

The status of the USB flash drive can be monitored on the touch screen:

-  USB-stick is connected and all data is sent to the USB drive. If this icon blinks, this means that data is being sent to the drive.
-  USB-stick is disconnected or there is a problem. If this icon blinks, the internal T-RAC II memory is full.

[WARNING] The USB flash drive may only be inserted or removed from T-RAC II when T-RAC II is switched off. Otherwise, donation data may get lost and the USB drive may get damaged.

[WARNING] Use only the USB flash drive which is delivered with the device. The system is not validated with other drives. Don't put other files on the flash drive or modify the flash drive. The device may become unusable and data may get lost.

[WARNING] Never switch off T-RAC II while data is being transferred. When data is being transferred, this is visualized on the touch screen by a green blinking USB-symbol.

[NOTE] Terumo BCT recommends not to keep the USB flash drives inserted in the T-RAC II devices when the devices are connected to the network. If the network fails, the data would be copied automatically to the USB drive and not to the network, making it harder to find back.

[NOTE] The donation data files are stored in the folder ~/T-RACII/YYYYMMDD. Where YYYYMMDD is the date of the blood collection.

[NOTE] The data-transfer options can only be used with TOMEs or TOMEs Mobile Collection Tool. Following optional accessories are needed:


- TOMEs hardware (2METOMES / 8METOMESHW) or TOMEs Installer (2METOMEsINST)
- TOMEs software license for T-RAC II (2METOMESTRA01 / 2METOMESTRA05 / 2METOMESTRA10)

[NOTE] Please refer to the help-file of the T-RAC II TOMEs software module for more information.

5.12 Transport Case

As an optional accessory, a T-RAC II Transport Case can be ordered. All standard and optional accessories of each version of the T-RAC II system can be placed inside the case. Dimensions of the case: 220x360x580mm (W x H x D).

In the lid, there is space for following components:

The Hand Unit of T-SEAL Mobile can be put under the T-RAC II or in the tray of T-RAC II.

The Battery Unit of T-SEAL Mobile can stay connected at the rear of T-RAC II also when the case is closed.


For easier transportation, the cases can be stacked on top of each other.

For easier cable-management, 3 additional Velcro bands are supplied with the case.


[WARNING] Do not stack more than 3 cases on top of each other.

5.12.1 Installation


The case is designed to make installation for mobile collections as easy and as fast as possible. Proceed as follows to install a device:


1. Put the case next to the place where T-RAC II should be positioned during donation.

2. Open the case and put the lid upside down where T-RAC II should be placed.

3. Take out the accessories (barcode reader, Multifunctional Pole / Desktop Remote Control, Hand Unit of T-SEAL Mobile) and put them at the place where they should be while donating. The barcode reader, the Desktop Remote Control and the Battery Unit of T-SEAL Mobile can stay connected to T-RAC II, even when the case is closed.


4. Put the bottom side of the case on top of the lid (which is still upside down). The lid will serve as a table for the device.

5.12.2 Battery charging

The case features a special function which makes it possible to charge the battery while T-RAC II is inside the case.

In order to do this, proceed as follows:

Connect the internal power cable of the case to the power connector of T-RAC II.

 A line drawing of the T-RAC II device inside the case, showing the power cable from the power supply being connected to the power connector of the case. The connection point is highlighted with an orange line.

Connect the power cable coming from the power supply (which is outside the case) to the power connector of the case.

The case can now be closed and the power supply can be connected to the mains.

5.13 T-RAC II Quad charger

T-RAC II Quad Charger is intended to charge up to 4 T-RAC II high-energy batteries simultaneously: the operator takes the battery out of a T-RAC II device and insert it into the T-RAC II Quad Charger.

It is a standalone accessory which is not connected with T-RAC II. Although it is a medical device, it cannot be placed or used in the patient environment.

The device can be transported for usage at different places.


More information can be found in the separate T-RAC II Quad Charger IFU (LCT-7017).

6 Installation and operation instructions


6.1 Installation – first time use

Unpack the device and place it at the desired position in the donation area. T-RAC II is to be placed beside a blood donation bed in which the donor will lie. The front of the device is to be placed towards the operator. *At first time use, remove the blue protection foil from the battery and insert the battery in the device (chapter 5.1).* Connect the device to the power supply and the power supply to the mains. Press the on/off button to switch on the device. Keep the device connected to the mains until the battery pictogram shows that the battery is fully charged.


When the battery is charged, the device is ready for use. It can be configured to your needs (chapter 7) and used for blood donation (chapter 6.2.3). For more information about the parts of the device and its additional accessories, see chapter 5.

Connectors at the back:

USB connector for the barcode reader (chapter 5.6)


RJ-45 connector for Ethernet communication (chapter 5.10)

RJ-45 connector for wifi communication (chapter 13, only on T-RAC II WiFi)

Connector for the Desktop Remote Control (chapter 5.7.2)

Power connector, to be connected to the power supply.

[WARNING] Do not place the device in a manner so that it is difficult to reach the power plug. The power plug is to be used for disconnecting the device from the mains.

[WARNING] Never connect the Multifunctional Pole to T-RAC II at the same time that a Desktop Remote Control is connected. Never connect more than one Multifunctional Pole to one T-RAC II device. This can permanently damage the device.

[NOTE] When using the emergency battery on T-RAC II Junior, the device should always be connected to the mains. Its battery is only meant to finish an ongoing donation. See chapter 5.1.1.

6.2 Operating instructions

6.2.1 Startup


Switch on T-RAC II by pressing the on/off button. The device will perform its self check.

When not using the T-RAC II for a longer time (ie overnight), switch off the device by pressing and holding the on/off button for 3 seconds.

6.2.2 Idle state

Once T-RAC II is switched on, the touch screen will act as the main interface with the device. The 'idle screen' is shown. Starting from this screen, a donation cycle can be started (see 6.2.3) or the settings can be reached (chapter 7).

Following fields are always present on the touch screen, regardless its state:

Main window

This area informs about the current stage of the donation cycle. By using pictures, the window indicates which action needs to be performed. It can also ask the user for input (e.g. to choose or enter the donation volume).

Current date and time

Device status icons

Maximum 5 icons inform the user about the status of following parts. From left to right:

- Network connection
- Battery
- Mains connection
- USB flash drive
- Clamp

Each button can be pressed in order to obtain more information about a part's status or to change the status. The meaning of the icons is explained throughout chapter 5 .

Navigation buttons

Use the right button to go to the next stage of the donation cycle and the left one to go to the previous stage. On the idle screen, the left button leads to the statistics of the previous donation (if available).

When the device is in the idle state, the green button shows the current weight in the tray. Pressing the button will tare the scale. For the donation process, it's never necessary to press this button.

If a button is dark like the 'statistics' button, this means it is disabled. This counts for all buttons.

6.2.3 Donation cycle

This chapter will explain a complete donation cycle. It is assumed that all settings are set to their default values.

The most simple donation process is described, as it will also run on T-RAC II Junior. *On T-RAC II Ethernet / Wifi, events give a lot more flexibility (e.g. scanning barcodes) and are explained in chapter 6.2.4 and 7.2.2.*

During the cycle, various (error) messages or alarms may be shown. These are summarized in chapter 8 .

1. Start of donation cycle


The donation cycle can be started from the idle screen by putting the blood bag system in the tray. Make sure that the primary blood bag (the one in which blood will be collected) lies on top of the other bags.

T-RAC II will detect the bag and proceeds automatically to the next step. Otherwise, 'next' can be pressed.

[WARNING] When placing the bag system, make sure that the system is inside the tray and does not touch any other part of the T-RAC II.

2. Tube insertion in the clamp

Insert the donor tube into the clamp. Make sure that the tube is pushed well against the back of the clamp. Push the clamp backwards and then release it or press 'LOCK' on the touch screen. The clamp will go to the 'locked' state. For more information about the clamp, see chapter 5.2.

[WARNING] At the moment of locking the clamp:

- Make sure that the tube is inserted completely to the back of the clamp.
- Make sure that the tube is not bent. It needs to be straight.


[WARNING] After locking the clamp, always make sure that the tube is correctly locked inside the clamp. It must be possible to move the tube from left to right and it must be impossible to take the tube out of the clamp. If the tube is stuck in any way:

1. Open the clamp by pressing 'Open clamp' and 'Yes'.
2. Insert the tube correctly.
3. Lock the clamp again.
4. Recheck that the tube is correctly locked.

3. Events – volume chooser

In this phase information about the donation can be entered via 'events' displayed on the screen. By default, only a volume chooser is depicted here. *This screen can also show other 'events' besides the volume chooser (see chapter 6.2.4).*

Depending on the program setting 'choose vol.' (see chapter 7.2.1), the volume chooser looks different:


Volume chooser

Depending on the settings, this field looks different:

- **Presets**

The volumes are shown as presets (like shown in the picture). Simply press a preset to select the donation volume. Up to 4 presets can be configured.

- **Keyboard**

Press the digits to set the volume.

Press this button to confirm the volume.

Press this button to remove unwanted input.

- **Barcode only**

No manual input possible, the volume is determined by the scanned barcode.


[NOTE] In the settings it is also possible to set a fixed volume. In that case, the volume chooser won't be shown.

[NOTE] It is impossible to select a higher volume than the maximum donation volume configured in the settings.

[NOTE] When a network connection is present, it is also possible to receive the donation volume from the BBIS. Refer to the help-file of the T-RAC II TOMEs software module for setting this up.

4. Venipuncture

After giving input on the events, the clamp will close. The device resets the weight measured by the balance, performs a self-check and asks to perform the venipuncture. Perform the venipuncture and press 'Start' in order to start the donation.

[WARNING] Pay attention not to hit any cables or the donor tube as long as the needle is connected to the donor.

5. Donation in progress


During donation, the clamp is in the 'locked' state so blood can flow through the tube but the tube cannot fall out of the clamp. The screen and the antenna inform about the progress of the donation. During donation, a number of alarms can warn the user. An explanation of these alarms can be found in chapter 8 .

[NOTE] At a certain moment during donation, the clamp will automatically close for 15 seconds. This is done in order to check whether the tube is correctly inserted into the clamp. If weight addition is measured while the clamp is closed, the 'No tube in clamp' alarm will be shown as explained in chapter 8 . During these 15 seconds, the low flow (pre-)alarm will not be activated.


6. End of donation

After reaching the collection volume or the maximum donation time, the clamp will close and an alarm will sound. Proceed by pushing 'Next'.


7. Needle removal

Remove the needle from the donor's arm. Press 'open clamp' to proceed.


8. Removal of the blood bag system from the tray

The clamp is opened so the tube can be removed from the clamp. The blood bag system can be removed from the tray. T-RAC II will detect this and proceeds automatically to the next step. Otherwise, 'Statistics' can be pressed.

9. Statistics

Statistics about the donation are shown. The user can continue to the idle state by pressing 'Finish'. Alternatively, the screen proceeds automatically after 8 seconds.

If barcodes were programmed, this field indicates how many. It is red when not all codes have been scanned.

The incident icons are only shown when an 'incident event' is programmed. See chapter 6.2.4.

6.2.4 Events


The volume chooser is only one of 8 types of events which can be used. Other events, for instance, allow to read barcodes or to communicate with a BBIS. The input of each event is recorded in the donation data.

The events themselves and their location in the donation cycle cannot be configured on T-RAC II, they are configured on TOMEs. The events can be displayed at following locations:

- *At startup: before the 'idle screen'. Results of these events are recorded into the donation data from all donations which take place, until the device is restarted or until another program is selected.*
- *Before donation: during step 3 of the donation cycle.*
- *During donation: during step 5.*
- *After donation (but before needle removal): between step 6 and 7.*
- *After needle removal: between step 7 and 8.*

[NOTE] Both 2METOMES / 8METOMESHW or 2METOMEsINST and (2METOMESTRA01 / 2METOMESTRA05 / 2METOMESTRA10 are needed in order to use events on T-RAC II. Information about the configuration of events can be found in the help-file of the T-RAC II TOMEs software module.

Events are always shown on the touch screen as follows:


This arrow allows to go to the previous event and change it. When doing this, only the latest result will be saved.

List of consecutive events. The current event is shown on a white background. The results of previous events are shown.

Window with the contents of the current event. What is shown depends on the type of event and the event's configuration.

This arrow allows to skip an event and move on to the next one. This button is only present at events which can be skipped (configured on TOMEs).

Depending on the event type and its settings, the user can manually give input on an event by using an on-screen keyboard. It is possible to enter the characters of a barcode or to type a small message. The number of characters is limited to 25. A stylus may facilitate the typing. To reach the keyboard, press **ABC**.

Following types of events exist. They are explained throughout this chapter:


1. Volume chooser
2. Standard barcode
3. Multi-barcode
4. Number-barcode
5. Information
6. Incident
7. Program #
8. Weight check

1. Volume chooser

The volume chooser allows to set the donation volume. For more information on the volume chooser, see chapter 6.2.3.


2. Standard barcode

The user is requested to scan a barcode. The name of a standard barcode (in this case 'Scan nurse') as well as its picture are configured on TOMEs.


3. Multi-barcode

With this event, multiple barcodes of the same kind are read. Depending on the configuration of this event, the number of barcodes is fixed or taken from a 'number-barcode' event. This barcode is especially made for reading multiple bags or tubes.


4. Number-barcode

This event defines the number that is used in a 'multi-barcode' event. The number (1 – 6) can be entered with user-defined barcodes or by using the on-screen keyboard. When the event is skipped, the number is '0'. In that case, the corresponding multi-barcode event will be skipped automatically.

5. Information

Depending on the settings of this event, user-defined information or information from the BBIS is shown here. The user gives feedback on this information by pressing the 'agree' or the 'disagree' button. When 'disagree' has been pressed, it depends on the event settings whether the donation can be continued.


[NOTE] An active network connection is required to receive information from the BBIS.

6. Incident

Depending on this event's settings, the screen which is shown will look different. When touch screen input is chosen, following screen is shown:

Standard incidents

: donor unwell

: problem with disposable

: donor fainting


: problem with venipuncture

: problem with T-RAC

Icon to record that everything is OK.

This button switches between standard incidents and user-defined incidents.

User-defined incidents

II

: donor feels pain in his arm

: coagulation

: contamination

: leakage

Selected incidents are shown on a white background.

Button to reach the manual input via keyboard.

T-RAC II provides ten standard icons and additionally nine user-defined incidents in the form of numbers. One or more incidents can be selected by simply pushing the icons and/or the numbers. It is also possible to give manual input via an on-screen keyboard (max 25 characters). When the correct incidents are selected, the user needs

to confirm them by pressing .

Incidents by barcode

When barcode input is chosen instead of keyboard input, following screen is shown:

Multiple user-defined barcodes can be scanned and they will be shown on the touch screen one by one. They need to be confirmed by pressing . When an incorrect barcode was read, the scanned barcodes can be erased by pressing .

7. Program

As can be read in chapter 7 , T-RAC II can run up to ten programs. T-RAC II will automatically run program 1 when no other program is selected using this event. Selecting another program is done by pressing a program number on the screen. Then, the corresponding program name will be displayed. To load it, confirm by pressing .

When loading another program is not needed, simply skip the event by pressing .

The program can also be changed by scanning a program barcode. These barcodes are defined by the user, see setting 'prog. barcode' in chapter 7.2.1.

[NOTE] It is best to put this event as early as possible in the program flow. By switching to another program, the program flow will change as well.

[NOTE] When this event is shown at startup of the device, the chosen program will be the default program for all donations until the device is restarted.

[NOTE] Results from events preceding a 'Program #', won't be saved into the donation statistics if another program was selected.

8. Weight check

Some donation center's procedures request a confirmation of the balance's calibration. This event doesn't calibrate the scale but offers a means of validation. When the event is shown, the user needs to place a 200g weight in the tray. If the 200g ($\pm 3g$) is correctly weighed during 2 seconds, the event will proceed automatically to the next event. Otherwise, the scale needs calibration. Contact Terumo BCT for this service. It is not possible to skip this event unless it occurs after donation.

[NOTE] Only qualified Terumo personnel are authorized to do the calibration of the scale. Normally, no calibration will be needed during the life time of T-RAC II.

7 Settings


In order to offer maximum flexibility, the configuration of T-RAC II is divided in a number of setting types. The device settings and the communication settings are inherent to each device and need to be configured separately on each T-RAC II device. The programs, which contain most settings, are interchangeable with other T-RAC II devices and with TOMEs. This can be done via a TOMEs network (chapter 5.10) or via a USB flash drive (chapter 7.2). Some parts of the programs can only be configured on TOMEs.

On T-RAC II Jr, there is only 1 program available which cannot be exchanged with other devices.


Setting type	Description	Configurable on			
		T-RAC II Junior	T-RAC II Ethernet /Wifi	TOMEs + transmission via network	TOMEs + transmission via USB flash drive
Device settings	Device-dependent settings and information (e.g. serial number, language).	Y	Y	Y	N
Communication settings	Needed when setting up a network with TOMEs. More information can be found in the help-file of the T-RAC II TOMEs software module.	NA	Y	N	N
Programs	Each program consists of following 4 parts:	Y, 1 program available	Partly, 10 programs available	Y	Y
Programs – program settings	These are settings which are part of programs.	Y	Y	Y	Y
Programs – program flow	This part of the program configures the place of the 'events' in the donation cycle.	NA	N	Y	Y
Programs – barcode rules	This part of the program defines how barcodes should be checked with each other.	NA	Y	Y	Y
Event settings	These are the definitions of the events. This is explained in the help-file of the T-RAC II TOMEs software module.	NA	N	Y	Y

(Y = yes, N = no, NA = not available)


In order to reach the settings which are available on T-RAC II, press 'Settings' on the idle screen. A password will be asked. The default password is '9876'.


The settings are presented as follows:

To change a setting, push the field that needs to be altered. The field will become selected and two extra buttons will be shown:

In case the selected field is a number or a text, a keyboard will be shown. For instance, when altering the 'Time' setting, following screen is shown:

Instead of using the alphanumeric keyboard, it is also possible to enter the setting by reading a barcode.

[NOTE] A stylus is handy for easier navigation in the settings screens.

7.1 Device settings

These are the device-dependent settings, meaning they cannot be copied to other devices.
When using a TOMEs network, it is possible to alter them via TOMEs.

Following device settings are available:

Serial number

Information

Serial number of the device. This setting cannot be changed by the user.

Firmware v

Information

Version of the firmware. This setting cannot be changed by the user.

Service (don)

Information

Number of donations which have been performed since last servicing. This setting cannot be changed by the user.

Password

0000 – 9999

Default: 9876

Allows to change the password for entering the settings.

Service (mon)

Information

Number of months since last servicing. This setting cannot be changed by the user.

Language

DAN / DEU / ENG / ESP / FRA / ITA / NED / NOR / POR / RUS / SUO / SVE / MAG /
POL / BRA / TUR

Default: ENG

Language of the firmware.

Time

00:00 – 23:59 (hh:mm)

Sets the time of the internal clock. When T-RAC II is connected to TOMEs via the network, the time is automatically taken from TOMEs. In that case, this setting cannot be changed.

Date

01-01-2010 – 31-12-2099 (DD-MM-YYYY)

Sets the date of the internal clock. When T-RAC II is connected to TOMEs via the network, the date is automatically taken from TOMEs.

Manage progs

Both / TOMs

Default: Both

Both: Programs can be managed on both T-RAC II and TOMEs.

TOMs: Programs can only be managed on TOMEs. On T-RAC II, it is not possible to edit / duplicate / delete programs or to change the order of the programs. This option provides greater certainty that all devices contain the same programs.

The programs can still be loaded/saved via the USB flash drive.

TOMEs short

1 – 999

Default: 1


Only useful when using a TOMEs network. On TOMEs, T-RAC II devices are referenced with their TOMEs short numbers instead of their serial numbers. It is a unique number for each T-RAC II device in the network. When two devices with the same number get connected, TOMEs will automatically give one of the devices a new TOMEs short number.

[NOTE] Changing this settings is disabled while T-RAC II is connected to TOMEs. Disconnect from TOMEs (e.g. by setting the program setting 'Enable LAN' to 'No') before changing this setting.

7.2 Programs

T-RAC II Jr has one program which cannot be transferred to other devices.


T-RAC II Ethernet and Wifi can contain up to ten programs. The programs are interchangeable with other devices and with TOMEs via the TOMEs network or via a USB flash drive. At startup, T-RAC II will always load the first program in the program list (program number 1). Switching to another program is possible when the event 'Program #' is provided (chapter 6.2.4).

[NOTE] When the device setting 'Manage progs' is set to 'TOMs', all options on this screen are disabled. In that case, it's only possible to view the programs or to transfer programs via the TOMEs network or via the USB flash drive.

Copying programs from/to the USB flash drive

From the programs screen, press . A new screen is displayed:

Press 'USB > T-RAC' in order to overwrite the internal programs of the current T-RAC II device with the programs of the USB flash drive. The batch of programs on the USB flash drive must have the name 'TRACPROG.RAP'.

Press 'T-RAC > USB' in order to write all programs currently present in T-RAC II to the USB flash drive. The programs are stored in one file carrying the name 'TRACPROG.RAP'. This file can be imported in other T-RAC II devices or in TOMEs.


[WARNING] The USB flash drive may only be inserted or removed from T-RAC II when T-RAC II is switched off. Otherwise, donation data may get lost and the USB drive may get damaged.

[WARNING] Never switch off T-RAC II while data is being transferred. When data is being transferred, this is visualized on the touch screen by a green blinking USB-symbol.

[NOTE] It is not possible to copy a single program. All ten programs will be exported / imported at once. Existing programs will be overwritten!

Contents of programs

The contents of the programs are explained in the following chapters. They can be reached by selecting a program and pressing 'edit'. A new screen will be displayed:

7.2.1 Program settings

Following program settings can be found on screen 1:

Name of program

20 characters

This name will be shown in the program overview and in the event 'Program #'.

Updated on

Information

Moment this program was edited last time.

Alarm volume

1 / 2 / 3 / 4

Default: 2

Audible volume of the alarms.

Button volume

0 / 1 / 2 / 3 / 4 Default: 1
 Audible volume when pressing a button on the touch screen or on the remote control. When set to 0, no audible signal is given.

All alarms aud

Yes / No	Default: Yes
Yes:	All alarms are audible, including the low flow rate pre-alarm and alarms that happen when not in donation.
No:	Only alarms that happen during donation are audible. The low flow rate pre-alarm is not audible.

End alarm

Off / Mel / MelR / Mon / MonR	Default: Mel
Off:	No audible alarm at the end of donation.
Mel:	Melody (1x)
MelR:	Melody (repeating until user shuts it down)
Mon:	Monotone (1x)
MonR:	Monotone (repeating)

Bcklight lev.

0 / 1 / 2	Default: 0
0:	Backlight level of the touch screen is regulated by T-RAC II. When connected to the mains, the level is at maximum and when only on battery power, the level is at minimum. <i>Also another power-saving feature is activated: when running only on battery power, the backlight of the touch screen will be dimmed after 30 seconds of inactivity on the idle screen. Touch the touch screen or the tray or press the buttons of the multi-functional pole to illuminate the screen again.</i>
1:	Backlight level: dimmed.
2:	Backlight level: bright.

Prog. barcode

Barcode (max 25 characters)

This barcode can be read during the event 'Program #'. T-RAC II will then automatically select the program which is linked to this barcode. Because this barcode selects the program, it must be different for each program.

Don. Organ.

20 characters

The content of this setting will be recorded in the data of each donation. It is meant to contain the name of the blood donation organization. In the dataview on TOMEs, it can be found with the name 'Donation organization'.

Following program settings can be found on screen 2 (configuration of event 'volume chooser'):

Minimum vol.

1 – 650ml Default: 1
 Minimum volume which can be chosen in the event 'volume chooser' or which can be received from BBIS/TOMEs.

Maximum vol.

2 – 650ml Default: 500
 Maximum volume which can be chosen in the event 'volume chooser' or which can be received from BBIS/TOMEs.

Choose vol.

Fix / BBIS / Barc / Prst / Keyb	Default: Keyb
Fix:	The donation volume of preset 1 is chosen automatically. No 'volume chooser' event will be displayed before donation.
BBIS:	<i>The donation volume is obtained from BBIS/TOMEs. A volume confirmation message will be shown on the screen instead of the volume chooser event.</i> <i>A network connection with TOMEs is required for this to work, otherwise the donation will need to be cancelled. Also, special configuration for bi-directional communication is required, please refer to the TOMEs manual for more information.</i>
Barc:	<i>The 'volume chooser' event is activated. The donation volume can only be chosen via one of the barcodes 'Barcode vol.'</i>
Prst:	<i>The 'volume chooser' event is activated. The donation volume can be selected on the touch screen using one of the preset buttons, or via the barcodes 'Barcode vol.' if they are defined.</i>
Keyb:	<i>The 'volume chooser' event is activated. The donation volume can be selected on the touch screen using one of the preset buttons, by using the keyboard (via the 123 button) or via the barcodes 'Barcode vol.' if they are defined.</i>

Preset vol. 1 – Preset vol. 4

0 – 650ml Default: 400, 450, 500 and 0ml
 Here the donation volume presets are defined which can be chosen in the event 'volume chooser'. When a preset is set to 0, it will not be shown in the event. The volumes which can be chosen here, are not limited to the limits set in the settings 'Minimum vol.' and 'Maximum vol.'. The first preset cannot be set to 0.

Barcode vol. 1 – Barcode vol. 4

Barcodes (max 25 characters)
A barcode can be assigned to each donation volume preset. When the setting 'Choose vol.' is set to Barc, Prst or Keyb, the barcode can be read to select the corresponding volume.

Following program settings can be found on screen 3:

Low flow limit

1 – 200ml/min Default: 40ml/min
 Minimum flow rate for triggering the low flow rate alarm and the low flow rate pre-alarm.

Low flow time

0 – 60s

Default: 15s

Time in seconds the low flow rate pre-alarm is activated. When the flow rate is still too low after this time, the low flow rate alarm will be activated. When set to 0, the pre-alarm is disabled and the low flow rate alarm will be activated immediately in case of a low flow rate.

High flow limit

10 – 300ml/min

Default: 200ml/min

When the flow rate passes above this limit during donation, the clamp will close and the high flow rate alarm will be triggered.

Cont high flow

Yes / No

Default: No

Yes: After a high flow alarm, the user will be able to continue the ongoing donation.

No: After a high flow alarm, the ongoing donation will be discontinued.

Max. don. time

00:01 – 30:00 (mm:ss)

Default: 15:00

Sets the maximum donation time. When the elapsed time passes this limit, the donation is stopped automatically, the clamp is closed and an alarm is triggered.

Warning timer

00:00 – 30:00 (mm:ss)

Default: 12:00

A special time warning can be given on the touch screen using this setting. When the elapsed time passes the limit, the 'time elapsed' as shown on the touch screen turns red.

When the estimated time to finish the donation passes the limit, the 'estimated time left' turns red.

Using this setting, it can be visualized during donation whether or not the blood will be usable for making platelet concentrate. When setting the timer to 00:00, the warning timer is disabled.

Time estim.

Left / Tot

Default: Left

Left: During donation, the estimated time before finishing the donation is displayed. The value is calculated based on the flow rate, the total donation volume and the time elapsed.

Tot: The estimated total donation time is shown.

Enable pause

Yes / No

Default: No

Yes: Special option for labeling of the bag during donation. Between 1% and 80% of the donation (based on donated volume), a pause button is displayed. When pressed, the clamp will close and the user has 60 seconds to put labels on the blood bag. After that time, or when pressing 'resume', the donation continues automatically. The pause function is only available once during a donation.

No: The pause option is not available. A donation can still be paused by pressing 'Stop' and 'Continue'.

Mix freq. start

10 – 30 oscillations per minute Default: 30
 Sets the agitating frequency of the tray at the start of donation.

Mix freq. end

10 – 30 oscillations per minute Default: 30
 Sets the agitating frequency of the tray at the end of donation. If this frequency is different from the frequency set in 'Mix freq. start', the oscillation frequency will change gradually during donation.

Start agitator

0 – 60 ml Default: 0
 Volume to be collected before the tray starts agitating. When set to 0, the agitator starts immediately.

Keep mixing

Yes / No	Default: Yes
Yes:	After donation, the tray keeps on agitating until the clamp goes open.
No:	After donation, the agitating is stopped.

Following program settings can be found on screen 4:

Stats. timer

Yes / No	Default: Yes
Yes:	The statistics screen proceeds automatically to the idle screen after 8 seconds.
No:	The user needs to press the 'finish' button on the statistics screen in order to proceed.

Blood density

1.05 – 1.06g/ml Default: 1.06
 The blood volume which is displayed on T-RAC II is based on the weight of the blood and the blood density. The blood density can be adjusted here.

Proceed weigh

0 – 300g Default: 20
 When the idle screen is displayed, the screen continues automatically to the next screen when a bag is placed in the tray. 'Proceed weigh' is the weight difference needed to trigger this. 0 disables auto-proceeding, the user needs to press 'next' in order to proceed.

Auto-clamp

Yes / No	Default: Yes
Yes:	When tube is inserted in the clamp before donation, the clamp locks after pushing the clamp head backwards and then releasing it.
No:	The user needs to press the button 'lock clamp' on the touch screen for locking the clamp.

PVFlow. delay

Default: 00:00

00:00 – 05:00 (mm:ss)

Sets the interval between venipuncture (moment of pressing 'start' button) and the moment blood should flow into the bag. If this time is exceeded, an alarm will sound. If set to 00:00, the post venipuncture flow delay alarm is disabled and no audible low flow alarm will be given until the first blood is collected.

Go PVFlow. de

Yes / No

Default: Yes

- Yes: When the post venipuncture flow delay time is exceeded, the user can still continue the donation.
- No: When the time set in PVFlow. delay is exceeded, the donation will be cancelled automatically.

Clamp interval

Yes / No

Default: No


- Yes: After donation, the clamp will go to the 'locked' state for 10 seconds each minute. In this way, more blood might be donated but the risk of coagulation in the area of the venipuncture is less.
- No: After donation, the clamp stays closed until the user presses 'open clamp'.

Time tubes

Yes / No

Default: No

- Yes: *When collecting blood in tubes before donation, this is normally done on the venipuncture screen. On that screen, an extra button will be shown which should be pressed between venipuncture and the collecting of the blood tubes. In this way, both actions can be timed. These durations are recorded into the donation data.*


- No: *No extra button or timer is activated.*

Re-scan codes

Yes / No

Default: No

- Yes: *On the statistics screen at the end of donation, the user can go back to view and edit the barcodes that have been read throughout the donation process.*
Once the idle screen is reached, the barcodes cannot be changed anymore.


- No: *The inputs of events cannot be edited anymore when the statistics screen is reached.*

BBIS Program

Yes / No

Default: No

- Yes: *The donation program can be chosen by the BBIS. The program which is chosen, will be displayed after the BBIS code has been read. The user needs to accept the chosen program, otherwise the donation will be cancelled automatically.*

- No: *The BBIS can still select the collection volume and other things, but not the donation program.*

- [NOTE] *This setting requires an active network connection to work. TOMEs and the BBIS should be set-up to communicate with each other.*

Enable LAN

Yes / No

Default: Yes

Yes: Ethernet is activated.

No: Ethernet is not activated. None of the features requiring a network connection (to TOMEs) will work. The network status icon is shown as

It is possible to quickly change this setting as well from the LAN status screen.

In order to do so, press the Ethernet status icon

On T-RAC II Wifi, disabling the LAN functionality will also switch off the internal wifi module, saving power.

Left / right

Off / Ask / Left / Rght

Default: Off

Off: The device doesn't send information towards TOMEs whether the device is set-up for donation with left or right arm.

Ask: After every startup of T-RAC II, the device will ask on the idle screen whether the device is set-up for donation with left or right arm. This information will be sent to TOMEs and can be viewed in the device status screen on TOMEs for easier assignation of donors to the donation beds.


Left: The device will send information to TOMEs that the device is set-up for donation with left arm.

Right: The device will send information to TOMEs that the device is set-up for donation with right arm.

[NOTE] This setting requires an active network connection to work.

7.2.2 Program flow

On screen 5, the order of the events can be viewed. This program flow can only be changed on TOMEs.

7.2.3 Barcode rules

On screen 6, the way barcodes are compared to each other can be configured. The checks are done after the reading of each code in a donation cycle. When an error towards the scheme is detected, an error message will be displayed. T-RAC II won't accept the wrong barcode, another one will need to be scanned. A maximum of 20 equations can be programmed.

Two comparators are available:

- <>
The barcodes at each side of the equation should be different.
- =
The barcodes should be equal.

In each equation, two barcodes can be chosen. Following codes can be chosen:

- **Codes present in the program**

The available codes are the ones that are present in the program which is currently being edited.

- **BBIS CODE 1 to BBIS CODE 5**

These are codes which can be configured through BBIS. These can only be used when the program is configured for bi-directional communication in TOMEs. Depending on the version of TOMEs the wording of these codes might differ from what is shown on T-RAC II:

T-RAC II	TOMEs
BBIS CODE 1	BBMS CODE 1 / BBIS IDENTIFIER
BBIS CODE 2	BBMS CODE 2 / BBIS CODE 1
BBIS CODE 3	BBMS CODE 3 / BBIS CODE 2
BBIS CODE 4	BBMS CODE 4 / BBIS CODE 3
BBIS CODE 5	BBMS CODE 5 / BBIS CODE 4

- * ('all barcodes which are scanned during the donation cycle')

When using '*' along with the comparator '<>', the code on the other side of the equation should be different from all other codes read during the donation cycle. Using this character along with the '=' comparator, is only useful when (a part of) all barcodes should be equal.

- ... ('empty')

When at least one of both sides of the equation is '...', the equation is not taken into account. Use this to delete a barcode comparison.

It is not always required to compare the complete result of the barcodes. Following parts can be chosen:

- **Full**

The complete barcodes will be compared.

- **Prfx**

Only the start flag of the codes (prefix) will be compared.


- **Sufx**

Only the end flag of the codes (suffix) will be compared.

- **Cntr**

Only the center part of the codes (without the start and end flag) will be compared.

The lengths of the flags are configured in the barcode event settings on TOMEs. Only parts of the same length can be compared to each other.

[NOTE] In the barcode event settings on TOMEs, it is possible to program additional barcode checks.

7.3 Communication settings

The communication settings define how T-RAC II should connect to TOMEs or to TOMEs Mobile Collection Tool. These settings can only be changed on each T-RAC II device individually. Following settings are available:

Static/dyn IP

- stat / dyn Default: dyn
- stat: T-RAC II will use a static IP-address. This IP-address, the subnet mask and the gateway IP-address can be defined in the settings 'Static IP', 'Subnet mask' and 'Gateway IP'.
- dyn: T-RAC II will receive an IP-address automatically from the router via DHCP. The settings 'Static IP', 'Subnet mask' and 'Gateway IP' are not available.

Static IP

000.000.000.000 – 255.255.255.255 Default: 192.168.001.002

This setting is only available when the option 'Static/dyn IP' is set to 'stat'. It defines the IP-address of the T-RAC II.

Subnet mask

000.000.000.000 – 255.255.255.255 Default: 255.255.255.000

This setting is only available when the option 'Static/dyn IP' is set to 'stat'. It defines the subnet mask of the T-RAC II.

Gateway IP

000.000.000.000 – 255.255.255.255 Default: 192.168.001.001

This setting is only available when the option 'Static/dyn IP' is set to 'stat'. It defines the IP-address of the gateway the T-RAC II will connect to.

TOMEs IP

000.000.000.000 – 255.255.255.255 Default: 192.168.1.100

This setting defines the IP-address of the TOMEs device or the PC running TOMEs Mobile Collection Tool that T-RAC II should communicate with.

Port

0 - 65534

Default: 21010

This setting defines the IP port number that T-RAC II should use for communication with TOMEs.

[NOTE] Depending on the changed communication settings, T-RAC II might switch off itself after pressing 'SAVE'. In that case, press the ON/OFF button to restart the device.

8 Alarms and error messages

Various alarms can be triggered by T-RAC II. Alarms that happen during donation are always audible. Whether the other alarms are audible, can be configured by the program setting 'All alarms aud'.

At some alarms, following button is shown: . This button will only silence an alarm temporary without disabling it. After 45 seconds, the alarm will sound again if the problem is not resolved.

The duration of alarms which take place between venipuncture and needle removal, are visualized on the screen *and recorded into the donation statistics*.


Medium priority alarms are accompanied by an acoustic signal (c-c-c), a low priority alarm is accompanied by an acoustic signal (e-c) (according to IEC 60601-1-8).

8.1 Donation alarms

Following alarms can only happen during donation. They are triggered based on limits set by the user in the program settings. All these alarms are medium priority alarms, except for the low flow rate pre-alarm.

Low flow rate pre-alarm

Indicated by: Touch screen: flow rate indication turns red.

Clamp: stays locked.

Antenna: color doesn't change, but starts flashing.

Audible alarm: only if the program setting 'All alarms aud.' is set to 'Yes'.

Cause: The flow rate is lower than the flow rate limit configured in the currently running program.

Solution: The flow rate should increase. The donor can react on the alarm given by the antenna by squeezing his hands to raise the blood flow.

Low flow rate alarm

Indicated by: Touch screen: red indication 'LOW FLOW'.

Clamp: stays locked.

Antenna: fast red flashing

Cause: The flow rate has been lower than the flow rate limit, during a longer time than the low flow time. Both settings are configured in the running program.

Solution: The flow rate should increase.

High flow rate alarm

Indicated by: Touch screen: message 'HIGH FLOW RATE'.

Clamp: closes automatically.

Antenna: fast red flashing

Cause: The flow rate was higher than the high flow limit, configured in the currently running program.

Solution: Possibly the venipuncture needle is inserted in an artery. Precautions need to be taken to prevent the donor from losing too much blood, fainting etc. Depending on the setting 'Cont high flow', the donation can be continued after a high flow rate alarm or not.

Donation finished: volume reached

Indicated by: Touch screen: flashing message 'DONATION FINISHED'. Clamp: closes automatically.
Antenna: red flashing.

Cause: The collection volume is reached.

Solution: Press 'Next' to continue.

Donation finished: time reached

Indicated by: Touch screen: flashing message 'MAXIMUM TIME REACHED'. Clamp: closes automatically.
Antenna: red flashing.

Cause: The maximum donation time is reached.

Solution: Press 'Next' to continue.

Post venipuncture flow delay exceeded

Indicated by: Touch screen: message 'PVFD. EXCEEDED'. Clamp: closes automatically.
Antenna: fast red flashing.

Cause: No blood flow was detected for a period of time longer than the period of time defined in the program setting 'PVFlow. delay'.

Solution: Depending on the program setting 'Go PVFlow. de', the donation can be continued or not.

No tube in clamp

Indicated by: Touch screen: message 'NO TUBE IN CLAMP'.

Clamp: closes automatically.

Cause: No tube was detected inside the clamp. This was checked by closing the clamp for 15 seconds during donation. During the check, weight was added to the tray.

Solution: Open the clamp by pressing the clamp status icon  and insert the tube correctly. Press 'Recheck' to check again whether the tube is inserted correctly. Press 'Continue' in order to proceed without rechecking.

8.2 System alarms

At following alarms, the antenna will flash red fast (medium priority alarm) when a donation is in progress. Otherwise, the antenna will flash white (low priority alarm).

Low battery level

Indicated by: Touch screen: message 'LOW BATTERY'.
Cause: If shown during donation:
The battery is almost empty. The clamp will close automatically.
If shown when not in donation with high-energy battery:
The battery level is too low, no donations can be started.
If shown when not in donation with emergency battery:
T-RAC II Jr is not connected to the mains.
Solution: Connect T-RAC II to the mains.
If a donation is ongoing and no mains power is available, finalize the donation.
If no donation is ongoing, an empty high-energy battery can be replaced by a charged one. Switch-off the device before replacing the battery, unless it is connected to the mains.

Battery error

Indicated by: Touch screen: message 'BATTERY ERROR'.
Cause: A problem with the battery has been detected.
Solution: In case of high-energy battery:
Connect T-RAC II to the mains or replace the battery. Switch-off the device before replacing the battery, unless it is connected to the mains.
In case of emergency battery (only on T-RAC II Jr):
Contact your Terumo BCT distributor.

Agitator error

Indicated by: Touch screen: message 'AGITATOR ERROR'.
Cause: T-RAC II detected an error when trying to swing the agitator. The agitator may be blocked or the weight on the scale may be too large.
Solution: Check the agitator. If the problem persists without obvious reason, contact your Terumo BCT distributor.

Scale error

Indicated by: Touch screen: message 'SCALE ERROR'.
Cause: This error can be caused by several reasons:

- It was not possible to tare the balance because its weight was too high, too low, or not constant.
- The scale is out of calibration.
- The balance indicated a negative weight.

Solution: Check the balance visually to see if there are no abnormalities. If the problem persists without obvious reason, contact your Terumo BCT distributor.

Clamp error

Indicated by: Touch screen: message 'CLAMP ERROR'.
Cause: T-RAC II detected an error when trying to move the clamp. The clamp may be blocked or broken.
Solution: Remove the clamp using the emergency exit feature (chapter 5.2) and check it visually to see if there are any abnormalities. If the problem persists without obvious reason, contact your Terumo BCT distributor.

System alarm

Indicated by: Touch screen: message 'SYSTEM ALARM' + code.
Cause: Cause depends on the code of the alarm.
Solution: Contact your Terumo BCT distributor for more information.

Following alarms won't happen when a donation is in progress. *The antenna will flash white at these alarms (low priority alarm).*

Tray abnormally underloaded / overloaded

Indicated by: Message on the touch screen during self-check.
Cause: At startup, the weight of the tray is checked. The device detects an abnormality of the tray weight.
Solution: Check whether the tray is present and whether something is laying on it.

No battery present at startup

Indicated by: Message on the touch screen.
Cause: This warning informs the user that no battery is present at startup. When the mains power is lost, the device will shutdown automatically without closing the clamp.
Solution: In case of high-energy battery:
Put a battery in place to prevent a loss of power due to a mains interruption.
In case of emergency battery (only on T-RAC II Jr):
Contact your Terumo BCT distributor.

Memory 80% full

Indicated by: Touch screen: message 'MEMORY 80% FULL'.
Cause: *The internal memory which is used for saving donation data, is almost full.*
Solution: *Connect T-RAC II to TOMEs via the network, or switch off the device, insert a USB flash drive in the USB slot at the front and switch on the device again. When the device is idle, it will automatically copy the donation data to TOMEs / USB flash drive. Data which is successfully copied, is automatically erased from the internal memory.*

Memory full

- Indicated by:* Touch screen: message 'MEMORY FULL'.
- Cause:* The internal memory which is used for saving donation data, is full. When a new donation would be started, the oldest donation data will get lost.
- Solution:* Connect T-RAC II to TOMEs via the network, or switch off the device, insert a USB flash drive in the USB slot at the front and switch on the device again. When the device is idle, it will automatically copy the donation data to TOMEs / USB flash drive. Data which is successfully copied, is automatically erased from the internal memory.

Device settings changed by TOMEs

- Indicated by:* Touch screen: message 'Device settings were changed by TOMEs'.
- Cause:* TOMEs sent new device settings to the device. They need to be approved before the device can continue.
- Solution:* Check whether the new settings are OK by pressing 'Settings'. When the settings are OK, approve them by pressing 'OK'.

Programs changed by TOMEs

- Indicated by:* Touch screen: message 'Programs were changed by TOMEs'.
- Cause:* TOMEs adapted the programs of the device. They need to be approved before the device can continue.
- Solution:* Check whether the new settings are OK by pressing 'Settings'. When the settings are OK, approve them by pressing 'OK'.

TOMEs short number changed

- Indicated by:* Touch screen: message 'The TOMEs short number of this device was changed by TOMEs to #XX'.
- Cause:* TOMEs adapted the TOMEs short number of the current device. The new short number is depicted on the touch screen.
- Solution:* Approve by pressing 'OK'.

USB flash drive error

- Indicated by:* Touch screen: message 'USB flash drive error'.

Status icon:  icon appears even when a USB drive is connected.

- Cause:* T-RAC II cannot handle the USB flash drive which was inserted. Donation data will not be written to the drive.
- Solution:* Use only the USB flash drive which was delivered with the device. If that doesn't work, contact your Terumo BCT distributor.

8.3 Event warnings

Barcode events can be configured to check the scanned barcodes. When a wrong code is read, an error will be given. Also, when communicating with the BBIS, special warnings may be given. The antenna will flash white when one of the following errors/warnings occurs.

Barcode doesn't match with the barcode rules

- Indicated by:* Message on the touch screen.
- Cause:* The scanned barcode is compared to another barcode which has been scanned earlier. The comparison is not correct according to the barcode rules (chapter 7.2.3).
When using bi-directional communication, it is possible that the barcode was compared to a barcode given by the BBIS.
- Solution:* Scan a correct barcode or adjust the barcode rules.
Depending on the setting of the barcode event as configured in TOMEs, it is possible to skip it and go to the next event. In that case, the faulty barcode won't be saved into the donation data.

Code flags do not match

- Indicated by:* Message on the touch screen.
- Cause:* On TOMEs, it is possible to define the first and last 10 characters of a barcode that have to be checked. When the scanned barcode doesn't match with this setting, the error will be given. The error is also given in case a multi-barcode event is configured to be completely equal or unequal and an error occurs against this rule.
- Solution:* Scan a correct barcode or adjust the barcode event settings in TOMEs.
Depending on the setting of the barcode event, it is possible to skip it and go to the next event. In that case, the faulty barcode won't be saved into the donation data.


Bi-directional communication: waiting for info from the BBIS

- Indicated by:* Message on the touch screen.
- Cause:* The event has been configured to receive information from the BBIS using bi-directional communication, but the link with the BBIS has not been established yet.
- Solution:* Wait until the connection is established. Make sure that the Ethernet connection to TOMEs is active. If the problem persists, check in TOMEs whether a donation input of the scanned BBIS Code 1 is present in Column BBMS Code 1 or BBIS Identifier in TOMEs (depending on which TOMEs version is installed).

Bi-directional communication: info from BBIS

Indicated by: Message on the touch screen.

Cause: Via the BBIS, a program and/or a volume was set to be used during the donation. The program and/or volume is depicted on the screen.

Solution: Check whether the program and volume are OK. If so, continue the ongoing donation by pressing . When pressing , the donation will automatically be cancelled.

Bi-directional communication: info from BBIS not plausible

Indicated by: Message on the touch screen.

Cause: Via the BBIS, a volume was set to be used during the donation. The volume which was selected on the BBIS, is not within the range allowed in the program on T-RAC II.

Solution: The information from the BBIS needs to be resent (to be performed on the BBIS). Otherwise, the donation will be cancelled.

Bi-directional communication: donor approval not OK


Indicated by: Message on the touch screen.

Cause: On the BBIS, the doctor has refused the donor to donate.

Solution: Press 'Continue', the donation will be cancelled automatically.

9 Troubleshooting

For alarms and error messages on the touch screen, please refer to chapter 8 . Other problems and their solutions are described below.

Symptom	Cause	Solution
The device doesn't start up.	1. The device is not connected to the mains or the internal battery is missing, broken or empty.	Connect the power supply to the mains and to T-RAC II. If the device is connected to the mains, the LED on the on/off button should be illuminated. If this is not the case, contact Terumo BCT for servicing. When working on a high-energy battery and no connection to the mains is available, replace the battery by a charged one.
	2. The device is powered-off.	Power-on the device by pressing the on/off button. The touch screen should become illuminated. If this is not the case, contact Terumo BCT for servicing.
The touch screen doesn't react correctly	The touch screen is out of calibration.	 <p>In case the touch screen doesn't respond correctly anymore, it needs to be calibrated. In order to do this, switch off the device. Then, press the on/off button and push the clamp backwards at the same time. This will enable the touch screen calibration mode. Calibration is done by pressing on the touch screen at the spots indicated by the arrows.</p>

10 Cleaning and maintenance

10.1 General cleaning

Always keep T-RAC II clean.

Clean the surface of T-RAC II with a soft cloth soaked in mild detergent or a chlorhexidine gluconate solution of less than 0.5% and well squeezed. When cleaning T-RAC II, be careful not to spill any cleaning solution inside the instrument. If this happens, do not operate the instrument. Professional cleaning and service is required if any liquid (blood, blood product, cleaning solution or other fluid) is spilled into the internal components of the machine.

[WARNING] Before cleaning the unit, disconnect the device from the mains by plugging out the power cable.

[WARNING] Abrasive solutions such as acetone or ammonia and organic solvents (except for following alcohols) shall not be used! Allowed alcohols: ethanol and isopropanol.

[WARNING] Do not use a moist cloth to clean the electrical contacts for the battery and do not spill any liquid on these contacts. This can create short circuits and damage the battery irreversibly.

10.2 Preventive maintenance

10.2.1 Check of the weighing function

Terumo BCT recommends to check the weighing function of the device weekly or every time the device is moved (e.g. for mobile collection). This can be done by:

1. Switching on T-RAC II until the idle state is reached (see chapter 6.2).
2. Pressing the 'TARE' button on the touch screen.
3. Placing a calibrated weight of minimum 200g and maximum 500g on the tray and checking its weight. The maximum allowable deviation from the weight is ±3g. If this is not the case, please have your T-RAC II device serviced.

10.2.2 Check of the alarm function

Regularly check the alarm function in following way:

1. Go to the settings menu and make sure that the program setting 'All alarms aud.' is set to 'Yes'. Refer to chapter 7 for more information.
2. Switch-off T-RAC II by pressing the on/off button and holding 3 seconds.
3. Pull-out the battery and make sure T-RAC II is still connected to the mains.
4. Switch-on T-RAC II by pressing the on/off button.
5. An alarm should sound and the antenna LED should be illuminated.

10.2.3 General maintenance

To assure safe and trouble free operation, please have your T-RAC II serviced at least every 12 months, depending on the frequency of use. If you experience a system failure, please call Terumo BCT or your Terumo BCT representative for servicing. When calling for servicing, please be able to describe the problem and have the serial number available so that faster handling of the repair will be possible. Do not attempt to dismantle or repair the unit. Only qualified personnel are authorized to service the machine.

Technical support - Medical Electronics
Call **free** number from

A	0800-293711	B	0800-94410	DK	808-80701
SF	0800-115226	IRL	1800-553224	CH	0800-563694
F	0800-908793	I	800-785891	UK	0800-9179659
D	0800-1808183	NL	0800-0222810	N	0800-12270
GR	00800-3212721	E	900-963251	S	020-791373

Other countries call +32 16 39 14 00 at international rates

E-mail: meservice@terumobct.com

11 High-energy battery shipping information

IMPORTANT: WHENEVER POSSIBLE, REMOVE THE HIGH-ENERGY BATTERY WHEN SHIPPING THE T-RAC II!

The T-RAC II high-energy battery (8ME05TRC226) and high-energy battery XLR (8ME05TRC320) are Li-Ion batteries with a nominal energy of 51Wh (2 cells, >20Wh/cell) and 49.6Wh (2 cells, >20Wh/cell). It has been tested successfully in accordance with IEC 62133 and UN38.3 (UN Manual of Tests and Criteria Part III Subsection 38.3), and is subject to the requirements of the IATA Dangerous Goods regulations (DGR). If contained in the T-RAC II, it is classified in Class 9 – Miscellaneous Dangerous Goods as UN3481 (Li-Ion batteries contained in equipment) and the requirements contain the following:

- The T-RAC II containing the high energy battery 8ME05TRC226 or high energy battery XLR 8ME05TRC320 must be packed in a strong outer packaging in accordance with Packing Instruction 967, Section I (DGR Section 5).
- The package displays at least following labels / markings / documents (DGR Section 7):
 1. Miscellaneous Dangerous Goods Class 9 label:

2. UN3481 shipping name label:

UN3481
 Lithium ion batteries contained
 in equipment, 0.5kg
3. Shipper's name and address marking.
4. Consignee's name and address marking
5. Dangerous Goods Declaration (DGR Section 8).
- Dangerous Goods training for the personnel involved in the shipping (DGR Section 1.5).

The original packaging of the T-RAC II complies with PI967 and contains the above mentioned labels, which are only necessary when the T-RAC II is shipped together with the high-energy battery 8ME05TRC226 or high-energy battery 8ME05TRC320. At Terumo BCT, the personnel involved in the shipping is trained in Dangerous Goods and is therefore allowed to make the Dangerous Goods Declaration. As this is usually not the case for other persons that might need to ship a T-RAC II:

Terumo BCT strongly recommends to ship the T-RAC II always WITHOUT the high energy batteries, in its original packaging WITHOUT the Class 9 label and the UN3481 label (i.e. labels 1 and 2 shown above need to be removed from the package or covered).

The Class 9 and UN3481 label need to be removed or covered because otherwise the shipper or customs will refuse or block the package because they will assume the package contains a Li-Ion battery but the necessary documents are missing.

Please also note that the shipping of damaged or defective batteries is strictly forbidden. Thus, if a battery is suspected to be defective, it should not be shipped under any circumstances.

If for some reason, the T-RAC II needs to be shipped with the high-energy battery, it is the sole responsibility of the shipper to comply with all requirements of the DGR, of which the most important ones are mentioned above. The most recent information on DGR and packing instructions of Li batteries can be found on:

<http://www.iata.org/whatwedo/cargo/dgr/Pages/lithium-batteries.aspx>

12 Disposal (end of life)

Electrical and electronic equipment (EEE) and batteries contain materials, components and substances which can be dangerous to the environment and harmful to human health if waste electrical and electronic equipment (WEEE) and batteries are not disposed of correctly.

Waste electrical and electronic equipment and batteries must not be disposed of with the remainder of unseparated waste, but should instead be collected separately. In this way, the environmental impact associated with the disposal of WEEE and batteries is reduced and there will be more opportunity for reusing, recycling and recovering WEEE and recycling batteries.

The Li-ion or the NiMH battery should be removed from the equipment.

At end of life, please dispose of this device/battery according to your local regulations. Contact your local distributor or municipality to know the available collection schemes. The embedded back-up coin-type lithium battery, collected together with this equipment at end of life as WEEE, is to be removed and treated by the recycling center.

T-RAC II has a life span of approx. 10 years.

[NOTE:] Warranty does not apply to any consumable part like a battery.

13 T-RAC II Wifi setup


Each T-RAC II Wifi device needs to be set-up separately in order to communicate with the wireless access point. To do so, proceed as follows

1. Put the IP-address of an Ethernet port of a PC on 192.168.1.100 and the subnet mask on 255.255.255.0.
2. Connect the Ethernet port of T-RAC II which is marked as follows directly to the Ethernet port of the PC.

[NOTE] If the PC doesn't support Auto-MDIX, a crossover cable is needed for this.

3. Switch-on T-RAC II Wifi.
4. On the PC, open a web browser and surf to 192.168.1.253. Following screen should be shown:

5. Enter 'dpac' in both the username and password fill-in box and press enter.

6. On screen 'Module Status' press 'Configuration'.

Module Status

Module Firmware Version:	4.00
Radio Firmware Version:	3.1.64.49
Link Status:	Connected
SSID:	linksys
MAC Address:	000B280D6963
BSSID:	C0C1C049F4AD
Transmit Rate (Mb/s):	54
Signal Level (dBm):	-41
Noise Level (dBm):	-96
IP Address:	192.168.1.103
Subnet Mask:	255.255.255.0
Default Gateway:	0.0.0.0
Primary DNS:	0.0.0.0
Secondary DNS:	0.0.0.0
Up Time (Sec):	34

7. A new screen is shown. Fill-in the SSID of the wireless access point. Afterwards press 'Commit':

WLAN Parameters

Radio Startup Mode:	On
WLAN Connection Type:	Infrastructure
SSID:	linksys
WLAN TX Power (dBm):	15 dBm
Maximum Wireless Data Rate:	Auto
Specify Specific Data Rates:	0x00000000
Use Fixed Data Rate:	Disabled
WLAN Region:	United States (US)
Preferred Radio Band:	Auto
WLAN Maximum Retries:	13
WLAN RTS Threshold:	0

Commit

8. After the configuration changes were committed successfully, press 'Reload' in this screen:

Configuration changes committed successfully

Reload

9. Press 'WLAN Security Settings' on the left of the screen:

The screenshot shows the B&B Electronics web interface. The top navigation bar includes links for Status, Configuration, Certificates, Network, and Maintenance, along with a note that the device is an Ethernet Router. The main content area is titled 'WLAN Parameters' and contains various configuration fields such as Radio Startup Mode (On), WLAN Connection Type (Infrastructure), SSID (linksys), WLAN TX Power (15 dBm), Maximum Wireless Data Rate (Auto), and WLAN Region (United States (US)). At the bottom of this section are 'Commit', 'Cancel', and 'Defaults' buttons. On the left sidebar, under the 'Configuration' section, the 'WLAN Security Settings' link is circled in red.

10. A new screen is shown. Fill in the security settings of the wireless access point and press 'Commit':

This screenshot shows the 'WLAN Security Parameters' configuration page. It includes fields for WLAN Security Type (WPA2-PSK-TKIP) and Authentication Type (Auto). Below these is a 'WPA / WPA2 Settings' section with a 'Pre Shared Key (PSK)' field. At the bottom of the page are 'Commit', 'Cancel', and 'Defaults' buttons. The 'Commit' button is circled in red.

11. After the configuration changes were committed successfully press 'Restart' in this screen:

The screenshot displays a success message: 'Configuration changes committed successfully'. Below this message are two buttons: 'Reload' and 'Restart'. The 'Restart' button is circled in red.

12. Wait 30 seconds while the module restarts:

The screenshot shows a web browser window with the URL <http://192.168.1.253/>. The page title is "B&B ELECTRONICS" and the sub-page title is "Make The Right Connections". The main content area displays the message "Restarting the Module" and "Please wait 30 seconds while the module restarts". Below this, it says "The page will refresh automatically". On the left sidebar, there is a list of configuration options including Express Setup, WLAN Settings, WLAN Security Settings, Network Settings, Serial Port Settings, Serial Port 2 Settings, Connection Settings, Ethernet Settings, Event Settings, Advanced Settings, Upload Configuration File, List Configuration Files, Delete Configuration File, Active Configuration, and User Configuration. At the bottom right, there are links for "B&B Electronics", "Contact Us", and "Technical Support".

13. After the module is rebooted, the screen 'Module Status' will show again. Check the 'link status' to see if T-RAC II is correctly connected to the access point:

The screenshot shows a web browser window with the URL <http://192.168.1.253/>. The page title is "B&B ELECTRONICS" and the sub-page title is "Make The Right Connections". The main content area displays the "Module Status" table. The "Link Status" row is circled in red. The table data is as follows:


Parameter	Value
Module Firmware Version:	4.00
Radio Firmware Version:	0.1.54.49
Link Status:	Connected
SSID:	linksys
MAC Address:	000B280D6963
BSSID:	C0C1C049F4AD
Transmit Rate (Mb/s):	54
Signal Level (dBm):	-41
Noise Level (dBm):	-96
IP Address:	192.168.1.103
Subnet Mask:	255.255.255.0
Default Gateway:	0.0.0.0
Primary DNS:	0.0.0.0
Secondary DNS:	0.0.0.0
Up Time (Sec):	34

At the bottom right, there are links for "B&B Electronics", "Contact Us", and "Technical Support".

14. Fill in the IP-address settings in the T-RAC II Communication settings window (see 7.3), according to your network.

15. Disconnect the Ethernet cable between T-RAC II and the PC.

16. Connect the small Ethernet cable between the two Ethernet ports of T-RAC II:


17. Restart T-RAC II.

18. The Ethernet icon on T-RAC II should now become green, indicating that T-RAC II


Wifi is correctly connected to TOMEs:

[WARNING] Make sure that T-RAC II WIFI is using the frequency band of 2.4 GHz and not 5GHz.


14 Symbol LS-2208 barcode reader set-up

Set All Defaults


Enable Codabar


Low Frequency

OR


*Medium Frequency
(Optimum Settings)


Low Volume


OR


Medium Volume

Baud Rate 2400

*None


*8-Bit

*1 Stop Bit

Scan Options

<DATA> <SUFFIX>

Enter

ACK/NAK

CE 2797

LC-7468en, Rev 8

TERUMO BCT

Terumo BCT, Inc.
10811 W. Collins Ave.,
Lakewood, CO 80215, USA

®: Registered Trademark

©: Copyright Terumo BCT, Inc. – All rights reserved