



Service Manual

Ultra-Low Temperature Freezer

MDF-DU502VH

MDF-DU702VH

Panasonic Healthcare Co., Ltd.
Biomedical Div.

SM0000002-03

Effective model

This service manual is effective for following models.

| Model | Voltage | Frequency |
|----------------|----------------|-----------|
| MDF-DU502VH-PA | 220V | 60Hz |
| MDF-DU502VH-PE | 220V/230V/240V | 50Hz |
| | | |
| MDF-DU702VH-PA | 220V | 60Hz |
| MDF-DU702VH-PE | 220V/230V/240V | 50Hz |

Contents

| | Page |
|---|------|
| Features, Caution ----- | 1 |
| Specifications ----- | 2 |
| -Structural specifications | |
| -Control specifications | |
| -Performance specifications | |
| Dimensions ----- | 5 |
| Cooling unit parts ----- | 7 |
| Refrigeration circuits ----- | 9 |
| Components on PCB ----- | 10 |
| Connections on PCB ----- | 11 |
| Electric parts ----- | 12 |
| Specification of sensor ----- | 13 |
| Wiring diagram ----- | 14 |
| Circuit diagram ----- | 15 |
| Operation ----- | 18 |
| Control specifications ----- | 25 |
| Inverter signal ----- | 33 |
| Parts layout ----- | 34 |
| Repairing unit/Gas charge ----- | 36 |
| Test data ----- | 52 |
| MDF-DU502VH,MDF-DU702VH | |
| -Pull-down & Pull-up temperature | |
| -Pull-down temperature | |
| -Pull-down pressure | |
| -Pull-down current-input | |
| -Temperature uniformity (17points measured) | |
| -Amount of power consumption | |
| -Cycle running | |
| Backup cooling kit installation and setting procedure ----- | 64 |

Features

- Energy saving
- Inverter control of compressor
- Improvement of storage efficiency
- Auto air intake port
- Enhancement of safety functions
- New control panel with LCD touch panel

Caution

*Parts replacement and option unit installation must be done by trained service engineer.
*Service engineer must refer to the section “Electric parts” and “Cooling unit parts”
about the parts for those operations.

Specifications

Structural specifications

| Item | MDF-DU502VH | MDF-DU702VH | | |
|------------------------|---|---|--|---|
| Name | Ultra-Low Temperature Freezer | | | |
| Exterior dimensions | W 790 × D 882 × H 1993 mm | | W 1030 × D 882 × H 1993 mm | |
| Interior dimensions | W 630 × D 600 × H 1400 mm | | W 870 × D 600 × H 1400 mm | |
| Effective capacity | 528 liters | | 729 liters | |
| Exterior | Painted steel | | | |
| Interior | Painted steel | | | |
| Outer door | Painted steel | | | |
| Inner door | 2 | | | |
| Shelf | 3 | | | |
| Outer door latch | 1 | | | |
| Outer door lock | 1 | | | |
| Insulation | Rigid polyurethane foam+VIP PLUS | | | |
| Access port | 3 (back x 1, bottom x 2) | | | |
| Auto air intake port | Inside of control panel (built-in anti freezing heater) | | | |
| Manual air intake port | Product left side lower part | | | |
| Compressor | VNEU213U (single-phase 220-240V 50/60Hz) x 2 | | | |
| Evaporator | High: Cascade, Low: Tube on sheet | | | |
| Condenser | High: Fin and tube, Low: Cascade | | | |
| Refrigerant | High Temp. side | Low Temp. side | High Temp. side | Low Temp. side |
| | R290 135 g ±5 g | R170 70 g ±3 g n-pentane 10 g (16cc) 0+3 g (12.5%wt) | R290 135 g ±5 g | R170 85 g ±3 g n-pentane 12 g (19cc) 0+3 g (12.5%wt) |
| Refrigerant oil | Ze-NIUSL22SA | | | |
| Battery | For power failure alarm, Lead storage battery, DC6 V 7200 mAh, Auto-recharge | | | |
| Weight | 246 kg | | 278 kg | |
| Accessories | 1 key, 1 scraper, 1 stick for air intake port cleaning | | | |
| Optional components | Temperature recorder: (MTR-85H MTR-G85C), Recorder fixing (MDF-S3085, MTR-85H) Recorder sensor Cover:(MTR-DU700SF) Backup cooling kit:(MDF-UB7);For Liquid CO ₂ | | | |
| | Small inner door (MDF-5ID4, MDF-5ID5) | | Small inner door (MDF-7ID4, MDF-7ID5) | |
| | Drawers (MDF-50R) | | | |
| | Storage rack (MDF-70SC) Inventory rack (IR-224U, IR-316U) Interface board (MTR-L03)*; For LAN Interface board (MTR-480)*; For RS-232C/RS-485 | | | |

* For the data acquisition system MTR-5000 user only.

Control specifications

| Item | MDF-DU502VH | MDF-DU702VH |
|--------------------------------|---|---|
| Temperature controller | Microprocessor control system: LCD touch panel input Settable range : -90 °C~ -50 °C (unit:1 °C) Memorized by Non-volatile memory | |
| Temperature sensor | Platinum resistance (Type: PT 1000 Ω) | |
| Temperature display | WVGA full color LCD Digital display (Unit: 1 °C) | |
| Alarms | High / Low temp. alarm | SV±5 °C~±40 °C, changeable (default:±10 °C) When the chamber temperature exceeds the alarm set temperature, "Alarm" is displayed alternately in normal characters and reverse video in the alarm display field and chamber temperature blinks. After 15 min. (0-15 min. changeable, default: 15°C), "Warning" is displayed in the message display field, buzzer beeps and remote alarm activates and chamber temperature blinks. |
| | Door alarm | When the door opens, "Open" is displayed alternately in normal characters and reverse video in the alarm display field. After 2 min. (0-15 min. changeable), buzzer beeps, no remote alarm activates. |
| | Fan lock alarm | No detection function. |
| | Power failure | "Warning" is displayed alternately in normal characters and reverse video in the alarm display field. And the notification is displayed in the message display field. Buzzer beeps and the remote alarm activates. |
| | Filter alarm | Filter alarm indicator is lit when the excessive dust is accumulated on the condenser filter. When this indicator is lit, clean the condenser filter. |
| | Remote alarm | Remote alarm terminal 3P; contact capacity DC 30 V, 2 A (Max) Remote alarm activates when temp. alarm or power failure occur. |
| Parts replacement notification | Battery age | Two type of battery accumulation time display in service mode. At replacement time, it is informed in the message display field. |
| | Fan motor age | Fan motor accumulation time display in service mode. Replacement notification is not performed. |
| Status | Ambient temperature abnormality | When the sensor is out of the range of 35 °C from 0 °C, the notification display in the message display field. |
| | Supply voltage drop | No detection function. |
| | High-load operation display | When low side compressor operation rate is over 95 %, the notification display in the message display field. |
| Self-diagnosis function | When a sensor is failed, the notification display in the message display field. Buzzer beeps and remote alarm activates. | |
| Compressor protection | Overload relay (internal) Detected value of filter sensor is 50 °C or more. "Warning" is displayed alternately in normal characters and reverse video in the alarm display field and the notification is displayed in the message display field. Buzzer beeps and remote alarm activates. | |
| Key lock | Selection of ON or OFF by slide in Key Lock mode. OFF: Unlocked ON: Locked | |

Performance specifications


| Item | MDF-DU502VH | MDF-DU702VH |
|---------------------------|--|--|
| Cooling performance | -86 °C (AT;30 °C, no load) | |
| Temperature control range | -86 °C~ -50 °C (ambient temperature; 30 °C, no load)* | |
| Rated voltage | AC 220/230/240 V | |
| Rated frequency | 50 Hz (PE) 60 Hz (PA) | |
| Rated power consumption | 430 W (Max. 820 W)(PA) 420 W (Max. 840 W)(PE) | 545 W (MAX 930W/ 945 W/ 960 W) (PE) 550 W (MAX 930 W) (PA) |
| Rated heat radiation | 1548 kJ/h (Max. 2952 kJ/h) 1512 kJ/h (Max. 3024 kJ/h) | 1962 kJ/h (Max.3348 kJ/h /3402 kJ/h /3456 kJ/h) (PE) 1980 kJ/h (Max. 3348 kJ/h))(PA) |
| Noise level | 52 dB [A] (background noise; 20 dB) | |
| Maximum pressure | 2200 kPa | 2200 kPa |

*Maximum cooling performance.


The chamber temperature can be reached at -86 °C at ambient temperature 30 °C with no load.

Dimension

MDF-DU502VH



MDF-DU702VH



Cooling unit parts

MDF-DU502VH

1Φ220-240V 50/60Hz


| Parts description | Specification | |
|-----------------------------|--|---|
| | Circuit H(Left side at rear view) | Circuit L(Right side at rear view) |
| Compressor | | |
| Type | VNEU213U | VNEU213U |
| Compressor cord | LDDG000200 | LDDG000200 |
| Rated power supply | Single phase, 220-240 V, 50/60 Hz | Single phase, 220-240 V, 50/60 Hz |
| Refrigeration oil | Ze-NIUSL22SA Q'ty:500 cc | Ze-NIUSL22SA Q'ty:500 cc |
| Cooling method | Forcible air circulation (partially) | Forcible air circulation (partially) |
| Inverter | | |
| Type | VENEU213U | VENEU213U |
| Inverter code | LDLB041800 | LDLB041800 |
| power supply | 220-240 V, 50/60 Hz | 220-240 V, 50/60 Hz |
| Condenser | | |
| Type | Finless tube | Finless tube |
| Condenser | 12 columns x 2 lines 5pitch Fin 50spc | Coil pipe Φ6.35 mm |
| Pre-condenser | W250 mm (within 1column x6lines pre-con) | |
| Frame pipe | Φ6.35mm | |
| Evaporator | Cascade condenser, Shell & Tube | Tube on sheet, Φ7.94mm (Sharing with interior) |
| Type | Φ80 mm | |
| Accumulator | Φ38 mm | |
| Accumulator | Φ38mm | |
| Capillary tube | | |
| Resistance | 78 PSI | 4.0 kg/cm ² |
| Length | 3000 mm | 2000 mm |
| Outer diameter | Φ2.4 mm | Φ1.8 mm |
| Inner diameter | Φ1.2 mm | Φ0.9 mm |
| Refrigerant | TYPE R-290 | Amount 135g ±5 g |
| | | TYPE R-170 n-pentane |
| | | Amount 70g ±3 g 10g (16cc)-0+3g |
| Dryer | 4AXH-9 18 g | 4AXH-9 55 g |
| Condensing fan | Material: ABS Blade: 4 pcs. 230 mm | Material: ABS Blade: 4 pcs. Φ230 mm |
| Condensing fan motor | SV4-11AB5P (440VAC,1.0μF running condenser) | SV4-11AB5P (440VAC,1.0μF running condenser) |
| Type | | |

MDF-DU702VH
1Φ220-240V 50/60Hz


| Parts description | Specification | | | |
|-----------------------------|--|---|---------------|----------------------------------|
| | Circuit H(Left side at rear view) | Circuit L(Right side at rear view) | | |
| Compressor | | | | |
| Type | VNEU213U | VNEU213U | | |
| Compressor cord | LDDG000200 | LDDG000200 | | |
| Rated power supply | Single phase, 220-240 V, 50/60 Hz | Single phase, 220-240 V, 50/60 Hz | | |
| Refrigeration oil | Ze-NIUSL22SA Q'ty:500 cc | Ze-NIUSL22SA Q'ty:500 cc | | |
| Cooling method | Forcible air circulation(partially) | Forcible air circulation(partially) | | |
| Inverter | | | | |
| Type | VENEU213U | VENEU213U | | |
| Inverter code | LDLB041800 | LDLB041800 | | |
| power supply | 220-240 V, 50/60 Hz | 220-240 V, 50/60 Hz | | |
| Condenser | | | | |
| Type | Finless tube | Finless tube | | |
| Condenser | 12 columns x 2 lines 5pitch Fin 50spc | Coil pipe Φ6.35mm | | |
| Pre-condenser | W250 mm (within 1column x6lines pre-con) | | | |
| Frame pipe | Φ6.35mm | | | |
| Evaporator | Cascade condenser, Shell & Tube Φ80mm | Tube on sheet, Φ7.94mm (Sharing with interior) | | |
| Accumulator | Φ38mm | | | |
| Capillary tube | | | | |
| Resistance | 78PSI | 4.0 kg/cm ² | | |
| Length | 3000mm | 2000mm | | |
| Outer diameter | Φ2.4mm | Φ1.8mm | | |
| Inner diameter | Φ1.2mm | Φ0.9mm | | |
| Refrigerant | TYPE R-290 | Amount 135g ±5 g | TYPE R-170 | Amount 85g ±3 g |
| | | n-pentane | | 12g (19cc)-0+3g |
| Dryer | 4AXH-9 18 g | | 4AXH-9 | 55 g |
| Condensing fan | Material: ABS Blade: 4 pcs. 230mm | Material: ABS Blade: 4 pcs. Φ230mm | | |
| Condensing fan motor | SV4-11AB5P (440VAC,1.0µF running condenser) | | SV4-11AB5P | |
| Type | | | | (440VAC,1.0µF running condenser) |

Refrigeration circuit

Binary refrigerating circuit <cascade system> <MDF-DU502VH/DU702VH>



Components on PCB



Connections on PCB

The following shows connections of connector on Main PCB.

| Connector | Connects to | Usage |
|-----------|--|---|
| CN1 | Emulator | To version up firmware |
| CN2 | #1~#6: MTR-480/MTR-L03 (Option) | To connect with MTR-480/MTR-L03. |
| CN3 | Remote alarm terminal #1: COM #2: N.O. #3: N.C. | Remote alarm contact outputs. |
| CN4 | #1~#2 Jumper | To upgrade firmware version |
| CN5 | #1,#7:USB PCB #8,#9:Door switch | To connect with LCD PCB To detect door open and shut |
| CN6 | #1~#2:Temp. sensor | To control chamber temperature |
| CN7 | #1~#2:Cascade sensor #3~#4:Filter sensor | To detect temperature of cascade To detect temperature of filter |
| CN8 | #1~#2:AT sensor | To detect ambient temperature |
| CN9 | #1~#4:Unused | |
| CN11 | #1,#2:Main battery switch | To supply power during power failure |
| CN12 | #1,#3:Battery switch | To supply power during power failure |
| CN14 | #1~#6:Back-up battery charging PCB | To control battery charging |
| CN15 | #1~#3:Inverter L #3~#6:Inverter H | To control inverter L To control inverter H |
| CN16 | #1~#4:Unused | |
| CN17 | #1,#2: Auto air intake port (AIP) heater | To control AIP heater |
| CN101 | #1,#4:Recorder, Power switch | To supply power to PCB |
| CN103 | #1,#3:Fan motor condenser | To control condensing fan motor |


Electric parts


| MDF-DU502VH MDF-DU702VH | | -PA | -PE |
|---|----------------------------|----------------------------------|----------------------------------|
| Compressor(A)(B) | Type Rating | VNEU213U 220-240 V, 50/60 Hz | VNEU213U 220-240 V, 50/60 Hz |
| Inverter | Type Rating | VCCHP2456 220-240 V, 50/60 Hz | VCCHP2456 220-240 V, 50/60 Hz |
| Condensing fan motor | Type Rating | SV4-11AB5P 230 V,10 W | SV4-11AB5P 230 V,10 W |
| Capi. Tube Heater(A),(B) MDF-DU702VH S/N:~18010001 only | Rating Resistance(25°C) | 230V,12W 4700Ω | 230V,12W 4700Ω |
| Temp. sensoor | Type Rating | Pt (THC-663) 1000Ω | Pt (THC-663) 1000Ω |
| Cascade sensor | Type Rating | 502AT-11 5 kΩ,25 °C | 502AT-11 5 kΩ,25 °C |
| Filter sensor | Type Rating | 502AT-11 5 kΩ,25 °C | 502AT-11 5 kΩ,25 °C |
| Door Switch | Type Rating | SDKNA20700 5 V 5 mA | SDKNA20700 5 V 5 mA |
| Battery | Type Rating | LC-P067R2J 6 V 7.2 AH | LC-P067R2J 6 V 7.2 AH |
| Battery Switch | Type Rating | SLE6A2-5 4 A 250 VAC | SLE6A2-5 4 A 250 VAC |
| Power Switch | Type Rating | JW-L21RRK-P4G 16A, 250VAC | JW-L21RRK-P4G 16A, 250VAC |
| Power Cord | Plug Length | U3 250VAC 15A 3.8m | K3 250VAC 16A 3.5m |

Specifications of sensor


The following shows the temperature in thermal sensor (502AT-11) and its resistance value.

| Temp. (°C) | Resistance Value (kΩ) | Temp. (°C) | Resistance Value (kΩ) | Temp. (°C) | Resistance Value (kΩ) |
|---------------|--------------------------|---------------|--------------------------|---------------|--------------------------|
| -50 | 154.6 | -7 | 17.93 | 12 | 8.167 |
| -45 | 116.5 | -6 | 17.16 | 13 | 7.853 |
| -40 | 88.91 | -5 | 16.43 | 14 | 7.553 |
| -35 | 68.19 | -4 | 15.74 | 15 | 7.267 |
| -30 | 52.87 | -3 | 15.08 | 16 | 6.993 |
| -25 | 41.21 | -2 | 14.46 | 17 | 6.731 |
| -20 | 32.44 | -1 | 13.86 | 18 | 6.481 |
| -19 | 30.93 | 0 | 13.29 | 19 | 6.242 |
| -18 | 29.51 | 1 | 12.74 | 20 | 6.013 |
| -17 | 28.16 | 2 | 12.22 | 25 | 5.000 |
| -16 | 26.88 | 3 | 11.73 | 30 | 4.179 |
| -15 | 25.66 | 4 | 11.25 | 35 | 3.508 |
| -14 | 24.52 | 5 | 10.80 | 40 | 2.961 |
| -13 | 23.43 | 6 | 10.37 | 45 | 2.509 |
| -12 | 22.39 | 7 | 9.960 | 50 | 2.137 |
| -11 | 21.41 | 8 | 9.569 | 55 | 1.826 |
| -10 | 20.48 | 9 | 9.196 | 60 | 1.567 |
| -9 | 19.59 | 10 | 8.840 | | |
| -8 | 18.74 | 11 | 8.496 | | |

The following shows the temperature in thermal control sensor (PT1000Ω) and its resistance value.


| Temp. (°C) | Resistance Value (Ω) | Temp. (°C) | Resistance Value (Ω) | Temp. (°C) | Resistance Value (Ω) |
|---------------|-------------------------|---------------|-------------------------|---------------|-------------------------|
| -140 | 452.8 | -70 | 730.3 | 0 | 1000.0 |
| -130 | 493.0 | -60 | 769.3 | 10 | 1038.0 |
| -120 | 533.1 | -50 | 808.1 | 20 | 1076.0 |
| -110 | 572.9 | -40 | 846.7 | 30 | 1113.8 |
| -100 | 612.6 | -30 | 885.2 | 40 | 1151.4 |
| -90 | 652.0 | -20 | 923.6 | 50 | 1189.0 |
| -80 | 691.3 | -10 | 961.9 | 60 | 1226.4 |


Wiring Diagram




Circuit Diagram

main







Operation


How to operate Service # 1 and Service # 2 on the service menu is described below. For other operations on the LCD, refer to the instruction manual.

| Term | Function |
|-------------------------|--|
| Runtime | Display the accumulate time of Main battery, Backup battery FAN A and FAN B, Reset them. |
| Version UP | Update the firmware of main board and touch panel. |
| Temp. Calibration | Chamber adjustment method. |
| System Configuration #1 | Setting Capi. Heater and its timer, compressors, Fan A/B . |
| Status | Monitoring of chamber temp. , ambient temp, comp. evaporator and fan. |
| System Configuration #2 | Turning On/Off of AIP heater, and Backup. Setting overload alarm delay.. |
| Reset Data | Initialization of data, and flash memory. Setting model code. |
| Error Log | Output of error log data. |
| Setup Data | Input and output the data. |
| Other Setting | Turning On/Off Screen Capture, Demo Mode, Self-diagnosis Function and FAN Alarm. |


CAUTION) Firmware update, compressor, FAN, and cap heater settings affect basic performance.

Be careful when operating.

(Runtime, Version Up)




Keep touching MENU KEY (5 sec)




Move to Menu display, touch service #1 key.

Service code "384"



- Main Battery
Main Battery accumulated time
- Backup Battery
Backup Battery accumulated time
- FAN
FAN accumulated time

When each part replace, reset corresponding data.




Move to Service Code input display, so input "384" or "335232"



Move to Service #1 display, touch Runtime or Version Up key

Service code "335232"

Input "384",
only can check
the firmware
version.



- Confirm the existence of file in predefined folder of USB memory.
- If file is not existence, message "File not found." is displayed in "Notice" screen and the process end.
- If file is a existence, if it is not compatible model, "Confirm firmware" is displayed and the process is terminated. If it is a compatible model, "Confirm firmware rewrite ..." is displayed and when "Yes" is selected, version upgrade of the touch panel embedded software is started, and when "No" is selected, the process is ended.
- During version up, wheel appear in the screen.
- When version upgrade is completed, "Firmware rewrite complete. ..." is displayed.
- If an error (other than Verify check) occurs during version upgrade, "Firmware rewrite failure ..." is displayed
- If a Verify error occurs in the Verify check, "Verify check NG. ..." is displayed.

※When some version up files are plurality, the first found file is used for version up
※During appearance "wheel", all operations are not able to receive.


※"OK" button in "Result" screen is disable, and all operations are not able to receive. (unless user turn off power)

※"OK" button in "Error" screen is disable, and all operations are not able to received. (unless user turn off power)


※During version up, no communication is not carried out and no log.is not output.

※If you cancel after or during upgrade, restart the main power supply.

(Temperature calibration procedure : chamber sensor)



Keep touching MENU KEY (5 sec)




Move to Service Code input display,
so input "335232"




Move to Menu display, touch Service #2 key.




Move to Service display, touch Temp.
Calibration key




Touch mark area.



Move to input value display.
Input measured value by thermometer



Cal key is vialed within setting range, touch
Cal key.



Touch TOP key, so return to Menu display.
Input value area return to zero.

NOTE:: Ambient Temp. sensor and condenser A&B sensor can calibrate in the same procedure.


(Operation System Configuration #1)



Move to Menu display, touch Service #2 key.



Move to Service display, touch System Configuration #1 key



- Capi. Heater Timer
Slide key ON(enable),
setting time(Default 8 min).
- Capi. Heater Start
Capi. heater control setting
(0) Heater on once in 18 hours
(1) Force on now
(2) Force off now(default)
- Comp. Control
Compressor control setting
(0)Off
(1)On
(2)Auto(default)
(3)Manual
- FAN control setting
(0) Off
(1) On
(2) Auto(default)

Refer to Control Specification for Capi. Heater Start,
Comp. Control, FAN Control.

(Operation Configuration #2, Status)

The diagram illustrates the navigation flow through the system configuration menus:

- Service #2** screen (Top Left):
 - Buttons: Temp. Calibration, System Configuration #1, **System Configuration #2**, Status, Reset Data, Setup Data, Error Log, Other Setting.
 - Red circles highlight "System Configuration #2" and "Status".
- System Configuration #2** screen (Top Right):
 - Buttons: **AIP Heater ON/OFF : ON**, Backup Function : OFF, Overload Alarm Delay : No. 0, Apply.
- Status** screen (Bottom Left):
 - Text: Compressor Control : ON
 - Table:

| Temp. SV | Temp. PV | Ambient Temp. | Running rate |
|----------|----------|---------------|--------------|
| -000.0°C | -000.0°C | -000.0°C | 000% |

| Comp. | PV | ON | OFF | Filter Temp. | FAN |
|-------|------|--------|--------|---------------|-------------|
| | | | | Cascade Temp. | Capi.Heater |
| H | 000% | 000min | 000min | -000.0°C | ON |
| L | 000% | 000min | 000min | -000.0°C | OFF |
- Note:** Use for validation

Move to Status display

Status screen (Bottom Left):


- Buttons: Page Select, Top, Back.
- Table:

| Eva. In Temp. | Eva. Out Temp. | Ref. Temp. | Cal. Voltage |
|---------------|----------------|------------|--------------|
| -000.0°C | -000.0°C | -000.0°C | 0000mV |


| Comp. | Speed | FAN (AD) | Main Battery | Air Intake Port |
|-------|-------|----------|--------------|-----------------|
| | | | Sub Battery | |
| H | 0000 | -0000 | 00.0v | ON |
| L | 0000 | -0000 | 00.0v | |

Move to next page.


(Operation Reset Data , Error Log)



Move to Reset Data display




Move to Error Log display



- Reset data
On(enable), OFF(disable)
- Function Code
925: Flash memory Initialization
- Model Code
Set model
(1) MDF-DU502VH
(2) MDF-DU702VH
Use for the version up of the firmware

- Check if USB memory is connected
- If not connected, "USB memory is disconnected." is displayed
- Communicate with ErrorLog command read request (0xA0)
(Suspends communication with Data command read request (0xFF))
- If communication fails, "Communication Failure." is displayed.
- If communication is successful, output error log data (2,048 bytes) to USB memory in binary format (file name: ErrorLog.bin)
- During communication and output, a wheel appears in the screen.
- It outputs log information of all seven periods including service (service mode #2) in the CSV format to the USB memory.
 - ① Actual Temp. (°C)
 - ② Filter Temp. (°C)
 - ③ Cascade Temp. (°C)
 - ④ Ambient Temp. (°C)
 - ⑤ Door Opening
- Output alarm information of all periods to the USB memory in CSV format.
- If an error occurs during output, display a warning message depending on the situation and interrupt processing
- When output is completed, "Export complete." is displayed in Result box.
- ※ If settings storage data already exists in USB memory, overwrite it.

(Operation Setup Data , Other Setting)



| ■ Other Setting | |
|--------------------------------------|--------------------------------------|
| | |
| Screen Capture | : <input type="button" value="ON"/> |
| Demo Mode | : <input type="button" value="OFF"/> |
| Self-diagnosis Function | : <input type="button" value="ON"/> |
| FAN Alarm | : <input type="button" value="ON"/> |
| <input type="button" value="Apply"/> | |

- Screen Capture
- Demo Mode
- Self-diagnosis Function
- FAN Alarm

USB memory data can be import and export
Import data: setting data
Export data: daily information

■ Data Export

- Confirm the connection of USB memory
- If it is not connected, message; "USB memory is disconnected" is displayed in "Notice screen".
- The communication carry out by Setup Data command read request(0x90). (By Data command read request(0xFF), the communication is temporary halt.)
- If communication is failure, "Communication Failure" is displayed .
- If communication is success, the saved setting data + LCD control value output to USB memory. by BIN format
- During communication and output, wheel appear in the screen.
- If error occurs during output, depending on any situations, warning message is displayed and the process break.
- When the output complete, "Export complete" is displayed in "Result" box.
- ※ If the previous saved setting data exit in USB memory, the data is overwritten.

■ Data Import

- Confirm the existence of file in predefined folder of the USB memory.
- If file is not existence or CRC does not match, "File not found." is displayed in "Notice" screen and the process end.
- If file is existence, "Yes" is displayed in confirm screen A.
And if "Yes" is selected, the communication is carried out by SetupData command wright request (0x91) and saved setting date.is sent.
The setting data is upgraded to saved setting value of LCD control
- During communication, wheel appears in the screen.
If the communication is successful, "Complete" is displayed in Result box.
- If the communication is failure, "Error" screen is displayed

Control specification

1. Temperature control

(1) Outline

This unit controls two comp.s on the H side and the L side in cooperation, and controls the temp. inside the refrigerator to the set temp. (two-way freezing).

Furthermore, since the two comp.s are of inverter type and the number of revolutions can be changed, PID control is used in the vicinity of the set temp. inside the compartment.

(2) Setting(user)

| Item | Set range | Default |
|------------------|-----------------|---------|
| Chamber temp. | -90 °C ~ -50 °C | -80 °C |
| Compressor delay | 3 min ~ 15 min | 3 min |

(3) Restrict of compressor starting

- Immediately after turning on the power supply or immediately after recovery from power failure, do not start the compressor for the compressor delay time.
- When stopping the compressor, do not start for 3 minutes.

(4) Compressor mode

| Name | Action | Condition |
|--------|--|---|
| NORMAL | Temperature control by NORMAL MODE | During normal operation |
| ON | Basically it is the same as NORMAL model, and both comp.s are always ON without temp. control. The rotational speed is maintained constant speed operation | During security operation (chamber temp. unknown) |
| OFF | Both comp.s are always OFF | Demo mode, Power failure Security operation (comp. protection temp. alarm) |

(5) Condition and Priority of comp. operation mode

- The condition and the priority of comp. operation mode as following bellow table.
When the state overlap, the comp. operation mode is determined by the priority.
For each state,
- For the power failure and demo mode, both are same comp. running operation mode. So, they are coedited as same level.
- Test Program and capirarily heater control are not run within being higher priority condition.
If higher priority condition starts, these are stopped.

| Priority | Condition |
|----------|---|
| High | Power failure, Demo mode |
| | Security operation (There is also priority in this) |
| | Test Program |
| | Capi heater control |
| Low | Other than those above |

Comp.Control in service menu

2. Compressor normal mode


(1) (Speed control)

Three type of comp. operation mode depending the chamber temp.

| Chamber temperature | > set temp.(SV) + Proportional band lower limit (EL) | \leq set temp.(SV) + Proportional band lower limit (EL) \geq set temp.(SV)- Proportional band upper limit (EH) | < set temp.(SV) - Proportional band upper limit (EH) |
|---------------------|---|---|---|
| Speed control | Constant speed control | PID control | stop |

(Constant speed operation)

The control method in the case of starting the comp. from the state where both the comp.s are stopped is shown below



① H side start condition (AND)

- After the H side stopped, the threshold (3 min) elapsed
- (Immediately after H side start failure, 10 sec elapsed)
- If the chamber temp. is higher than the set temp. plus the threshold (0.5°C)

② L side start condition (AND)

- After the L side stopped, the threshold (4 minutes) elapsed
- (Immediately after L side start failure, 10 min elapsed)
- A time longer than the threshold(※1) has elapsed after the H side start.

threshold (※1) condition

short(60 sec) Cascade temp. is below threshold (-20°C) and within 2 hours from the L side stop
long (480 sec) Other than those above

<After security operation>

- After the L side stop, the threshold(4 min) elapsed.
- After the H side start, 1min elapsed
- The cascade temp. is less than -34°C or 15 min elapse from the above

③ H side stop condition (OR)

- the chamber temperture is lower than set temp.-threshold (0.5°C).
- (security operation)

④ L side stop condition(OR)

- H side stop
- Outside the initial startup period and the cascade temp. is higher than the threshold (-10°C)

⑤ Condition for L side restart (AND)

- After the L side stopped, the threshold (4 min) elapsed
- (Immediately after H side start failure, 10 seconds elapsed)
- The cascade temperature is -34°C or less, or even 15 minutes elapse

Number of revolutions at L side restart

First at 2500 rpm,
 Return to 3000 rpm with cascade temperature below -20 °C
 ⑥ Transition conditions from the initial startup period to the constant-speed initial period (OR)
 • Initial startup period has passed 120 sec or more
 • Other than the conditions below
 Extended for 30 sec when the cascade temp. is higher than the threshold (-10 °C) and the cascade temp. drop detection.
 During the startup period, if the ambient temp. is 32 °C or more, increase the H rotation to 4000 rpm
 ⑦ Transition conditions to constant speed initial period or constant speed period
 • Detect maximum peak value of the cascade temp. and it is -30 °C or less.
【Error handling】
 If the following conditions occur, reset the microcomputer once
 • Command failure occurred 20 times consecutively
 If the above condition occurs even after resetting the microcomputer, it outputs a communication failure warning.
 If the communication failure is canceled, the alarm is canceled.

【Addition】

If the following condition (AND) under the constant speed period and the forced-ON mode, change the L side rotational speed from 3000 to 3150.
 • Set temp. is -86 °C or less
 • Chamber temp. is -83 °C or less (return to 3000 rpm at -82 °C over)
 If the following conditions, change the rotation speed during PID control from 3000 to 3150.
 • Set temp. is -83 °C or less

3. Fan control

(1) Control

The fan is operated in the following mode against the compressor operation.
 table FAN operation mode

| Name | Operation | Condition |
|-------|---|---|
| A_ON | Fan is always On | When On is set by FAN control※1 |
| A_OFF | Fan is always Off | When Off is set by FAN control※1 |
| AUTO | Change control according to ambient temp. | Normal When Auto is set by FAN control |

※1 Temporary setting in manufacturing process and service work.

Then, the state transition of AUTO mode is shown in the table below.

- FAN_OFF : 31°C (Ambient temp. measurement value)
- FAN_ON : 32°C (Ambient temp. measurement value)

- ① AUTO mode
 FAN_ON \leq T both A and B, ON
 FAN_OFF < T < FAN_ON refer to ※2 T \leq FAN_OFF Link with comp.

Table AUTO mode

| State | | Chamber temp.(T) | | |
|-------|----------------|----------------------|----------------------|-----------------|
| No. | FAN | T \leq FAN_OFF | FAN_OFF < T < FAN_ON | FAN_ON \leq T |
| C11 | Stop | Link with H comp. ※2 | ※2 | ON C12 |
| C12 | Run | Link with H comp. ※2 | ※2 | (none) — |
| | Control policy | Link with H comp. ※2 | none | Always ON |

※2 If it does not match either state of link with comp. or AB_ON, it is compatible with comp.

4. Capi. heater control

(1) Outline

In order to prevent oil clogging in the capillary, the compressor is periodically stopped, and the heater attached to the capillary is energized. Even if the compressor is in operation, the compressor is forcibly stopped for a certain period depending on the condition.

(2) Set value(service work)

| Item | Range | Default |
|-------------------|----------------------------|---------|
| Capi. Heater time | 0 min, 6~15 min | 8 min |
| Capi heater mode | 0: AUTO / 1 : ON / 2 : OFF | 2 : OFF |

(3) Control

If AUTO is set, the following operation is performed after 18 hours of product operation start (every 18 hours thereafter). If OFF is set, the same operation as AUTO is performed but only the energization of the capi. heater is canceled.

- ① Energize the capillary every 18 hours and do not start the compressor during this time.
- ② If the internal temperature has not reached the OFF1 point, set the compressor to the forced ON state and wait until OFF1 or less.
- ③ If it does not reach OFF1 or less even after 60 min elapse, forcibly stop the compressor and energize the capillary.

When On is set, the forcible compressor is stopped and the capillary is energized. After that, it switches to AUTO.




Fig. 4 Capi. energization timing

5. AIP heater control

This heater is always turned ON, but for energy saving or protection of abnormal heating, is turned OFF under the following conditions

Heater OFF condition

- When the chamber temp is more than -39°C, turn off. (when it is less than -40°C, turn ON)
- Opening the door
- On power failure or demo mode

6. Buzzer control

Refer to the following to the types and specifications or buzzer sound and the beep conditions.

Buzzer sound spec.

| Spec.(On ms/Off ms/On time) | Type | Condition |
|-----------------------------|--------------------|--|
| (120ms/120ms/twice) | Starting sound | Starting sound(normal) |
| (780ms/220ms/continue) | Alarm sound(long) | Alarming |
| (60ms/1600ms/continue) | Alarm sound(short) | Power failure, alarm history, ALARM test |

7. Alarm type

Table Alarm type

| Code | Name | Class | Judgment condition |
|------------|---------------------------------------|---------------------|---|
| W04 A04 | High temp. alarm | Temp. alarm | the chamber temp. exceeds the warning temp. |
| W05 A05 | Low temp. alarm | Temp. alarm | The chamber temp. falls below the warning temp. |
| W01 | Power failure alarm | Power failure alarm | During detecting power failure |
| - | Door alarm | Door alarm | During opening door |
| - | Filter alarm | Alarm Notice | Filter sensor temp. is 45 °C or more Release condition: Filter sensor temp. is 40 °C or less, or condenser temp. abnormality occurred |
| W02 | Condenser temp. abnormal | Warning alarm | Filter sensor detected value is 50°C or more Release condition: filter sensor is AT + 5°C or less |
| W09 | Chamber temp. control sensor open | Warning alarm | After 5 V voltage conversion, 4950 mV or more |
| W10 | Chamber temp. control sensor short | Warning alarm | After 5 V voltage conversion 25 mV or less |
| S10 | Cascade sensor open | Status alarm | After 5 V voltage conversion, 4950 mV or more |
| S11 | Cascade sensor short | Status alarm | After 5 V voltage conversion, 25 mV or less |
| S12 | Filter sensor open | Status alarm | After 5 V voltage conversion, 4950 mV or more |
| S13 | Filter sensor short | Status alarm | After 5 V voltage conversion, 25 mV or less |
| S14 | AT sensor open | Status alarm | After 5 V voltage conversion, 4950 mV or more |
| S15 | AT sensor short | Status alarm | After 5 V voltage conversion, 25 mV or less |
| S20 | Battery switch | Status alarm | Detected 7.2 V or more |
| W06 | H side inverter communication failure | Warning alarm | Communication failure (more than 20 consecutive occurrences, more than 20 consecutive occurrences after reset), timeout (more than 10 consecutive times), circuit board short |
| W07 | L side inverter communication failure | Warning alarm | Communication failure (more than 20 consecutive occurrences, more than 20 consecutive occurrences after reset), timeout (more than 10 consecutive times), circuit board short |
| - | Communication abnormal | - | Detected on the touch panel side |
| S02 | Over range of ambient AT | Status notice | The acquired value is 0 °C or less, 35 °C or more |
| S01 | Overload running notice | Status notice | It does not reach the set temp. for more than 5 days (120 hours) |
| S18 | Main battery change notice | Status Notice | The number of working days of the product has passed 1022 days (about 2.8 years) |
| S19 | Backup battery change notice | Status Notice | 1022 days (about 2.8 years) passed while the Backup function was on |
| S03 | AIP port heater abnormal | Status Notice | Heater control state and acquisition voltage do not match at ON control, less than 0.2 V At OFF control, error judgment is made at 0.1 V or more |
| S16 | Battery re-charge error | Status Notice | Voltage of the main battery is 6 V or less and there is no voltage rise of 0.3 V or more even after 6 hours |
| S17 | Backup battery re-charge error | Status Notice | Voltage of the main battery is 6 V or less and there is no voltage rise of 0.3 V or more even after 6 hours |
| S21 | Main microcomputer flash error | - | Detected on the touch panel side |
| S22 | Communication protocol error | - | Detected on the touch panel side |

8. Alarm function

(1) Alarm delay

- Even if the judgment is alarm, buzzer sound and remote alarm is not output within the set time.
(after passing set time, the status become alarm.)
- Set range of alarm delay time:※ restrict input by touch panel side.
Temp. alarm : 0 min (Immediately: no alarm delay) , 1min~15min (default: 15min)
Door alarm 0min (Immediately: no alarm delay) , 1min~15min (default: 2min)

(2) Alarm buzzer output and erase

- The alarm buzzer output on alarm status, and if push "Buzzer" key during output of alarm buzzer, mute the buzzer sound.
- Set range of alarm delay time:※ restrict input by touch panel side.
Temp. alarm : 0 min (Immediately: no alarm delay) , 1min~15min (default: 15min)
Door alarm 0min (Immediately: no alarm delay) , 1min~15min (default: 2min)

(3) Ring back

- After passing set time from alarm buzzer mute, buzzer sound again.
- Set range of ring back, ※restrict input by touch panel side.
Ring back function OFF , 1min~99min (default: 30min)

(4) Remote alarm

- This is linkage with alarm buzzer output. And restrict the output from remote alarm port.
(if the alarm buzzer output, the remote alarm output too.)
However, when the buzzer alarm is no-linkage, even if the alarm is mute status, the remote alarm keep output.
- Remote alarm set
Buzzer no-linkage(ON)、buzzer linkage(OFF) (default:buzzer no-linkage)

(5) Alarm history

If the judge of alarm from alarm status is released, the status transitions alarm history, and that buzzer has outputted and alarm has occurred are displayed.

(6) Alarm history release

If put the buzzer button within alarm history state, buzzer erase and move to normal state.

Table function vs alarm

| Function Type | (1)Alarm delay | (2)Alarm buzzer output / erase (4)Remote alarm | (3)Ring back | (5)Alarm history (6)History release |
|---------------------|----------------|---|--------------|--|
| Temp. alarm | ○ | ○ | ○ | ○ |
| Power failure alarm | — | ○ | ○ | ○ |
| Door alarm | ○ | ○(Remote alarm —) | — | — |
| Warning | — | ○ | ○ | — |
| Alarm Status | — | ○(Remote alarm —) | — | — |
| Notice | — | — | — | — |

10. Model

- (1) Model code
 1: MDF-DU502VH
 2: MDF-DU702VH

- (2) Parameter

| Module | Term | Default | | Range |
|---------------------|------------------------------|---------------|---------------|---------------|
| | | MDF-DU502VH | MDF-DU702VH | |
| Comp. control | Set temp. | -80°C | -80°C | -50°C ~ -90°C |
| | Control mode | Auto | Auto | fix |
| | Comp. first start delay time | 3 min | 3 min | 3 min~15 min |
| | Comp. restart delay time | 3 min | 3 min | fix |
| | ON point | +0.5°C | +0.5°C | fix |
| | OFF point | -0.5°C | -0.5°C | fix |
| | Offset value | +1.0°C | +1.0°C | fix |
| | Ambient temp.correct | -1.0°C | -1.0°C | fix |
| | ON/OFF adjustment value | existance | existance | fix |
| FAN control | Control mode | AUTO | AUTO | Control modes |
| | ON point | 27°C | 27°C | -20°C~40°C |
| | OFF point | 26°C | 26°C | fix |
| Capi heater control | Control ON/OFF | ON | ON | ON, OFF |
| | Control mode | OFF | OFF | ON, OFF, AUTO |
| | Maximum cooling time | 60 min | 60 min | fix |
| | Action interval | Every 18 hour | Every 18 hour | fix |
| AIP heater control | Control ON/OFF | ON | ON | ON, OFF |
| | ON point | -40°C | -40°C | -50°C~10°C |
| | OFF point | ON point+1°C | ON point+1°C | fix |

※1 Configurable on service menu

| Module | Item | Default | | Range |
|--------------------|--|-------------------------------|-------------------------------|-------------------|
| | | MDF-DU502VH | MDF-DU702VH | |
| Alarm | High temp. alarm temp. | +10°C | +10°C | +5°C ~ +40°C |
| | Low temp. alarm temp. | -10°C | -10°C | -5°C ~ -40°C |
| | Temp. alarm delay time | 15 min | 15 min | 0 min~15 min |
| | Door alarm delay time | 2 min | 2 min | 0 min ~15min |
| | Ring back function | ON | ON | ON, OFF |
| | Ring back time | 30 min | 30 min | 1 min ~99min |
| | Remote alarm linkage | OFF | OFF | ON, OFF |
| | Filter alarm temp. | 45°C | 45°C | ※1 |
| | Filter alarm release temp. | Alarm temp. -5°C or less | Alarm temp. -5°C or less | fix |
| | Comp protection alarm temp. | 50°C | 50°C | ※1 |
| | Comp protection release alarm temp. | Ambient temp. +5°C or less | Ambient temp. +5°C or less | fix |
| | Sensors open alarm | 4950mV | 4950mV | fix |
| | Sensors short alarm | 25mV | 25mV | fix |
| | Battery SW judgment | 100mV | 100mV | fix |
| | Port heater voltage judgement | 0.2V | 0.2V | fix |
| | Range over of ambient temp. | 0~35°C | 0~35°C | fix |
| | Overload running notice | 5 day | 5 day | 6hour,2day,5day※1 |
| | Battery change notice | 1022 day | 1022 day | fix |
| | Touch panel communication error | 130 sec | 130 sec | fix |
| | Backup battery changing failure | 0.3V / 6 hour | 0.3V / 6 hour | fix |
| Calibration value | Chamber temp. sensor PCB calibration value | 0°C | 0°C | ※1 |
| | Chamber temp. sensor calibration value | 0°C | 0°C | ※1 |
| | Ambient temp. sensor calibration | 0°C | 0°C | ※1 |
| | Cascade temp. sensor calibration value | 0°C | 0°C | ※1 |
| | Filter sensor calibration value | 0°C | 0°C | ※1 |
| Security operation | Touch panel OFF | 4.45V | 4.45V | fix |

※1 Variable in service menu



Inverter signal

(1) Outline

In this model, an inverter is used for compressor control. In order to investigate the status of the compressor and the inverter, signal verification of the inverter side is necessary.

(2) The error signal on the inverter side is as follows.

| Failure Mode | LED operation |
|--------------------|-------------------|
| Normal operation | 1 blink each 15 s |
| Thermostat failure | 2 blinks each 5 s |
| Inverter failure | 3 blinks each 5 s |
| Compressor failure | 4 blinks each 5 s |

Parts layout


【Front/Control panel】



【Opposite side of control panel】 【Opposite side of front door /auto air intake port】



【Left side】



Door latch support

Handle

Air intake port


Inner door latch support

Handle

Inner door latch

Handle

【Lower Right side】



Battery switch

Power switch

Remote alarm terminal


【Lower front】

【Chamber】



PT sensor

【Rear side】



Air intake vent(Grill)

Space for Temperature recorder
(Option)

【Lower right side/Electric Box】




Main PCB

【Lower front/Condenser】



【Rear bottom/Cooling UNIT】



Inverter L

Inverter H

Compressor L Compressor H

Repairing unit / Gas charge

■Preparations (H/L side common)

- Refrigerant cylinder (type, amount required to be confirmed), refrigerant cylinder hose
- Alarm
- Detector (gas leak tester for HC refrigerant)
- Dedicated gauge manifold (dedicated GM) (combined with R290 · R170)
- General purpose gauge manifold (GM)
- Long hose (10 m or more: length to go outdoors) Field check required in advance
- Charge hose (4 or more), branching instrument
- Tools (Driver, pliers, pliers, etc.)
- Replacement parts (depending on repair contents)
- Piping
- Explosion-proof vacuum pump
- Powerful vacuum pump
- Piercing tool (2 to 4)
- Service valve
- Service port (check valve) * (See work contents in 2.1)
- Valve driver
- Scale (electronic weighing scale)
- Coupler
- Pinch pliers (2)
- Nitrogen cylinder, pressure regulator
- Leak Check Spray (Liquid detergent may be used)
- Gas welding machine
- Silver braze, flux, exclusive lighter
- Plumbing tools (cutter, reamer, bender, flare tool, expander, wrench, mirror, light etc.)
- Welding cooling tool (water container, waste etc.)
- Welding protective equipment (leather gloves, goggles, flameproof sheets, steel plates, etc.)
- Refrigerant protection (leather gloves, goggles)

■Low side only

- Teflon hose (for pentane)
- Pentane Charging Cylinder
- N-Pentane (Containing amount required)
- Pentane dryer

■Caution

- Use of fire is strictly prohibited during HC gas release and its charge. Especially, pay shut attention for product around and the release port.
- Confirm in advance about the power supply of the work place, weld work availability, carry-in / out route.
- Due to constraints of work environment, welding work may be carried out by moving the product outdoors. In that case, confirm the securing of the power supply at the work place.
- Be careful of burns due to touching the heated part during welding operation.
- Since necessary preparations change depending on the work place and repair content, please check in advance.
- Do not insert the power plug easily without confirming to customer such as emergency outlet.
- Do not use a hose that is cracked or is hard for connector connection, as it may cause gas leakage.

■Common work for high temperature side and low temperature side


1. Release gas

1.1 Securing work space

Preparations : Alarm (Note: **Alarm should be operated at all times during work**)

Procedure: Move to the space where you can work the product, open the window / door and set the alarm on the floor within 50 cm radius from the unit.


Note: Alarms are generally for one type of gas and there is no mode switching function according to multiple gases. Use a detector for detailed measurements such as exploration of gas leakage location.



1.2 Tool connection

Preparations : Detector, dedicated GM, piercing tool, long hose, charge hose, branching instruments, tools, leak check spray (liquid detergent can also be substituted)

Procedure : Connect piercing tool, charge hose, dedicated GM, long hose, as shown below, and lead the long hose tip outdoors.



Note: **Make sure connection is not leaky.**

Although connection to Dedicated GM can be done by hand tightening, confirm that there is no leakage from the connection part when starting gas flow. Make sure there is no break in the middle of the hose (especially long hose).

1.3 Refrigerant release


Preparation: Refrigerant protection equipment (leather gloves, goggles), tools

Procedure: ~~※~~Wearing leather gloves and goggles because it handles refrigerant.

~~※~~Confirm whether there are ignitable ones like electric appliances outdoors.

Install a piercing tool in the high / low pressure service pipe of the product, open the valve of the dedicated GM (however, VR is not opened), and release the refrigerant from the long hose to the outdoor.

(Valve open: VC⇒VG⇒VL⇒VH)



Note: Piercing tool is unnecessary if service valve is attached.

Place the detector near the tip of the outdoor long hose. Be sure to use a detector that can detect HC refrigerant.


Make sure that the release refrigerant is almost gone by the guidelines and detectors of the dedicated GM.

1.4 Vacuum evacuation of cooling circuits

preparation : Explosion-proof vacuum pump, detector for HC refrigerant

procedure: Since the refrigerant is dissolved in the compressor oil, connect the explosion-proof vacuum pump to the tip of the long hose led outdoors, and perform vacuum evacuation from the refrigeration circuit for at least 10 minutes. (Place vacuum pump outdoors)

When no gas is detected, shut all the valves and stop the pump.




Note: Apply a detector to the exhaust port of the vacuum pump to confirm that the refrigerant has been released.

1.5 Vacuum evacuation before nitrogen blowing

Preparation : Nitrogen cylinder, pressure regulator, charge hose, branching instruments, tools, leak check spray (liquid detergent can also be substituted)

Procedure: Connect the nitrogen cylinder as shown below and operate the explosion-proof vacuum pump. Open the valves VC, VG, VR and perform vacuum evacuation between the nitrogen cylinder and the dedicated GM.(About 1 minute)



Note: make sure connection is not leaky.

1.6 Nitrogen charge before nitrogen blowdown

procedure : After 1.5 , VC, VG are shut, VR, VH, VL are opened and nitrogen is charged (about 30 seconds) (however, the secondary pressure is 0.5-0.7 MPa)


note: In order to prevent air from entering the freezing circuit, fill with nitrogen.

| | Valve state | |
|----|-------------|------|
| | open | shut |
| VC | | ○ |
| VG | | ○ |
| VL | ○ | |
| VH | ○ | |
| VR | ○ | |

1.7 Nitrogen blow

Procedure: Once all the valves are shut, reconnect as shown below and blow nitrogen for about 30 seconds into the circuit. (Nitrogen replacement: provided that the secondary pressure is 0.5 to 0.7 MPa)

Caution! Since the refrigerant may burn with the flame of the welding machine, always carry out the above procedure.



Note: Nitrogen blow is performed to completely release the refrigerant. The flow rate of nitrogen may be small (to the extent that it feels a little with the skin).

2. Repair

2.1 Replacement parts

Preparation: Replacement parts, tools, service valves, check valves, piping, piping tools
(Cutter, reamer, bender, flare tool, expander, wrench, mirror etc.)

If you do not use a check valve, you need a coupler, piercing tool, pinch off tool, etc.

Procedure: Exchange necessary parts according to repair contents.

In this system, install a service valve on the low pressure side and a check valve on the high pressure side.

When installing the service valve, replace the low-pressure side service pipe. After that, install the service valve on the low pressure side service pipe.

(Instead of installing a check valve, a coupler or a piercing tool may be used, after vacuum evacuation, a pinch-off tool and welding may be used for sealing)

Also, when repairing the refrigeration circuit, be sure to replace the dehydrator.


※ When replacing parts, follow each parts replacement procedure.

※ Place the removed refrigeration circuit parts in a place away from the work place without fire.



In the PA model, the service pipe label (red tag) is attached to the pipe for service valve installation.

Remove the label of the pipe to be welded.



NOTE : Check valve is also called "access connector".

When the check valve is welded, it will be hot, so remove the valve before doing.

If the check valve cannot be installed on the high pressure side, in the process of "13.

Vacuum exhaust ①"

After vacuuming for 2 to 3 hours, once the high pressure side is pinched and welded,

Vacuum is carried out for 1 to 2 hours from the low pressure side.

2.2 Preparation for gas welding


Preparation: Flameproof sheet, tools, nitrogen cylinder, charge hose, branching instrument

Procedure : Pull the welded part as far as possible from the product and cover the periphery with a flameproof sheet as necessary.

In order to prevent oxide film of pipe inside during welding, a nitrogen cylinder is connected to the low pressure side service valve (SVL)

Replace the nitrogen in the piping. (Nitrogen substitution: provided secondary pressure is 0.5 to 0.7 MPa) At this time,

Remove the valve in the check valve. (Outlet of nitrogen)



Note: Insulation materials are easy to ignite, be sure to cover with fireproof sheets.

The flow rate of nitrogen may be small (to the extent that it feels a little with the skin).

2.3 Piping connection by gas welding

Preparation: gas welding machine, exclusive lighter, silver braze, flux piping tool

(Cutter, reamer, bender, flare tool, expander, wrench, mirror, light etc.)

Welding cooling product (water-filled container, waste etc.),

Weld protective product (leather gloves, goggles, flameproof sheets, steel plates, etc.)

Procedure: Take care not melt the pipe, heat it so that a braze melts around the connection, and then braze to it. Check the brazed state with a mirror, light, etc.

When using silver wax, apply flux to the place where you usually braze.

Note: Forced cooling immediately after brazing heating adversely affects the structure of copper metal. Please leave it for about 2 minutes and then apply wet cloth.

2.4 Refrigeration circuit leak check

Preparations: Leak check spray (liquid detergent can also be substituted)

Procedure: Continue with nitrogen replacement as it is, install a valve in the check valve, and perform a leak check of the entire unit circuit.

(Nitrogen replacement: provided that the secondary pressure is 0.5 to 0.7 MPa)

Note: Use a leak check spray to check the gas leak at the replacement (welding) point.

3. Evacuation

3.1 Leak check of connection point

Preparation: Leak check spray (liquid detergent may be used), dedicated GM, general purpose GM *, powerful vacuum pump, Nitrogen cylinder, charge hose, branching equipment, tools


* In case of evacuating simultaneously on the low temperature side and the high temperature side

Procedure: Connect as shown below and prepare for vacuum evacuation. When repairing on the low temperature side was also done, add general purpose GM and connect. ← Vacuum exhaust at the same time at high temperature and low temperature for shortening the time

(In the case of only the high temperature side, connect the hose from the nitrogen cylinder directly to the refrigerant port)

Open the unopened valves of nitrogen cylinders, VR, VL, VH and general purpose GM,

Perform a leak check of the connection point. (Nitrogen replacement: provided that the secondary pressure is 0.5 to 0.7 MPa)



Note: Use a leak check spray to check the gas leak at the connection point.

Generally, the explosion-proof vacuum pump has a small capacity, so prepare a separate powerful vacuum pump.

(All flammable HC refrigerants are released, so there is no problem even if they are not explosion-proof type)

Connection should be made surely without leakage.

3.2 Evacuation①

Procedure: Remove the hose of the general purpose GM (or coolant port) that was connected from the nitrogen cylinder, open the nitrogen, shut the valve (or VR) of the removed hose, and perform vacuum evacuation to fill the refrigerant. After operating the powerful vacuum pump, open the valve of VC. Vacuuming time: about 4 hours or more with a vacuum pump with a capacity of 300 L / min (When using a vacuum pump with a small capacity value, vacuuming is performed until the degree of vacuum is 3 Torr or less as a guide)

(For high temperature side only, VR shut →)

| | Valve state | |
|----|-------------|------|
| | open | shut |
| VC | ○ | |
| VG | ○ | |
| VL | ○ | |
| VH | ○ | |
| VR | ○ | |

3.3 Complete evacuation

Procedure : Specified time (about 4 hours or more) Shut each valve and stopper temporarily after evacuation and stop the pump.

If vacuum evacuation on the low temperature side is carried out together, disconnect from the refrigerant port.

4. Gas charge

4.1 Evacuation②

Preparations: Dedicated GM, Refrigerant Cylinder, Bottle Cap, Powerful Vacuum Pump, Scale (Electronic Weighing Scale), Charge hose, branching tool, tool


Procedure : After operating the powerful vacuum pump with the refrigerant cylinder on the scale,

Open the valve in the following order and evacuate to the refrigerant cylinder.

VS⇒VR⇒VG⇒VC.

Vacuuming time: about 1 minute with a vacuum pump with capacity of 300 L / min

(Standard: vacuum degree 3 Torr or less)



4.2 Completion of evacuation

Procedure : Specified time (about 1 minute) after evacuating,

shut the valve in the following order and stop the pump.

VC⇒VG⇒VR⇒VS

Remove the vacuum pump from the hose.

Valve state

| | open | shut |
|------|------|------|
| VC | | ○ |
| VG | | ○ |
| VL | | ○ |
| VH | | ○ |
| VR | | ○ |
| VS | | ○ |
| main | | ○ |
| SVL | | ○ |

4.3 Preparation for charge

Preparation: Leak check spray

Procedure: Re-check the type of refrigerant and the amount sealed.

Open the valve slowly in the order of the main valve of the refrigerant cylinder → VS → VR → VL, the indication of scale.

After stabilizing (calm down), reset the display of scale. (0 display)

Perform a leak check to the dedicated GM.

Note : Use a leak check spray (also liquid detergent is acceptable) to check the gas leak at the connection point.

Pay attention to the posture of the refrigerant cylinder when it is enclosed, depending on the structure of the cylinder, because it may be inverted or sideways.

Valve state

| | open | shut |
|------|------|------|
| VC | | ○ |
| VG | | ○ |
| VL | ○ | |
| VH | | ○ |
| VR | ○ | |
| VS | ○ | |
| main | ○ | |
| SVL | | ○ |

4.4 Starting gas charge

Procedure: Open SVL slowly and check that the display value of scale slowly decreases.

When the specified amount approaches, start slowly closing the SVL, and when reaching the specified amount, shut the SVL. (Please note that there is a time lag between SVL closing operation and display of scale responding to it) Record the final charge amount and shut the main valve of the refrigerant cylinder → VS → VR → VL.

Note: Instead of SVL, VS can be opened and shut to adjust the gas filling amount, but the accuracy of filling quantity is bad. In that case, it is necessary to operate the body compressor and suck in the refrigerant from SVL to VS

valve state


| | open | shut |
|------|------------|------|
| VC | | ○ |
| VG | | ○ |
| VL | ○ | |
| VH | | ○ |
| VR | ○ | |
| VS | ○ | |
| main | ○ | |
| SVL | adjustment | |

4.5 Stop gas charge

Preparations : Tools, Explosion-proof vacuum pump

Procedure : As shown in the figure below, let the long hose outdoors and connect the explosion-proof vacuum pump. Activate the pump, open the valves (VC, VG, VL, VR, VS) and suck out the remaining refrigerant of the hose / general purpose GM.

When release of the refrigerant is completed, shut all the valves and stop the pump. Remove the hose from the low pressure service valve.




Note: Be sure to use explosion-proof type when using vacuum pump to suck refrigerant.

Do not use the powerful vacuum pump used for evacuation by mistake.

4.6 Gas leak check

Preparation : Detector

Procedure : Leak check low pressure service valve with charge hose removed using detector. Subsequently, the high pressure side service port and the part welded and connected at the time of repair are likewise leak checked. If there is no abnormality, attach the cap to the low-pressure side service valve and the high-pressure side service port



Note: Do not apply excessive force to the welded portion of piping when installing the rear cover after the work (prevent occurrence of cracks).

■L (low temperature) cooling circuit

The work content on the low temperature side is work on the high temperature side 3.3 After completing vacuum evacuation, work to enclose pentane is added. After that, work from the work of vacuum evacuation ② in the same procedure as the high temperature side.

Describe additional pentane inclusion work below.

4.7 Enclosed in pentane


4.7.1 Pentane collection

Preparation : Pentane charging cylinder, n-pentane, pentane dryer, explosion-proof vacuum pump, charge hose, branching instrument, tool

Procedure1 : Prepare the pentane, reconfirm the enclosed amount and convert the "weight ⇒ scale" using the conversion table.

Attach a pentane dryer to the tip of the valve VCL of the pentane cylinder and connect the tip of the valve VCU to the vacuum pump with a charge hose. (At this time, VCU is "open", VCL is "shut" state)

Operate the vacuum pump and evacuate it to a state where pentane solution can be injected. (About 1 minute)




n-pentane

Note1 : Because pentane vapor is generated, use an explosion-proof vacuum pump.
The pentane dryer has the function of removing the same moisture as the dehydrator.

Procedure2: Dip the tip of the pentane dryer in pentane solution.

With the VCU opened slightly, opening the VCL slowly will cause the pentane to begin to be sucked into the cylinder.

When the specified amount of pentane is sucked into the cylinder, pull out of pentane ends when VCL is "shut". Shut the VCU, shut down the vacuum pump, remove the charge hose from the VCU, and remove the pentane dryer from the VCL.



Note2: When sucking up pentane, in order to prevent over suction, VCL slowly gently opens and shut and adjusts the suction amount.


Remove the pentane solution from the tip of the pentane dryer for a while, so remove the pentane dryer from the pentane cylinder and remove it from the pentane container after all the pentane has fallen.

4.7.2 Charge pentane to refrigeration circuit

Preparations : Powerful vacuum pump, Teflon hose (for pentane)

Procedure: With the VR shut, connect the valve VCL of the pentane cylinder with the refrigerant port of the dedicated GM with a Teflon hose.

Operate the vacuum pump, open the valves VC, VG, VR and evacuate the Teflon hose.



Note 1: In order to prevent corrosion of hose by pentane and flow into service valve, be sure to use Teflon hose for connection between pentane cylinder and dedicated GM.

Procedure 2: Shut the valves VC and VG, then open the valve VL.

Holding the pentane cylinder vertically, slowly open the VCL and aspirate the pentane into the circuit.

At this time, pull in while watching the guide of the cylinder and stop the VCL when it reaches the specified amount (stop position).

Be sure to stop with the pentane remaining in the cylinder.

Confirm that there is no pentane left in Teflon hose and shut all valves.

Remove the Teflon hose and return the remaining pentane to the container.

Valve Operation

| | open | shut |
|-----|-----------------------|-----------------------|
| VC | | <input type="radio"/> |
| VG | | <input type="radio"/> |
| VL | <input type="radio"/> | <input type="radio"/> |
| VH | | <input type="radio"/> |
| VR | <input type="radio"/> | |
| VCU | | <input type="radio"/> |
| VCL | adjustment | |
| SVL | <input type="radio"/> | |

Note2: If the cylinder becomes empty beyond the specified amount, the air will enter the refrigeration circuit, so be sure to stop with the pentane remaining.

4.7.3 After pentane circuit enclosure, just like high temperature side work, 4. Perform gas charge work.

Note 1 : Paste the new service pipe label (red tag) again near the original pasting position.
Leave the service valve without removing it.
This operation is PA only. It is unnecessary for other destinations.

Note 2 : Do not apply excessive force to the welded part of piping (prevention of cracking) when installing the rear cover after the work is finished.





Confirm operation as final confirmation

Turn on the freezer and check the operation







Test Data





MDF-DU502VH(1φ230V 50Hz)







MDF-DU502VH(1φ220V 60Hz)




MDF-DU702VH(1φ230V 50Hz)



MDF-DU702VH(1φ220V 60Hz)



Temperature uniformity - 17points measuring
MDF-DU502VH



Internal Temperature Uniformity (Reference Data)

MDF-DU502VH

<Conditions>

Ambient temperature: 23/30°C

Load: Unloaded

<Distribution data>

Temperature of the cycle in each area (SV=-80°C, air temperature)

Unit:°C

| | | | Ambient temperature 23°C | | | | | | | |
|---------|-----------------------|-------------|--------------------------|---------|-----------------|--------------|-----------|---------|-----------------|--------------|
| | | | 230V,50Hz | | | | 220V,60Hz | | | |
| | | | Maximum | Minimum | Middle of cycle | Differential | Maximum | Minimum | Middle of cycle | Differential |
| (1) | Upper area | Left front | -77.0 | -79.6 | -78.3 | ±1.3 | -76.9 | -79.6 | -78.3 | ±1.4 |
| | | Left back | -77.2 | -80.1 | -78.7 | ±1.5 | -77.1 | -80.0 | -78.6 | ±1.5 |
| | | Right front | -77.4 | -80.0 | -78.7 | ±1.3 | -77.4 | -79.9 | -78.7 | ±1.3 |
| | | Right back | -77.3 | -79.9 | -78.6 | ±1.3 | -77.2 | -79.8 | -78.5 | ±1.3 |
| | | Center | -77.2 | -79.5 | -78.4 | ±1.2 | -77.1 | -79.5 | -78.3 | ±1.2 |
| (6) | Middle area | Left front | -79.8 | -80.9 | -80.4 | ±0.5 | -79.8 | -80.8 | -80.3 | ±0.5 |
| | | Left back | -80.3 | -81.6 | -81.0 | ±0.6 | -80.3 | -81.6 | -81.0 | ±0.6 |
| | | Right front | -79.7 | -81.0 | -80.4 | ±0.6 | -79.8 | -80.8 | -80.3 | ±0.5 |
| | | Right back | -80.2 | -81.7 | -81.0 | ±0.8 | -80.2 | -81.5 | -80.9 | ±0.6 |
| | | Center | -80.4 | -81.1 | -80.8 | ±0.4 | -80.4 | -81.0 | -80.7 | ±0.3 |
| (11) | Bottom area | Left front | -78.7 | -80.6 | -79.7 | ±1.0 | -78.7 | -80.4 | -79.6 | ±0.9 |
| | | Left back | -79.1 | -81.2 | -80.2 | ±1.1 | -79.1 | -81.1 | -80.1 | ±1.0 |
| | | Right front | -79.2 | -81.2 | -80.2 | ±1.0 | -79.2 | -81.1 | -80.2 | ±1.0 |
| | | Right back | -79.2 | -81.5 | -80.4 | ±1.2 | -79.3 | -81.4 | -80.4 | ±1.1 |
| | | Center | -79.1 | -81.3 | -80.2 | ±1.1 | -79.1 | -81.1 | -80.1 | ±1.0 |
| (16) | Center of Upper shelf | -79.3 | -80.1 | -79.7 | ±0.4 | -79.3 | -80.0 | -79.7 | ±0.3 | |
| (17) | Center of Lower shelf | -80.5 | -81.3 | -80.9 | ±0.4 | -80.5 | -81.1 | -80.8 | ±0.3 | |
| Average | | | - | - | -79.8 | - | - | - | -79.8 | - |

Unit:°C

| | | | Ambient temperature 30°C | | | | | | | |
|---------|-----------------------|-------------|--------------------------|---------|-----------------|--------------|-----------|---------|-----------------|--------------|
| | | | 230V,50Hz | | | | 220V,60Hz | | | |
| | | | Maximum | Minimum | Middle of cycle | Differential | Maximum | Minimum | Middle of cycle | Differential |
| (1) | Upper area | Left front | -76.7 | -79.4 | -78.1 | ±1.4 | -76.8 | -79.5 | -78.2 | ±1.3 |
| | | Left back | -76.9 | -79.9 | -78.4 | ±1.5 | -77.0 | -80.0 | -78.5 | ±1.5 |
| | | Right front | -77.2 | -79.8 | -78.5 | ±1.3 | -77.2 | -79.8 | -78.5 | ±1.3 |
| | | Right back | -77.0 | -79.7 | -78.4 | ±1.4 | -77.1 | -79.8 | -78.5 | ±1.4 |
| | | Center | -76.9 | -79.3 | -78.1 | ±1.2 | -77.1 | -79.4 | -78.3 | ±1.2 |
| (6) | Middle area | Left front | -79.8 | -81.1 | -80.5 | ±0.6 | -79.8 | -81.1 | -80.5 | ±0.6 |
| | | Left back | -80.4 | -81.7 | -81.1 | ±0.6 | -80.3 | -81.8 | -81.1 | ±0.8 |
| | | Right front | -79.8 | -81.0 | -80.4 | ±0.6 | -79.7 | -80.9 | -80.3 | ±0.6 |
| | | Right back | -80.3 | -81.7 | -81.0 | ±0.7 | -80.3 | -81.7 | -81.0 | ±0.7 |
| | | Center | -80.5 | -81.2 | -80.9 | ±0.4 | -80.5 | -81.2 | -80.9 | ±0.4 |
| (11) | Bottom area | Left front | -79.1 | -81.2 | -80.2 | ±1.1 | -79.2 | -81.2 | -80.2 | ±1.0 |
| | | Left back | -79.4 | -81.5 | -80.5 | ±1.1 | -79.5 | -81.5 | -80.5 | ±1.0 |
| | | Right front | -79.5 | -81.4 | -80.5 | ±1.0 | -79.6 | -81.4 | -80.5 | ±0.9 |
| | | Right back | -79.6 | -81.8 | -80.7 | ±1.1 | -79.6 | -81.8 | -80.7 | ±1.1 |
| | | Center | -79.5 | -81.6 | -80.6 | ±1.1 | -79.5 | -81.6 | -80.6 | ±1.1 |
| (16) | Center of Upper shelf | -79.3 | -80.0 | -79.7 | ±0.3 | -79.3 | -80.0 | -79.7 | ±0.3 | |
| (17) | Center of Lower shelf | -80.7 | -81.5 | -81.1 | ±0.4 | -80.8 | -81.5 | -81.2 | ±0.3 | |
| Average | | | - | - | -79.9 | - | - | - | -79.9 | - |





Unit:°C





<Distribution data>

Temperature of the cycle in each area (SV=-70°C, air temperature)


Unit:°C

| | | | Ambient temperature 23°C | | | | | | | |
|------|-----------------------|-------------|--------------------------|---------|-----------------|--------------|-----------|---------|-----------------|--------------|
| | | | 230V,50Hz | | | | 220V,60Hz | | | |
| | | | Maximum | Minimum | Middle of cycle | Differential | Maximum | Minimum | Middle of cycle | Differential |
| (1) | Upper area | Left front | -67.3 | -70.4 | -68.9 | ±1.6 | -67.2 | -70.5 | -68.9 | ±1.7 |
| | | Left back | -67.5 | -71.1 | -69.3 | ±1.8 | -67.5 | -71.1 | -69.3 | ±1.8 |
| | | Right front | -67.4 | -70.6 | -69.0 | ±1.6 | -67.5 | -70.6 | -69.1 | ±1.6 |
| | | Right back | -67.4 | -70.6 | -69.0 | ±1.6 | -67.5 | -70.6 | -69.1 | ±1.6 |
| | | Center | -67.3 | -70.6 | -69.0 | ±1.7 | -67.3 | -70.6 | -69.0 | ±1.7 |
| (6) | Middle area | Left front | -69.3 | -71.0 | -70.2 | ±0.9 | -69.3 | -70.9 | -70.1 | ±0.8 |
| | | Left back | -70.0 | -71.7 | -70.9 | ±0.9 | -69.9 | -71.8 | -70.9 | ±0.9 |
| | | Right front | -69.2 | -70.8 | -70.0 | ±0.8 | -69.2 | -70.8 | -70.0 | ±0.8 |
| | | Right back | -70.0 | -71.9 | -71.0 | ±1.0 | -69.9 | -71.9 | -70.9 | ±1.0 |
| | | Center | -70.1 | -71.1 | -70.6 | ±0.5 | -70.2 | -71.1 | -70.7 | ±0.5 |
| (11) | Bottom area | Left front | -67.4 | -69.6 | -68.5 | ±1.1 | -67.4 | -69.5 | -68.5 | ±1.1 |
| | | Left back | -68.0 | -70.3 | -69.2 | ±1.2 | -68.0 | -70.2 | -69.1 | ±1.1 |
| | | Right front | -67.7 | -69.8 | -68.8 | ±1.1 | -67.6 | -69.6 | -68.6 | ±1.0 |
| | | Right back | -68.1 | -70.6 | -69.4 | ±1.3 | -68.1 | -70.5 | -69.3 | ±1.2 |
| | | Center | -67.6 | -69.8 | -68.7 | ±1.1 | -67.6 | -69.6 | -68.6 | ±1.0 |
| (16) | Center of Upper shelf | -69.5 | -70.6 | -70.1 | ±0.6 | -69.5 | -70.5 | -70.0 | ±0.5 | |
| (17) | Center | | | | | | | | | |





Temperature uniformity - 17points measuring
MDF-DU702VH



Internal Temperature Uniformity (Reference Data)

MDF-DU702VH

<Conditions>

Ambient temperature: 23/30°C

Load: Unloaded

<Distribution data>

Temperature of the cycle in each area (SV=-80°C, air temperature)

Unit:°C

| | | | Ambient temperature 23°C | | | | | | | |
|---------|-----------------------|-------------|--------------------------|---------|-----------------|--------------|-----------|---------|-----------------|--------------|
| | | | 230V 50Hz | | | | 220V 60Hz | | | |
| | | | Maximum | Minimum | Middle of cycle | Differential | Maximum | Minimum | Middle of cycle | Differential |
| ① | Upper area | Left front | -76.1 | -79.2 | -77.7 | ±1.6 | -76.1 | -79.2 | -77.7 | ±1.6 |
| ② | | Left back | -76.7 | -79.9 | -78.3 | ±1.6 | -76.7 | -79.9 | -78.3 | ±1.6 |
| ③ | | Right front | -76.7 | -79.7 | -78.2 | ±1.5 | -76.7 | -79.7 | -78.2 | ±1.5 |
| ④ | | Right back | -77.0 | -80.2 | -78.6 | ±1.6 | -77.1 | -80.1 | -78.6 | ±1.5 |
| ⑤ | | Center | -76.5 | -79.5 | -78.0 | ±1.5 | -76.5 | -79.5 | -78.0 | ±1.5 |
| ⑥ | Middle area | Left front | -78.9 | -80.2 | -79.6 | ±0.6 | -78.9 | -80.2 | -79.6 | ±0.6 |
| ⑦ | | Left back | -79.1 | -80.6 | -79.9 | ±0.8 | -79.0 | -80.7 | -79.9 | ±0.9 |
| ⑧ | | Right front | -78.9 | -80.0 | -79.5 | ±0.5 | -78.9 | -80.0 | -79.5 | ±0.5 |
| ⑨ | | Right back | -79.2 | -80.9 | -80.1 | ±0.9 | -79.3 | -81.0 | -80.2 | ±0.8 |
| ⑩ | | Center | -79.4 | -80.5 | -80.0 | ±0.5 | -79.4 | -80.5 | -80.0 | ±0.5 |
| ⑪ | Bottom area | Left front | -78.0 | -80.6 | -79.3 | ±1.3 | -78.0 | -80.6 | -79.3 | ±1.3 |
| ⑫ | | Left back | -78.3 | -80.8 | -79.6 | ±1.3 | -78.3 | -80.8 | -79.6 | ±1.3 |
| ⑬ | | Right front | -78.1 | -80.5 | -79.3 | ±1.2 | -78.0 | -80.5 | -79.3 | ±1.3 |
| ⑭ | | Right back | -78.6 | -81.3 | -80.0 | ±1.4 | -78.6 | -81.3 | -80.0 | ±1.4 |
| ⑮ | | Center | -78.5 | -81.0 | -79.8 | ±1.3 | -78.5 | -81.0 | -79.8 | ±1.3 |
| ⑯ | Center of Upper shelf | | -78.9 | -79.8 | -79.4 | ±0.5 | -78.9 | -79.8 | -79.4 | ±0.5 |
| ⑰ | Center of Lower shelf | | -79.8 | -80.6 | -80.2 | ±0.4 | -79.9 | -80.6 | -80.3 | ±0.4 |
| Average | | | - | - | -79.2 | - | - | - | -79.2 | - |

Unit:°C

| | | | Ambient temperature 30°C | | | | | | | |
|---------|-----------------------|-------------|--------------------------|---------|-----------------|--------------|-----------|---------|-----------------|--------------|
| | | | 230V 50Hz | | | | 220V 60Hz | | | |
| | | | Maximum | Minimum | Middle of cycle | Differential | Maximum | Minimum | Middle of cycle | Differential |
| ① | Upper area | Left front | -75.9 | -78.9 | -77.4 | ±1.5 | -76.0 | -78.9 | -77.5 | ±1.5 |
| ② | | Left back | -76.6 | -79.8 | -78.2 | ±1.6 | -76.6 | -79.7 | -78.2 | ±1.6 |
| ③ | | Right front | -76.6 | -79.4 | -78.0 | ±1.4 | -76.6 | -79.4 | -78.0 | ±1.4 |
| ④ | | Right back | -77.0 | -80.0 | -78.5 | ±1.5 | -77.0 | -79.9 | -78.5 | ±1.5 |
| ⑤ | | Center | -76.4 | -79.1 | -77.8 | ±1.4 | -76.3 | -79.1 | -77.7 | ±1.4 |
| ⑥ | Middle area | Left front | -79.0 | -80.2 | -79.6 | ±0.6 | -78.9 | -80.2 | -79.6 | ±0.6 |
| ⑦ | | Left back | -79.1 | -80.6 | -79.9 | ±0.8 | -79.1 | -80.5 | -79.8 | ±0.7 |
| ⑧ | | Right front | -78.9 | -80.0 | -79.5 | ±0.5 | -78.9 | -80.0 | -79.5 | ±0.5 |
| ⑨ | | Right back | -79.3 | -80.8 | -80.1 | ±0.8 | -79.3 | -80.9 | -80.1 | ±0.8 |
| ⑩ | | Center | -79.5 | -80.4 | -80.0 | ±0.5 | -79.5 | -80.5 | -80.0 | ±0.5 |
| ⑪ | Bottom area | Left front | -78.1 | -80.6 | -79.4 | ±1.3 | -78.1 | -80.6 | -79.4 | ±1.3 |
| ⑫ | | Left back | -78.4 | -80.7 | -79.6 | ±1.2 | -78.4 | -80.7 | -79.6 | ±1.2 |
| ⑬ | | Right front | -78.2 | -80.6 | -79.4 | ±1.2 | -78.2 | -80.6 | -79.4 | ±1.2 |
| ⑭ | | Right back | -78.8 | -81.3 | -80.1 | ±1.3 | -78.8 | -81.3 | -80.1 | ±1.3 |
| ⑮ | | Center | -78.7 | -81.1 | -79.9 | ±1.2 | -78.7 | -81.1 | -79.9 | ±1.2 |
| ⑯ | Center of Upper shelf | | -78.9 | -79.7 | -79.3 | ±0.4 | -78.9 | -79.7 | -79.3 | ±0.4 |
| ⑰ | Center of Lower shelf | | -80.0 | -80.7 | -80.4 | ±0.4 | -79.9 | -80.7 | -80.3 | ±0.4 |
| Average | | | - | - | -79.2 | - | - | - | -79.2 | - |

Unit:°C





<Distribution data>





Temperature of the cycle in each area (SV=-70°C, air temperature)

Unit:°C

| | | | Ambient temperature 23°C | | | | | | | |
|---|-----------------------|-------------|--------------------------|---------|-----------------|--------------|-----------|---------|-----------------|--------------|
| | | | 230V 50Hz | | | | 220V 60Hz | | | |
| | | | Maximum | Minimum | Middle of cycle | Differential | Maximum | Minimum | Middle of cycle | Differential |
| ① | Upper area | Left front | -66.3 | -70.4 | -68.4 | ±2.1 | -66.3 | -70.3 | -68.3 | ±2.0 |
| ② | | Left back | -67.2 | -71.5 | -69.4 | ±2.2 | -67.2 | -71.5 | -69.4 | ±2.2 |
| ③ | | Right front | -66.8 | -70.5 | -68.7 | ±1.9 | -66.8 | -70.4 | -68.6 | ±1.8 |
| ④ | | Right back | -67.4 | -71.6 | -69.5 | ±2.1 | -67.5 | -71.5 | -69.5 | ±2.0 |
| ⑤ | | Center | -67.0 | -70.6 | -68.8 | ±1.8 | -67.0 | -70.5 | -68.8 | ±1.8 |
| ⑥ | Middle area | Left front | -68.4 | -69.8 | -69.1 | ±0.7 | -68.4 | -69.8 | -69.1 | ±0.7 |
| ⑦ | | Left back | -69.0 | -71.1 | -70.1 | ±1.1 | -69.0 | -71.1 | -70.1 | ±1.1 |
| ⑧ | | Right front | -68.3 | -69.7 | -69.0 | ±0.7 | -68.3 | -69.7 | -69.0 | ±0.7 |
| ⑨ | | Right back | -69.1 | -71.4 | -70.3 | ±1.2 | -69.1 | -71.3 | -70.2 | ±1.1 |
| ⑩ | | Center | -69.4 | -70.6 | -70.0 | ±0.6 | -69.4 | -70.6 | -70.0 | ±0.6 |
| ⑪ | Bottom area | Left front | -66.6 | -69.5 | -68.1 | ±1.5 | -66.6 | -69.4 | -68.0 | ±1.4 |
| ⑫ | | Left back | -67.5 | -70.3 | -68.9 | ±1.4 | -67.5 | -70.2 | -68.9 | ±1.4 |
| ⑬ | | Right front | -66.5 | -69.3 | -67.9 | ±1.4 | -66.5 | -69.3 | -67.9 | ±1.4 |
| ⑭ | | Right back | -67.4 | -70.8 | -69.1 | ±1.7 | -67.3 | -70.7 | -69.0 | ±1.7 |
| ⑮ | | Center | -67.2 | -69.9 | -68.6 | ±1.4 | -67.2 | -69.8 | -68.5 | ±1.3 |
| ⑯ | Center of Upper shelf | | -69.3 | -70.8 | -70.1 | ±0.8 | -69.3 | -70.7 | -70.0 | ±0.7 |
| ⑰ | Center of Lower shelf | | -69.1 | | | | | | | |

Running cycle
MDF-DU702VH














Backup cooling kit installation and setting procedure

Confirmation of packing items

① Before installing the optional component, check the accessories listed in table 1 below are enclosed.

Table 1

| No. | PARTS NAME | Q'TY. | APPEARANCE | THE EXPLANATION OF THE USE |
|-----|---|-------|---|--|
| 1 | CO2 solenoid valve Assy. | 1 |  | Electromagnetic valve for injecting the LCO ₂ gas into the freezer chamber (put together at the installation of MDF-UB7). |
| 2 | Back-up cooling kit body | 1 |  | For controlling the LCO ₂ gas injection to prevent the temperature rise in the event of power failure. |
| 3 | Battery | 1 |  | Power source for backup board. ※Never touch the terminal portion. That may cause electric shock. DC 6 V x 4 pcs. = DC 24 V |
| 4 | Connecting pipe and joint packing | 1 |  | Pipe for connecting the LCO ₂ gas cylinder with pipe joint. ※Metric screw type |
| 5 | Connecting pipe and joint packing | 1 |  | Pipe for connecting the LCO ₂ gas cylinder with pipe joint. ※Inch screw type |
| 6 | Nozzle cover attachment | 1 |  | For fixing the nozzle cover. |
| 7 | Nozzle cover | 1 |  | To protect the nozzle in the chamber. |
| 8 | CO ₂ injection nozzle and thermal insulation | 1 |  | Component for CO ₂ solenoid valve assy. |
| 9 | Solenoid cover | 1 |  | To protect the CO ₂ solenoid valve assy. |

| | | | | |
|----|---|----|---|---|
| 10 | Solenoid valve cover attachment | 1 |  | For fixing the solenoid valve cover. |
| 11 | Binder | 2 |  | For fixing the thermal insulation. |
| 12 | Truss screw M4x10 stainless | 12 |  | For fixing the back-up cooling kit body(4). For fixing the battery(4). For fixing the nozzle cover(2). For fixing the solenoid valve cover(2). |
| 13 | Truss tapping screw M4x12 stainless | 16 |  | For fixing the CO ₂ solenoid valve assy(4). For fixing the solenoid valve cover(4). For fixing the valve cover attachment(2). For fixing the nylon clip(6). |
| 14 | Nylon clip 6N | 1 |  | For fixing the thermostat sensor of the back-up cooling kit body. |
| 15 | Nylon clip 2N | 5 |  | For fixing the thermostat sensor of the back-up cooling kit body. |

BACKUP COOLING KIT INSTALLATION

- ② Unravel the bundled harness

Harness UB_BATTRY_1

R,BL 3P-3P

Harness MAIN_UBK

BL,R,OR.BU.Y.W 6P-6P,3P

Harness UB_VALVE_1

Y,G 2P-2P



Fig.2-1

- ③ Remove the product unit cover R and the switch box cover



Fig.3-1





Fig.3-2

- ④ While lightly inserting the back-up cooling kit body into the installation position of the product, wire harnesses through the wire saddle and wire as shown in the figure.



Wire harnesses through wire saddles as shown in the figure.

- Wire harness UB_BATTRY_1 through wire saddle ② through wire saddle ①.
- Wire harness MAIN_UBK is passed through the wire saddles ⑤ to ⑥ in order, and connected to the CN 14 of the main board.



Wire harness UB_VALVE_1 in the order of wire saddle ③ → ④ → ⑤ → ⑥ → ⑦, and connect it to the harness UB_VALVE_2 that comes out from the bottom of the cabinet.

- ⑤ • Wire the gas sensor as shown by the dotted line and set it free toward the upper part of the power supply box.

Fix the four marking parts in the right figure with the truss head screws M4x6.




Fig.5-1

⑥ Wire the gas sensor cable using the clamps on the marking section on the bellow.

Remove CAP A attached to the base and wire the cable of the gas sensor from the hole of the base to the back of the body.



Fig.6-1






Fig.6-2



⑦ Connect the harness UB_BATTRY_2 to the position marked with the square in the right figure.

In the figure on the right, connect the red wire to the solid square marker portion and the black wire to the dotted square marker portion respectively.


Fix the four marking parts in the right figure with the truss head screws M4x6.



⑧ Route UB_BATTRY_2 through the two clamps on the right figure.



⑨ Connect UB_BATTRY_1 and UB_BATTRY_2.



- ⑩ Remove the sensor outlet cover on the back and the temperature control sensor cover inside the chamber, and remove the silicone in the backup sensor insertion hole.

Insert the gas thermo into the backup sensor insertion hole from which silicon has been removed and fix the gas thermo using the nylon clip attached to the sensor mounting plate.

Fill the gap between sensor insertion hole and gas thermo with silicon.

Reinstall the temperature control sensor cover and the sensor outlet cover.




Fig.10-1 Fig.10-2




Fig.10-3 Fig.10-4

seal with silicon

- ⑪ Attach the nylon clips SL - 2 N and SL - 6 N to the holes for those on the back and secure the gas thermo piping with them. Insert a notch in the cap A that you removed, insert a gas thermo piping into it, and push A into the hole in the base so as to close the hole.




Fig.11-1 Fig.11-2

- ⑫ Remove the seal on the back side and take out the stored harness PT_AIP_UB_VALVE_ 2.





Fig.12-1 PT_AIP_UB_VALVE_2
Fig.12-2 solenoid valve power

Fig.12-3 insert hole for gas tube


- ⑬ Insert the valve pipe from the back and secure it with the mounting plate.
 Bend the valve tip downward at 90 degree.
 After inserting the gas pipe, fill the gap with silicon.
 When the cover is attached, the tip of the valve pipe does not interfere with the cover.




- ⑭ Fix the solenoid valve cover to the position shown in the figure with 2 truss tapping screws M4x12.



- ⑮ Cover the backup cover so that the various harnesses are not caught in, and fix the four places with screws.





WIRING



PRECAUTIONS

Before installation


- Liquid CO₂ cylinder should be a siphon type cylinder
 - Only liquid CO₂ can be used for this system. Never use other kinds of gases.
 - Liquid CO₂ loses its cooling capacity when the ambient temperature is above +31°C. Do not install the liquid CO₂ cylinder in a place exposed to direct sunlight or near heat-emitting appliances such as stoves and heaters. The back-up time per 1 liquid CO₂ cylinder depends on the ambient temperature. Install a liquid CO₂ cylinder in the cool location.
 - Liquid CO₂ cylinder should be installed in a special cylinder stand (see figure). Never install the cylinder without the stand. If you have difficulty in obtaining the stand, contact your dealer.
 - The length of time in which a 30kg liquid CO₂ cylinder can maintain a temperature of -70°C is shown in the another sheet "Ambient Temperature vs. -70°C Holding Time of Liquid CO₂ cylinder (30kg)".
- Install the required number of cylinders according to the customer's desired back-up time.



- Before installation, tell the customer about cylinder storage and exchange procedures at the time of cylinder inspection according to high-pressure gas regulations.

At installation

- The connecting pipe between the liquid CO₂ cylinder and the backup cooling kit should not be installed in a place where it may be bumped in a passage.
- Use the connecting pipe provided with the backup cooling kit. Run the connecting pipe in a coil or U-bend formation to act as a cushion (see figure).
- Do not let foreign particles or water get into the pipe or joints.






After installation

- Check that liquid CO₂ is injected into the freezer chamber. Refer to the "Injection Test" for an explanation of the testing procedure.
- Conduct a leakage test of the liquid CO₂ cylinder and connecting area.

INSTALLATION

1. Put the liquid CO₂ cylinder into the CO₂ cylinder stand. Secure each cylinder to the stand with chains.
2. Connect the valve on the CO₂ cylinder stand with the CO₂ outlet of the cylinder with a flexible hose.
3. Close the valve on the CO₂ cylinder stand completely and open the cylinder valve. Then check for CO₂ leakage at the cylinder valve, cylinder stand valve and the connecting area of the flexible hose.
4. When no leakage is found, close the cylinder valve.
5. Cut the connecting pipe (total length; 2.5m) enclosed with the backup cooling kit into an appropriate length. The cut end should be flared.
6. Remove the cap covering the connecting port on the backup cooling kit. Place the flared nut, ring, and sleeve onto the end of the pipe.
7. Insert the cut side of the connecting pipe all the way into the connecting port on the backup cooling kit and tighten the flared nut.
8. Connect the other side of the connecting pipe with the cylinder valve (when 1 cylinder) or the cylinder stand valve (when more than 1 cylinder). Tighten the flared nut all the way.
- Note: Run the connecting pipe in a coil or U-bend formation to act as a cushion (see Fig. 2).
9. Open the cylinder valve or cylinder stand valve (whichever applicable) and check for leakage at the both ends of the connecting pipe. The leakage can be checked by applying soapy water to the connecting area as shown in Fig. 3.
10. Close all of the cylinder valve and also the cylinder stand valve.
11. Check the injection by following the procedure in the "Injection Test" on the next page.



INJECTION TEST

1. Run the freezer until the chamber temperature reaches the appropriate level.
2. With the power switch of the backup cooling kit turned off, open the cylinder valve and the cylinder stand valve.
3. Set the temperature of the backup cooling kit to a temperature 10°C or more higher than the freezer chamber temperature.

Note:

The temperature setting range of the backup cooling kit is between -50 and -70°C. The consumption of liquid CO₂ gas is increased and the backup cooling time is decreased considerably if the set temperature is lower than -70°C.


4. Turn on the power switch of the backup cooling kit with the freezer door closed.
For the up-right freezer, check the indicator lamp (green) on the door switch is ON.
5. Press the test switch for the backup cooling kit to check the liquid CO₂ gas is injected into the freezer chamber.
With the back-up system set in this way, the liquid CO₂ gas is injected into the freezer chamber automatically to ensure freezing when the chamber temperature reaches the set temperature of the backup cooling kit.
6. If the set temperature of the backup cooling kit is below -60°C, turn off the power switch of the freezer.
7. Check the liquid CO₂ gas is injected into the freezer chamber automatically when the chamber temperature rises.
Note the backup cooling time is decrease if the liquid CO₂ gas is injected continuously.
8. Ensure the injection of the liquid CO₂ gas automatically stops when the chamber temperature goes down.

●Refer to the "Operating Instructions" for an explanation of the switches on the freezer.

SETTING



①Press Menu key at lower right of TOP display for 5 second




②Move to Service Code display. Input “335232” with number key at right of the display. And press OK



③Move to Menu display, press Service #2 key at the lower middle of the display.




④Move to Service #2 display, press System Configuration #2 key at upper right of the display.



⑤Move to System Configuration #2. Press and slide the slide bar at right of Backup Function item at upper of the display, (Then the slide bar change “ON”). Make sure that the slide bar change and press Apply key.





⑥Move to Service display, press TOP key at upper right of the display.



⑦Move to TOP display, make sure that display at Backup. Now this operation is completed.

MDF-DU500VX

AT23°C no load SV=70°C



| No. | thermo couple | Max. temp (°C) | Min. temp. (°C) | middle value (°C) |
|-----|---------------|-------------------|--------------------|----------------------|
| 1 | top stage | -55.1 | -66.6 | -60.9 |
| 2 | 2nd stage | -58.2 | -71.5 | -64.9 |
| 3 | 3rd stage | -57.8 | -72.4 | -65.1 |
| 4 | 4th stage | -58.2 | -72.1 | -65.2 |
| 5 | gas sensor | -58 | -66.8 | -62.4 |

***Revision history**

| Revision # | Date of revision | Reason of revision |
|--------------|------------------|--|
| SM0000002-01 | Oct.17. 2017 | Update Service manual for released catalogue. |
| SM0000002-02 | Oct.27. 2017 | Change firmware version to V1.04A for main PCB to fix the malfunction for optional backup system, MDF-UB7. |
| SM0000002-03 | Feb.7. 2018 | Update PL for MDF-DU502VH release. |
| | | |
| | | |
| | | |
| | | |